ACCENT-300 ALAT

Nr kat. **7-316** (PL)

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Cat. No **7-316** (EN)

INTENDED USE


Diagnostic kit for determination of alanine aminotransferase activity, intended to use in automatic analyzer ACCENT-300. The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Alanine aminotransferase (ALAT, ALT, GPT) is an enzyme participated in amino acids metabolism. ALAT is present in all tissues but the highest level is found in liver and kidney cells. Damage of hepatocytes or kidney cells causes significant release of ALAT into the circulation. Measurement of ALT activity in serum is valuable in the diagnosis of liver diseases: jaundice, mononucleosis or hepatic cirrhosis.

METHOD PRINCIPLE

Optimized, modified method according to International Federation of Clinical Chemistry (IFCC), without pyridoxal phosphate.

The rate of absorbance changing at $\lambda=340$ nm is directly proportional to alanine aminotransferase activity.

REAGENTS

Package	
1-Reagent	4 x 40 ml
2-Reagent	2 x 20.5 ml

The reagents, stored at 2-8°C are stable up to expiry date printed on the package. The reagents stored on board of the analyzer at 2-10°C are stable for 12 weeks.

Concentrations in the test

1-Reagent	
Tris buffer, pH 7,4	≤ 150 mmol/l
L-alanine	≤ 750 mmol/l
LDH	≤ 4 U/ml
stabilizer, preservative	
2-Reagent	
2-oxoglutarate	≤ 74 mmol/l
NADH	≤ 1.7 mmol/l
buffer	
preservatives	

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN

Serum, heparinized or EDTA plasma free from hemolysis. Hemolysis should be avoided, since ALAT activity in erythrocytes is 3 to 5 times higher than in normal serum. Do not freeze the samples. ALAT activity remains stable in specimen up to 3 days at 15-25°C or up to 7 days at 2-8°C. Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.
Deionised water is recommended as a reagent blank.

REFERENCE VALUES ⁶

serum / plasma	37°C	
women	up to 31 U/l	up to 0.517 μ kat/l
men	up to 41 U/l	up to 0.683 μ kat/l

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use with each batch of samples the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173). For the calibration of automatic analysers the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended. Deionised water should be used as a calibrator 0. The calibration curve should be prepared every 12 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analyser ACCENT-300. Results may vary if a different instrument or a manual procedure is used.

- LOQ (Limit of Quantitation):** 7 U/l (0.12 μ kat/l)
- Linearity:** up to 1000 U/l (16.7 μ kat/l)
- Specificity / Interferences**
Haemoglobin up to 0.31 g/dl, ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl and triglycerides up to 1000 mg/dl do not interfere with the test.

Precision

Repeatability (run to run) n=20	Mean	SD	CV
	[U/l]	[U/l]	[%]
level 1	31.9	1.7	5.4
level 2	107	1.3	1.2
Reproducibility (day to day) n=80	Mean	SD	CV
	[U/l]	[U/l]	[%]
level 1	33.1	1.5	4.7
level 2	104	2.3	2.2

Method comparison

A comparison between ALAT values determined at ACCENT-300 (y) and at BECKMAN COULTER AU680 (x) using 120 serum samples gave following results:
 $y = 0.9953x + 0.8962$ U/l;
 $R = 0.999$ (R – correlation coefficient)

WASTE MANAGEMENT

Please refer to local legal requirements.

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