## ACCENT-300 TOTAL IgE

Nr kat. 7-331 (PL)

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(EN)

### INTENDED USE

Diagnostic kit for determination of IgE levels, used in automatic analyser ACCENT-300.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

IgE is an immunoglobulin with a molecular weight of approximately 190 kD normally present in the blood in trace amounts. Continual production of IgE antibodies in response to common naturally occurring allergens, however, often results in elevated serum levels and in the development of such clinically important Type I allergic reactions as asthma, hay fever, dermatitis and food allergies. Elevated IgE levels are also seen in parasitic (helminth) diseases, IgE myeloma, and in hepatitis. The measurement of IgE in human serum is thus considered to be useful in the diagnosis, treatment, assessment of disease progression, or postoperative prognosis for such conditions.

### METHOD PRINCIPLE

When an antigen-antibody reaction occurs between IgE in a sample and anti-IgE antibody which has been sensitized to latex particles, agglutination results. This agglutination is detected as an absorbance change (572 nm), with the magnitude of the change being proportional to the quantity of IgE in the sample. The actual concentration is then determined by interpolation from a calibration curve prepared from calibrators of known concentration.

### REAGENTS

#### Package

1-Reagent 1 x 33 ml  
2-Reagent 1 x 17 ml

The reagents, stored at 2-10°C are stable up to expiry date printed on the package. The reagents stored on board of the analyser at 2-10°C are stable for 5 weeks.

### Concentrations in the test

suspension of latex particles sensitized with (mouse) anti-IgE antibodies (pH 7.3) 0.125 w/v%  
glycine buffer solution (pH 8.3)  
preservative

### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- After measurements are taken, reagent bottles should be capped and kept at 2-10°C.
- Do not interchange the caps of reagent bottles.
- Reagents with different lot numbers should not be interchanged or mixed.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

### SPECIMEN

Serum, plasma collected on heparin (sodium / lithium salt), EDTA (sodium / potassium salt) or citric acid.

If the test cannot be done immediately, the sample should be stored at -20°C. Avoid repeated freezing and thawing. Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-Reagent and 2-Reagent are ready to use.  
0.9% NaCl is recommended as a reagent blank.

### REFERENCE VALUES <sup>3</sup>

serum, plasma	< 358 IU/ml
---------------	-------------

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of sample the CORMAY IMMUNO-CONTROL II (Cat. No 4-290) s.

For the calibration of automatic analysers systems the CORMAY IgE CALIBRATORS kit (Cat. No 4-280) is recommended.

The calibration curve should be prepared every 5 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

These metrological characteristics have been obtained using an automatic analyser TBA-30R. Results may vary if a different instrument is used.

- Analytical range:** 25 – 1000 IU/ml.

For higher concentration, dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

- Specificity / Interferences**

Haemoglobin up to 0.5 g/dl, bilirubin up to 30 mg/dl, triglycerides up to 1500 mg/dl, RF up to 500 IU/ml do not interfere with the test.

- Precision**

Repeatability (run to run) n = 10	Mean [IU/ml]	SD [IU/ml]	CV [%]
level 1	40.5	2.7	6.57
level 2	427.4	7.7	1.80

- Method comparison**

A comparison between CORMAY reagent (y) and another commercially available assay (x) using 55 samples gave following results:

$$y = 1.01x + 11.7 \text{ IU/ml;}$$

$$R = 0.9967 \quad (R - \text{correlation coefficient})$$

### WASTE MANAGEMENT

Please refer to local legal requirements.

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