## ACCENT-300 RF

Nr kat. 7-337 (PL)

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Cat. No 7-337 (EN)

### INTENDED USE

Diagnostic kit for determination of rheumatoid factor levels, used in automatic analyser ACCENT-300.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Rheumatoid factor (RF) is an autoantibody against human IgG commonly seen in sera at a high concentration in some conditions, particularly in patients with rheumatoid arthritis. The measurement of RF value is useful in evaluating the diagnosis, effects of therapy and prognosis of RA, systemic lupus erythematosus, chronic hepatopathy, etc.

### METHOD PRINCIPLE

When an antigen-antibody reaction occurs between RF in a sample and denatured human IgG which has been sensitized to latex particles, agglutination results. This agglutination is detected as an absorbance change (572 nm), with the magnitude of the change being proportional to the quantity of RF in the sample. The actual concentration is then determined by interpolation from a calibration curve prepared from calibrators of known concentration.

### REAGENTS

#### Package

1-Reagent 1 x 40 ml  
 2-Reagent 1 x 12.5 ml

The reagents when stored at 2-10°C are stable up to expiry date printed on the package. The reagents stored on board of the analyser at 2-10°C are stable for 12 weeks.

### Concentrations in the test

suspension of latex particles sensitized with denatured human IgG (pH 7.3) 0.17 w/v%  
 glycine buffer solution (pH 8.3)  
 preservative

### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- After measurements are taken, reagent bottles should be capped and kept at 2-10°C.
- Do not interchange the caps of reagent bottles.
- Reagents with different lot numbers should not be interchanged or mixed.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

### SPECIMEN

Serum or plasma (Na-EDTA, K-EDTA, Na-Heparin, Li-Heparin, citric acid).

If the test cannot be done immediately, the sample should be placed in a tightly sealable container and stored at -20°C. Repeated freezing and thawing should be avoided.

Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-Reagent and 2-Reagent are ready to use.

0.9% NaCl is recommended as a reagent blank.

#### Actions required:

When performing assays in analyser ACCENT-300, there is a probability of **cross-contamination** affecting the tests results: RF – CALCIUM. To avoid this effect follow the recommendations contained in the instruction 51\_03\_24\_007\_ACCENT-300\_CARRYOVER.

### REFERENCE VALUES <sup>4</sup>

serum, plasma	< 18 IU/ml
---------------	------------

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples the CORMAY IMMUNO-CONTROL I (Cat. No 4-288).

For the calibration of automatic analysers systems the CORMAY RF CALBRATORS kit (Cat. No 4-277) is recommended. As a 0 calibrator 0.9% NaCl should be used.

The calibration curve should be prepared every 5 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using an automatic analysers HITACHI 917 and ACCENT-300. Results may vary if a different instrument is used.

- Sensitivity:** 12 IU/ml.
- Linearity:** up to 120 IU/ml.

For higher concentration, dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

#### Specificity / Interferences

Haemoglobin up to 1.0 g/dl, bilirubin up to 66 mg/dl, intralipid up to 20 g/l do not interfere with the test.

#### Precision

Repeatability (run to run) n = 21	Mean [IU/ml]	SD [IU/ml]	CV [%]
level 1	22.0	0.19	0.85
level 2	49.3	0.26	0.52
Reproducibility (day to day) n = 63	Mean [IU/ml]	SD [IU/ml]	CV [%]
level 1	20.0	1.03	5.12
level 2	48.0	0.89	1.86

#### Method comparison

A comparison between CORMAY reagent (y) and commercially available assay (x) using 60 samples gave following results:

$$y = 0.946x - 0.87 \text{ IU/ml};$$

$$R = 0.993 \quad (R - \text{correlation coefficient})$$

### WASTE MANAGEMENT

Please refer to local legal requirements.

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### PROGRAM NA ANALIZATOR / APPLICATION / АДАПТАЦИЯ:

#### Parameters

No.	34	Prim.Wave.	578
Test	RF	Sec.Wave.	
Method	Fixed-time	Sample Vol.	7.5
Direction	Ascend	R1 Vol.	240
Unit	IU/ml	R2 Vol.	70
Decimals	1	Line. Limit	20
Incubation	20	Antigen Check	
Reaction	5   11	Substrat	0
<b>R1 Blank</b>		<b>Mix. R Blank</b>	
Lower	0	Lower	0
Upper	0	Upper	0
<b>Response</b>		<b>Linearity</b>	
Lower	-2.5	Lower	12
Upper	2.5	Upper	120
Sample Vol.	45	Full Name	RF
Dilution	5	Print No.	34

#### Calibration

Rule	Logistic 5P
K Factor	0
Replicates	1
Interval	35
Sensitivity	0
Correlation	0
Difference	2.5
Blank Response	0   2.5
Coefficient	0
Difference	
Non-linear SD	0