## ACCENT-300 LDH

Nr kat. 7-339

(PL)

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Cat. No **7-339** (EN)

### INTENDED USE

Diagnostic kit for determination of lactate dehydrogenase activity, used in automatic analyser ACCENT-300.


The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Lactate dehydrogenase (LDH, LD) is intracellular enzyme occurred in all tissues. LDH catalyzes the reversible conversion of lactate to pyruvate using NAD<sup>+</sup> as a cofactor. LD is a tetramer containing two possible forms of subunits: H and M. The result is five isoenzymes termed LD-1 (H<sub>4</sub>) through LD-5 (M<sub>4</sub>). The isoenzymes are present in different proportion in each tissue and have different electrophoresis mobility, what is very useful for diagnostic.

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The rate of absorbance changing at  $\lambda=340$  nm is directly proportional to lactate dehydrogenase activity.

### REAGENTS

<b>Package</b>	
1-Reagent	2 x 40 ml
2-Reagent	1 x 20.5 ml

The reagents, stored at 2-8°C are stable up to expiry date printed on the package. Stability on board of the analyser at 2-10°C: 10 weeks.

### Concentrations in the test

phosphate buffer (pH 7.5)	50 mmol/l
pyruvate	0.6 mmol/l
NADH	0.25 mmol/l
preservative	

### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

### SPECIMEN

Serum, heparinized plasma free from hemolysis. Do not use hemolyzed blood or serum because erythrocytes contain 150 times more LDH activity than serum.

As an anticoagulant for plasma preparation use heparin lithium or ammonium salt.

LDH activity is unstable and is rapidly lost during storage. Specimens can be stored up to 4 hours at 15-25°C or 1-2 days at 2-8°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-Reagent and 2-Reagent are ready to use. Deionised water is recommended as a reagent blank.

### REFERENCE VALUES <sup>4</sup>

serum / plasma	37°C	37°C
adults	225 – 450 U/l	3.75 – 7.50 $\mu$ kat/l

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) with each batch of samples.

For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and/or LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

The calibration curve should be prepared every 8 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

These metrological characteristics have been obtained using automatic analyser ACCENT-300. Results may vary if a different instrument or a manual procedure is used.

- Sensitivity:** 38 U/l (0.633  $\mu$ kat/l).
- Linearity:** up to 2000 U/l (33.3  $\mu$ kat/l).

If LDH activity in tested sample exceeds 2000 U/l dilute the sample with 0.9% NaCl in the ratio of 1 to 9 and repeat the assay, multiply the result by 10.

- Specificity / Interferences**

Bilirubin up to 20 mg/dl, haemoglobin up to 5 g/dl, ascorbate up to 62 mg/l and triglycerides up to 1000 mg/dl do not interfere with the test.

- Precision**

Repeatability (run to run) n = 20	Mean	SD	CV
	[U/l]	[U/l]	[%]
level 1	380.70	6.21	1.63
level 2	762.55	3.97	0.52
Reproducibility (day to day) n = 80	Mean	SD	CV
	[U/l]	[U/l]	[%]
level 1	434.14	10.88	2.51
level 2	719.59	12.44	1.73

- Method comparison**

A comparison between LDH values determined at ACCENT-300 (y) and at COBAS INTEGRA 400 (x) using 96 samples gave following results:

$$y = 0.8658x + 58.87 \text{ U/l};$$

$$R = 0.9961 \quad (R - \text{correlation coefficient})$$

### WASTE MANAGEMENT

Please refer to local legal requirements.

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### PROGRAM NA ANALIZATORY / APPLICATION / АДАПТАЦИЯ:

#### Parameters

No.	25	Prim.Wave.	340
Test	LDH	Sec.Wave.	450
Method	Kinetic	Sample Vol.	4
Direction	Descend	R1 Vol.	200
Unit	U/l	R2 Vol.	50
Decimals	1	Line. Limit	20
Incubation	8	Antigen Check	
Reaction	3   15	Substrat	0
<b>R1 Blank</b>		<b>Mix. R Blank</b>	
Lower	0	Lower	0
Upper	0	Upper	0
<b>Response</b>		<b>Linearity</b>	
Lower	-2.5	Lower	38
Upper	2.5	Upper	2000
Sample Vol.	45	Full Name	LDH
Dilution	5	Print No.	25

#### Calibration

Rule	Two Point Linear / One Point Linear
K Factor	0
Replicates	3
Interval	56
Sensitivity	0
Correlation	0
Difference	2.5
Blank Response	0   2.5
Coefficient	0
Difference	
Non-linear SD	0