## ACCENT-300 ASO

Nr kat. **7-340** (PL)

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Cat. No **7-340** (EN)

### INTENDED USE

Diagnostic kit for determination of anti-streptolysin O, used in automatic analyser ACCENT-300.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### METHOD PRINCIPLE

When an antigen-antibody reaction occurs between ASO in a sample and SLO which has been sensitized to latex particles, agglutination results. This agglutination is detected as an absorbance change (572 nm), with the magnitude of the change being proportional to the quantity of ASO in the sample. The actual concentration is then determined by interpolation from a calibration curve prepared from calibrators of known concentration.

### REAGENTS

**Package**  
 1-Reagent 1 x 40 ml  
 2-Reagent 1x 28.5 ml

The reagents, stored at 2-10°C are stable up to expiry date printed on the package. The reagents stored on board of the analyser at 2-10°C are stable for 8 weeks.

### Concentrations in the test

suspension of latex particles sensitized with SLO (pH 8.2) 0.17 w/v%  
 glycine buffer solution (pH 8.3)  
 preservative

### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Reagents with different lot numbers should not be interchanged or mixed.
- Do not interchange the caps of reagent bottles.
- Reagents should be mixed before use by gentle inverting the bottle several times
- Immediately after use, recap the bottles and store at 2-10°C.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

### SPECIMEN

Serum or plasma collected on heparine (sodium or lithium salt), EDTA (sodium or potassium salt) or citric acid. After the blood is completely clotted, the sample should be centrifuged and separated from the cells and fibrinogen. If the test cannot be done immediately, the sample should be placed in a tightly sealable container and stored at -20°C. Avoid repeated freezing and thawing. Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-Reagent and 2-Reagent are ready to use.  
 0.9% NaCl is recommended as a reagent blank

### REFERENCE VALUES <sup>3</sup>

serum, plasma	< 160 IU/ml
---------------	-------------

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY IMMUNO-CONTROL I (Cat. No 4-288) with each batch of samples.

For the calibration of automatic analysers systems the CORMAY ASO CALIBRATOR kit (Cat. No 4-278) is recommended.

The calibration curve should be prepared every 8 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

These metrological characteristics have been obtained using automatic analysers Accent-300 and Hitachi 917. Results may vary if a different instrument is used.

- Analytical range:** 20 – 1200 IU/ml.

If ASO level exceeds analytical range, dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by the dilution factor.

### Specificity / Interferences

Haemoglobin up to 0.5 g/dl, bilirubin up to 20 mg/dl and triglycerides up to 500 mg/dl do not interfere with the test.

### Precision

Repeatability (run to run) n = 20	Mean [IU/ml]	SD [IU/ml]	CV [%]
level 1	46.8	1.07	2.29
level 2	80.2	1.31	1.63
level 3	221.7	2.62	1.18

Reproducibility (day to day) n = 12	Mean [IU/ml]	SD [IU/ml]	CV [%]
level 1	48.8	2.83	5.81
level 2	78.8	2.60	3.30
level 3	219.8	5.24	2.38

### Method comparison

A comparison between CORMAY reagent (y) and another commercially available assay (x) using 50 samples gave following results:

$$y = 1.11 x - 44 \text{ IU/ml};$$

$$R = 0.945 \quad (R - \text{correlation coefficient})$$

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