ACCENT-300 URINE PROTEINS

Nr kat. 7-342 (PL)

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(EN)

INTENDED USE

Diagnostic kit for determination of total proteins in urine and cerebrospinal fluid, intended to use in automatic analyser ACCENT-300.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

In healthy people with properly functioning kidneys, proteins are actively reabsorbed in the proximal tubules and only small amounts of proteins (several mg per day) are excreted in urine. The measurement of total proteins concentration in urine is used in the diagnosis and treatment of heart and thyroid diseases, which are characterized by proteinuria or albuminuria.

The measurement of total proteins concentration in cerebrospinal fluid (CSF) is especially useful in detecting increased permeability of the blood-brain barrier and in detecting increased intrathecal synthesis of immunoglobulins. Increased concentration of protein in CSF may indicate brain tumors, intracerebral hemorrhage, brain injury, bacterial and viral encephalitis and multiple sclerosis.

METHOD PRINCIPLE

Direct, colorimetric method with pyrogallol red. At an acidic pH the protein amino acid groups, with the pyrogallol red-molybdate complex, form a coloured compound which shows a maximum absorbance at a wavelength of 600 nm. Colour intensity is proportional to the concentration of proteins in the sample.

REAGENTS

Package

1-Reagent 2 x 27.1 ml

The reagent when stored at 15-25°C is stable up to expiry date printed on the package. The reagents stored on board of the analyser at 2-10°C are stable for 10 weeks.

Concentrations in the reagent

Succinate buffer ≤ 60 mmol/l
Pyrogallol red ≤ 0.07 mmol/l
Sodium molybdate ≤ 0.05 mmol/l
Stabilizers, preservative

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Do not freeze the reagent.
- Reagent should be mixed before use by gentle inverting the bottle several times.
- The appearance of turbidity or control urine values outside the manufacturer's acceptable range may indicate of reagent instability.

SPECIMEN

Urine: Urine used for analysis may come from the first morning sample, random sample or timed collection sample according to the classification by ECLM².

In order to collect and prepare the samples only dedicated tubes and containers should be used.

Do not use preservatives.

Freshly collected urine should be kept at room temperature for about an hour and then cooled to 4°C. Direct reduction of the temperature on freshly collected urine can cause precipitation of minerals.

Samples with visible turbidity should be centrifuged before analysis. Determination of uncentrifuged samples may give increased results.

Urine samples are stable for 2 days at 2-8°C⁶.

Nevertheless it is recommended to perform the assay with freshly collected samples!

Cerebrospinal fluid: Cerebrospinal fluid should be centrifuged before analysis. The presence of blood in samples of cerebrospinal fluid may result in false results of protein determination. For proper interpretation of the results, cerebrospinal fluid must be determined simultaneously with a sample of blood taken from a patient at the same time.

The stability of the cerebrospinal fluid is 3 days at 2-8°C, 6 months at the temperature -20°C⁶.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent is ready to use.

0.9% NaCl is recommended as a reagent blank.

Actions required:

When performing assays in analyser ACCENT-300, there is a probability of **cross-contamination** affecting the tests results: URINE PROTEINS – CREATININE, URINE PROTEINS – UA, URINE PROTEINS – UA PLUS, CREA ENZYMATYC – URINE PROTEINS, FERRUM – URINE PROTEINS, GGT – URINE PROTEINS, UA – URINE PROTEINS, UA PLUS – URINE PROTEINS. To avoid this effect follow the recommendations contained in the instruction 51_03_24_007_ACCENT-300_CARRYOVER.

Calculation

For the calculation of proteins excreted over 24 hours, multiply the concentration (mg/dl) by the volume (dl) of the 24 hours urine.

REFERENCE VALUES^{7,10}

urine (adults)	< 15 mg/dl (0.15 g/l)	
urine 24-h (adults)	< 100 mg (0.10 g)	
cerebrospinal fluid	mg/dl	g/l
children 0 - 4 weeks	20 – 80	0.20 – 0.80
children >4 weeks, adults	15 – 45	0.15 – 0.45

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY URINE CONTROL LEVEL 1 (Cat. No 5-161) and LEVEL 2 (Cat. No 5-162).

For the calibration of automatic analysers the CORMAY URINE PROTEINS CALIBRATORS (Cat. No 5-181) is recommended. 0.9% NaCl should be used as a calibrator 0.

The calibration curve should be prepared every 10 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analyser ACCENT-300. Results may vary if a different instrument or a manual procedure is used.

- Sensitivity:** 1.8 mg/dl (0.018 g/l).

- Linearity:** up to 153 mg/dl (1.53 g/l).

For higher concentration dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

Specificity / Interferences

Hemoglobin up to 0.004 g/dl, ascorbic acid up to 20 mg/dl, creatinine up to 6 g/l, bilirubin up to 5 mg/dl, conjugated bilirubin up to 60 mg/dl, uric acid up to 85 mg/dl, glucose up to 35 g/l, citrates up to 250 mg/dl, oxalates up to 90 mg/dl, calcium ions up to 130 mg/dl, magnesium ions up to 1.8 g/l, phosphate ions up to 1.2 g/l, urea up to 50 g/l do not affect the results of the determination.

High concentration of iron (II) ions in the test sample may cause interferences¹².

Acetaminophen and some antibiotics from penicillins and aminoglycosides can interfere^{4, 11-13}.

Precision

Repeatability (run to run) n=10	Mean	SD	CV
	[mg/dl]	[mg/dl]	[%]
level 1	20.54	0.15	0.73
level 2	55.54	0.86	1.54
Reproducibility (day to day) n=20	Mean	SD	CV
	[mg/dl]	[mg/dl]	[%]
level 1	20.50	0.47	2.30
level 2	55.90	1.14	2.05

Method comparison

A comparison between total proteins values determined at ACCENT-300 (y) and at ADVIA 1650 (x) using 75 urine samples gave following results:

$$y = 0.9737x + 16.068 \text{ mg/dl};$$

$$R = 0.994 \quad (R - \text{correlation coefficient})$$

A comparison between total proteins values determined at ACCENT-300 (y) and at ADVIA 1650 (x) using 30 CSF samples gave following results:

$$y = 1.0737x + 1.5928 \text{ mg/dl};$$

$$R = 0.999 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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PROGRAM NA ANALIZATOR / APPLICATION / АДАПТАЦИЯ :

Parameters

No.	56	Prim.Wave.	630
Test	PROT U	Sec.Wave.	0
Method	Endpoint	Sample Vol.	4
Direction	Ascend	R1 Vol.	180
Unit	mg/dl	R2 Vol.	0
Decimals	1	Line. Limit	20
Incubation	0	Antigen Check	
Reaction	0 8	Substrat	0
R1 Blank		Mix. R Blank	
Lower	0	Lower	0
Upper	0	Upper	0
Response		Linearity	
Lower	-2.5	Lower	1.8
Upper	2.5	Upper	153
Sample Vol.	45	Full Name	PROT U
Dilution	5	Print No.	56

Calibration

Method	Two – Point Linear
K Factor	0
Replicates	3
Interval	70
Sensitivity	0
Correlation	0
Difference	2.5
Blank Response	0 2.5
Coefficient	0
Difference	
Non-linear SD	0

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