ACCENT-300 PHOSPHORUS

Nr kat. 7-343 (PL)

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(EN)

INTENDED USE

Diagnostic kit for determination of inorganic phosphorus concentration, used in automatic analyser ACCENT-300.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Phosphorus is present in all body cells as a component of nucleic acids, phospholipids and phosphoproteins. Phosphorus is essential for intracellular storage and conversion of energy (ATP, creatine phosphate) and participates in carbohydrates metabolism. In the blood phosphorus is present as a mixture of inorganic phosphates HPO_4^{2-} and H_2PO_4^- . Besides phosphorus and calcium constitute mineral portion of bone. Continuous flux of phosphorus in organism is controlled by parathyroid hormone (PTH), vitamin D and calcitonin. Phosphorus serum level abnormalities are caused usually by disorders of vitamin D metabolism or parathyroid and kidney diseases.

METHOD PRINCIPLE

Direct phosphomolybdate reaction without deproteinization. Phosphate ions form with molybdate ions in acid solution proportional amounts of unreduced phosphomolybdate complex. The concentration of the complex formed is determined by measuring its absorbance.

REAGENTS

Package

1-Reagent 2 x 40 ml

The reagent, stored at 2-8°C is stable up to expiry date printed on the package. The reagent stored on board of the analyser at 2-10°C is stable for 12 weeks.

Concentrations in the test

ammonium molybdate	0.4 mmol/l
sulphuric acid	150 mmol/l
hydrochloric acid	100 mmol/l
detergents	

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Contaminated glassware is the greatest source of error. Disposable plastic ware is recommended for the test.
- 1-REAGENT meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Ingredients:

1- REAGENT contains sulfuric acid (VI) and hydrochloric acid.

Danger.

H314 Causes severe skin burns and eye damage.
P280 Wear protective gloves/protective clothing/eye protection/face protection.

P301+P330+P331 IF SWALLOWED: rinse mouth. Do NOT induce vomiting.

P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing.

P305 +P351 +P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P310 - Immediately call a POISON CENTER or doctor.

SPECIMEN

Serum, heparinized plasma (recommended: heparine lithium, sodium or ammonium salt) free from hemolysis, 24-hours urine.

Serum is the preferred specimen. Level of inorganic phosphate in heparinized plasma is about 0.2 to 0.3 mg/dl (0.06-0.10 mmol/l) lower than in serum.

Serum should be separated from red blood cells as soon as possible after blood collection, because erythrocytes contain several times higher phosphate concentration than normal serum.

Urine preparation: To prevent phosphate precipitation in urine, specimens should be collected in HCl, 20-30 ml of 6 mol/L for 24-h specimen. Then dilute 1 part of acidified urine with 10 parts of distilled water. Multiply the result by the dilution factor.

Serum and plasma can be stored up to 7 days at 2 – 8°C. For longer storage samples should be frozen at -20°C.

24-hours urine samples can be stored up to 7 days at 2-8°C. Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent is ready to use.

Deionized water is recommended as a reagent blank.

REFERENCE VALUES ⁷

serum / plasma	mg/dl	mmol/l
age: 0 – 10 d	4.5 – 9.0	1.45 – 2.91
10 d – 24 mo	4.5 – 6.7	1.45 – 2.16
24 mo – 12 y	4.5 – 5.5	1.45 – 1.78
12 – 60 y	2.7 – 4.5	0.87 – 1.45
> 60 y males	2.3 – 3.7	0.74 – 1.20
> 60 y females	2.8 – 4.1	0.90 – 1.32
24-hours urine	g/24h	mmol/24h
	0.4 – 1.3	12.9 – 42.0

It is recommended for each laboratory to establish its own reference ranges for local population.

Phosphorus concentration in 24-hours urine – calculation

phosphorus concentration in 24-hours urine [g/24h]	=	phosphorus concentration in sample of 24-hours urine [mg/dl]	x	urine volume of 24-hours urine [dl/24h]	÷	1000
--	---	--	---	---	---	------

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) for determination in serum or CORMAY URINE CONTROL LEVEL 1 (Cat. No 5-161) and LEVEL 2 (Cat. No 5-162) for determination in urine.

For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended. Calibrator and deionised water should be used for calibration.

The calibration curve should be prepared every 12 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analyser ACCENT-300. Results may vary if a different instrument or a manual procedure is used.

- Sensitivity:** 0.23 mg/dl (0.074 mmol/l).
- Linearity:** up to 18.6 mg/dl (6.01 mmol/l).

For higher concentration dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

Specificity / Interferences

Haemoglobin up to 0.16 g/dl, ascorbate up to 62 mg/l, bilirubin up to 15 mg/dl and triglycerides up to 1000 mg/dl do not interfere with the test.

Precision

Repeatability (run to run) n = 10	Mean	SD	CV
	[mg/dl]	[mg/dl]	[%]
level 1	3.86	0.02	0.63
level 2	6.76	0.04	0.52
Reproducibility (day to day) n = 20	Mean	SD	CV
	[mg/dl]	[mg/dl]	[%]
level 1	4.39	0.04	0.89
level 2	7.20	0.09	1.24

Method comparison

A comparison between phosphorus values determined at ACCENT-300 (y) and at ADVIA 1650 (x) using 23 samples gave following results:

$y = 1.0189x - 0.0609$ mg/dl;

$R = 0.974$ (R – correlation coefficient)

WASTE MANAGEMENT

Please refer to local legal requirements.

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PROGRAM NA ANALIZATOR / APPLICATION / АДАПТАЦИЯ:

Parameters

No.	15	Prim.Wave.	340
Test	PHOS	Sec.Wave.	700
Method	Endpoint	Sample Vol.	3
Direction	Ascend	R1 Vol.	300
Unit	mg/dl	R2 Vol.	0
Decimals	2	Line. Limit	
Incubation	0	Antigen Check	
Reaction	0 35	Substrat	0
R1 Blank		Mix. R Blank	
Lower	0	Lower	0
Upper	0	Upper	0
Response		Linearity	
Lower	-2.5	Lower	0.23
Upper	2.5	Upper	18.6
Sample Vol.	45	Full Name	Phosphorus
Dilution	5	Print No.	15

Calibration

Rule	Multi- Point Linear
K Factor	0
Replicates	3
Interval	84
Sensitivity	0
Correlation	0
Difference	2.5
Blank Response	0 2.5
Coefficient	0
Difference	
Non-linear SD	0

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