## ACCENT-300 ACP

Nr kat. 7-349

(PL)

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(EN)

### INTENDED USE

Diagnostic kit for determination of total, non prostatic and prostatic acid phosphatase activity, used in automatic analyser ACCENT-300.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Acid phosphatases (ACP) are a group of enzymes that can be found in liver, spleen, erythrocytes, platelets, bone marrow and prostate gland. The ratio between total acid phosphatase and prostatic phosphatase is healthy males is about 1:1. Activity of total acid phosphatase increases in case of Paget's disease, hyperparathyroidism, bone cancer, Gaucher'a disease, Niemann-Pick disease and myelocytic leukaemia. Prostatic acid phosphatase levels increase in case of prostatic cancer.

### METHOD PRINCIPLE

Acid phosphatases catalyses the hydrolysis, in acid medium, of 1-naphthylphosphate (1-NP) into 1-naphtol and phosphate. 1-naphtol reacts with diazo-2-chloro-5-toluene (Fast Red TR salt), forming an azo dye compound with absorbance increase is proportional to the total acid phosphatase activity (ACP). Prostatic acid phosphatase (ACP-P) is inhibited by tartrate and is detected with an indirect method by subtraction between ACP and non-prostatic acid phosphatase (ACP-NP).

### REAGENTS

#### Package

REAGENT 1a	6 x 20 ml
REAGENT 1b	6 vials (powder)
TARTRATE REAGENT	1 vial (powder)
ACETIC ACID SOLUTION (acetic acid 0.1 mol/l)	1 x 3 ml
small spoon	1 pcs.

The reagents, stored at 2-8°C are stable up to expiry date printed on the package.

#### Working reagents stability

Stability of working reagents ACP and ACP-NP after reconstitution: 7 days at 2-8°C, if contamination avoided and bottle recapped immediately after use.

#### Concentrations in the working reagent ACP

citrate buffer (pH 5.4)	0.1 mol/l
1-NP	10 mmol/l
Fast Red TR salt	> 1.0 mmol/l

#### Concentrations in the working reagent ACP-NP

citrate buffer (pH 5.4)	0.1 mol/l
1-NP	10 mmol/l
Fast Red TR salt	> 1.0 mmol/l
sodium tartrate	0.18 mol/l

### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- It is recommended to carry out the prostatic acid phosphatase test only on samples with total acid phosphatase activity > 3.5 U/l as it is possible to obtain not reliable with a total acid phosphatase activity lower than the above-mentioned limit.
- Controls and calibrators used in determination of ACP should be used according to the guidelines contained in their instruction for use.
- The reagents must be limpid; do not use if turbid.
- A slight pink coloration of working reagents ACP and ACP-NP does not affect the reagents performance.
- ACETIC ACID SOLUTION meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

#### Ingredients:

ACETIC ACID SOLUTION contains acetic acid.

#### Danger

H314 Causes severe skin burns and eye damage.

H318 Causes eye damage.

P280: Wear protective gloves/protective clothing/eye protection/face protection.

P302+P352: IF ON SKIN: Wash with plenty

of water.

P305+P351+P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P321: Specific treatment is urgently needed (go to see a doctor with the Safety data sheet for this product).

P501: Dispose of the contents/containers in accordance with the current legislation on waste treatment.

### SPECIMEN

Fresh serum.

Do not use haemolized, lipemic and icteric samples. Since the sample activity decreases of about 50% within an hour at room temperature, acidify the specimen immediately after collection by adding 20µl ACETIC ACID SOLUTION for each ml of fresh serum.

Stability of acidified sample: 5 days at 2-8°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

Deionised water is recommended as a reagent blank.

#### Actions required:

When performing assays in analyser ACCENT-300, there is a probability of **cross-contamination** affecting the tests results: CALCIUM - ACP. To avoid this effect follow the recommendations contained in the instruction 51\_03\_24\_007\_ACCENT-300\_CARRYOVER.

### ACP ASSAY

**Working reagent ACP:** Dissolve the contents of one bottle of Reagent 1b with the contents of one bottle of Reagent 1a. Let stand for about 10 minutes and gently mix. The contents pour into bottle of Reagent 1a, put in analyser and use application for ACP.

### ACP-NP ASSAY

**Working reagent ACP:** Dissolve the contents of one bottle of Reagent 1b with the contents of one bottle of Reagent 1a. Let stand for about 10 minutes and gently mix. In order to determine ACP-NP activity add one small spoon contents (corresponding to 35-40 mg) of Tartrate Reagent for each 5 ml working reagent. Gently mix until complete solution. Such solution put in analyser and use application for ACP-NP.

### CALCULATION OF ACP-P (prostatic acid phosphatase)

ACP-P [U/l] = ACP – ACP-NP

### REFERENCE VALUES <sup>4</sup>

ACP		
serum	U/l	µkat/l
male	2.5 – 11.7	0.042 – 0.195
female	0.3 – 9.2	0.005 – 0.154
ACP-P		
serum (tartrate inhibition)	U/l	µkat/l
male	0.2 – 3.5	0.003 – 0.058
female	0.0 – 0.8	0.000 – 0.013

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control of ACP it is recommended to use, with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173). For internal quality control of ACP-NP it is recommended to use, with each batch of samples, the control sera.

For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) or LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

The calibration curve should be prepared with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analyser Hitachi 704. Results may vary if a different instrument is used.

#### Sensitivity:

ACP: 0.35 U/l (0.006 µkat/l)

ACP-NP: 0.39 U/l (0.007 µkat/l).

- Linearity ACP/ACP-NP:** up to 75 U/l (1.25 µkat/l).

For higher activity dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

#### Specificity / Interferences

Haemoglobin up to 0.15 g/dl, triglycerides up to 400 mg/dl, bilirubin up to 20 mg/dl and do not interfere with the test

#### Precision

ACP			
Repeatability (run to run) n = 20	Mean [U/l]	SD [U/l]	CV [%]
level 1	6.7	0.26	3.87
level 2	13.8	0.37	2.70
level 3	28.0	0.64	2.30

Reproducibility (day to day) n = 20			
	Mean [U/l]	SD [U/l]	CV [%]
level 1	6.64	0.27	4.06
level 2	13.84	0.37	2.64
level 3	28.10	0.62	2.22
ACP-NP			
Repeatability (run to run) n = 20			
	Mean [U/l]	SD [U/l]	CV [%]
level 1	2.8	0.19	6.73
level 2	3.9	0.18	4.50
level 3	6.0	0.23	3.86
Reproducibility (day to day) n = 20			
	Mean [U/l]	SD [U/l]	CV [%]
level 1	2.93	0.17	5.88
level 2	3.90	0.19	4.95
level 3	5.96	0.23	3.82

#### Method comparison

A comparison between CORMAY reagent (y) and another commercially available assay (x) using 60 samples gave following results:

#### ACP

y = 1.0162 x + 0.10151 U/l;

R = 0.997 (R – correlation coefficient)

#### ACP-NP

y = 1.1325 x - 0.09887 U/l;

R = 0.998 (R – correlation coefficient)

### WASTE MANAGEMENT

Please refer to local legal requirements.

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