## ACCENT-300 GLUCOSE HEX

Nr kat. 7-352 (PL)

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### LITERATURA

- Barham P., Trinder P.: Analyst 97, 142-145 (1972).
- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics 4th ed., PA: WB Saunders, 868-869, 2006.
- McPherson Richard A., Pincus Matthew R.: Henry's Clinical Diagnosis and Management by Laboratory Methods, 22nd ed. Philadelphia, PA: WB Saunders, 2011.
- Dujmovic and F. Deisenhammer, Stability of cerebrospinal fluid/serum glucose ratio and cerebrospinal fluid lactate concentrations over 24 h: analysis of repeated measurements, Clinical Chemistry and Laboratory Medicine, Volume 48, Issue 2, 2010, pp. 209-212.
- Zalecenia kliniczne dotyczące postępowania u chorych na cukrzycę 2014, Diabetologia Kliniczna, tom 3, suplement A, 2014.
- Sacks, David B., et al. Guidelines and recommendations for laboratory analysis in the diagnosis and management of diabetes mellitus. Clinical chemistry 48.3 (2002): 436-472.
- Miles RR, Roberts RF, Putnam AR, Roberts WL. Comparison of serum and heparinized plasma samples for measurement of chemistry analytes. Clin Chem 2004;50:1704-1706.
- Alan H.B. Wu: Tietz Clinical Guide to Laboratory Tests, 4th ed. WB Saunders, 444-450 (2006).
- Young DS., Effects of drugs on clinical laboratory tests, 5. Vol. 2. Washington DC, USA: AACC Press (2000).

Data wydania: 10.2020

## ACCENT-300 GLUCOSE HEX

Cat. No **7-352** (EN)

### INTENDED USE

Diagnostic kit for determination of glucose concentration, used in automatic analyser ACCENT-300.


The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Glucose is a simple six-carbon sugar. Oxidative metabolism of glucose provides the energy for most cellular processes. Glucose level in the blood is tightly controlled by several hormones. Elevated glucose level is the classic sign of diabetes mellitus. Glucose level abnormalities (hyper- or hypoglycemia) might be caused also by pancreas tumors and diseases of liver, thyroid gland or adrenal glands.

### METHOD PRINCIPLE

Enzymatic method with hexokinase and glucose-6-phosphate dehydrogenase (G6P-DH).

The rate of NADH formation is directly proportional to the glucose concentration in the sample.

### REAGENTS

#### Package

1-Reagent 4 x 31.5 ml  
 2-Reagent 4 x 7.1 ml

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. The reagents are stable for 12 weeks on board the analyser at 2-10°C

### Concentrations in the tests

#### 1-Reagent

PIPES buffer (pH 7.5) 80 mmol/l  
 Mg<sup>2+</sup> 10 mmol/l  
 ATP 4 mmol/l  
 NAD 3 mmol/l

#### preservative

#### 2-Reagent

hexokinase ≥ 4500 U/l  
 glucose-6-phosphate dehydrogenase (G6P-DH) ≥ 14000 U/l  
 preservative

### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Do not use after expiry date.
- Do not freeze reagents.

- Do not interchange caps.
- Reagents should be mixed before use by gentle inverting the bottles several times.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

### SPECIMEN

EDTA or heparinized plasma / serum, free from hemolysis, cerebrospinal fluid, urine.

**Plasma / Serum.** Serum and plasma specimens should be separated from cells within 30 minutes after collection. Plasma specimen which is not assayed immediately after collection should be kept in tubes containing sodium fluoride or sodium iodoacetate. These compounds adding prevent glycolysis and stabilize glucose level.

Serum and plasma can be stored up to 2 days at 4°C.<sup>3</sup>

Plasma is the specimen recommended for the glucose determination in the blood.<sup>5</sup>

**Cerebrospinal fluid.** Glucose concentration in cerebrospinal fluid should be measured directly after specimen collection. Cerebrospinal fluid must be analysed simultaneously with a blood sample.

After centrifuge CSF sample can be stored up to 24 hours at 4°C.<sup>4</sup>

**Urine.** Collect 24-hour sample in dark bottle and keep on ice. Preserve sample by adding 5 ml of glacial acetic acid to the container before starting the collection. The final pH of the sample should be between 4 and 5. Centrifuge samples with visible turbidity or precipitates before analysis.

Urine can be stored up to 24 hour at 4°C. Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-Reagent and 2-Reagent are ready to use.

For reagent blank deionized water is recommended.

#### Actions required:

When performing assays in analyser ACCENT-300, there is a probability of **cross-contamination** affecting the tests results: GLUCOSE HEX - MG. To avoid this effect follow the recommendations contained in the instruction 51\_03\_24\_007\_ACCENT-300\_CARRYOVER.

### REFERENCE VALUES

	mg/dl	mmol/l
plasma, serum <sup>5,6,7</sup>	70 – 99	3.9 – 5.5
urine (24h) <sup>8</sup>	1 – 15	0.1 – 0.8
cerebrospinal fluid <sup>8</sup>	40 – 70	2.2 – 3.9

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use with each batch of samples, the following controls: CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) for determination in serum; CORMAY

URINE CONTROL LEVEL 1 (Cat. No 5-161) and LEVEL 2 (Cat. No 5-162) for determination in urine.

For the calibration of automatic analyser the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

The calibration curve should be prepared every 12 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analyser ACCENT-300. Results may vary if a different instrument or a manual procedure is used.

- Sensitivity:** 4.5 mg/dl (0.25 mmol/l).

- Linearity:** up to 1000 mg/dl (55.5 mmol/l).

For higher concentration dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

#### Specificity / Interferences

Haemoglobin up to 1.25 g/dl, bilirubin up to 40 mg/dl, ascorbate up to 62 mg/L and triglycerides up to 1000 mg/dl do not interfere with the test. Some medicines can interfere.<sup>9</sup>

#### Precision

Repeatability (run to run) n = 10	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	87.92	2.01	2.29
level 2	288.83	1.92	0.67
Reproducibility (day to day) n = 20	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	89.24	1.00	1.13
level 2	283.48	3.33	1.18

#### Method comparison

A comparison between glucose values determined at **ACCENT-300** (y) and at **COBAS INTEGRA 400** (x) using 100 samples gave following results:

$$y = 0.999x + 1.1334 \text{ mg/dl;}$$

$$R = 0.996 \quad (R - \text{correlation coefficient})$$

### WASTE MANAGEMENT

Please refer to local legal requirements.

### LITERATURE

- Barham P., Trinder P.: Analyst 97, 142-145 (1972).
- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics 4th ed., PA: WB Saunders, 868-869, 2006.
- McPherson Richard A., Pincus Matthew R.: Henry's Clinical Diagnosis and Management by Laboratory Methods, 22nd ed. Philadelphia, PA: WB Saunders, 2011.
- Dujmovic and F. Deisenhammer, Stability of cerebrospinal fluid/serum glucose ratio and cerebrospinal fluid lactate concentrations over 24 h: analysis of repeated measurements, Clinical Chemistry and Laboratory Medicine, Volume 48, Issue 2, 2010, pp. 209-212.
- Zalecenia kliniczne dotyczące postępowania u chorych na cukrzycę 2014, Diabetologia Kliniczna, tom 3, suplement A, 2014.
- Sacks, David B., et al. Guidelines and recommendations for laboratory analysis in the diagnosis and management of diabetes mellitus. Clinical chemistry 48.3 (2002): 436-472.
- Miles RR, Roberts RF, Putnam AR, Roberts WL. Comparison of serum and heparinized plasma samples for measurement of chemistry analytes. Clin Chem 2004;50:1704-1706.
- Alan H.B. Wu: Tietz Clinical Guide to Laboratory Tests, 4th ed. WB Saunders, 444-450 (2006).
- Young DS., Effects of drugs on clinical laboratory tests, 5. Vol. 2. Washington DC, USA: AACC Press (2000).

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### PROGRAM NA ANALIZATOR / APPLICATION / АДАПТАЦИЯ:

#### Parameters

No.	1	Prim.Wave.	340
Test	GLUC HEX	Sec.Wave.	
Method	Endpoint	Sample Vol.	3
Direction	Ascend	R1 Vol.	200
Unit	mg/dl	R2 Vol.	40
Decimals	1	Line. Limit	
Incubation	10	Antigen Check	
Reaction	-2   25	Substrat	
<b>R1 Blank</b>		<b>Mix. R Blank</b>	
Lower		Lower	
Upper		Upper	
<b>Response</b>		<b>Linearity</b>	
Lower	-2.5	Lower	4.5
Upper	2.5	Upper	1000
Sample Vol.	45	Full Name	GLUCOSE HEX
Dilution	5	Print No.	1

#### Calibration

Method	Two Point Linear
K-Factor	
Replicates	3
Interval	84
Sensitivity	0
Correlation	0
Difference	2.5
Blank Response	0   2.5
Coefficient	0
Difference	
Non-linear SD	0

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