## ACCENT-300 AMYLASE

Nr kat. 7-355 (PL)

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Cat.No 7-355 (EN)

### INTENDED USE

Diagnostic kit for determination of  $\alpha$ -amylase activity intended to use in automatic analyzer ACCENT-300.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

$\alpha$ -Amylase is a digestive enzyme secreted by salivary glands and pancreas. Low level of amylase is also found in skeletal muscle, adipose tissue and fallopian tubes.  $\alpha$ -Amylase is measured generally in pancreas diseases. Elevation of amylase activity is observed also due to inflammation of abdominal cavity or salivary glands.

### METHOD PRINCIPLE

2-Chloro-4-nitrophenyl- $\alpha$ -maltotrioidide (CNP-G3) is a direct substrate for determination of  $\alpha$ -amylase activity, which does not require the presence of ancillary enzymes.

10 CNP-G3  $\alpha$ -amylase  $\rightarrow$  9 CNP + CNP-G2 + 9 maltotriose + glucose

The rate of 2-chloro-4-nitrophenol formation can be monitored at 405 nm and is proportional to the  $\alpha$ -amylase activity.

### REAGENTS

**Package**  
 1-Reagent 2 x 40 ml

The reagent, stored at 2-8°C is stable up to expiry date printed on the package. The reagent is stable for 11 weeks on board the analyzer at 2-10°C.

### Concentrations in the test

MES buffer < 120 mmol/l  
 calcium acetate < 7 mmol/l  
 potassium hydroxide < 40 mmol/l  
 potassium thiocyanate < 1100 mmol/l  
 2-chloro-4-nitrophenyl- $\alpha$ -maltotrioidide < 2 mmol/l  
 preservative, stabilizer

### Warnings and notes

- Protect from direct sunlight.
- Prevent the reagent from microbiological contamination and from saliva and sweat  $\alpha$ -amylase! Saliva and sweat contain  $\alpha$ -amylase. Do not pipette by the mouth, avoid skin contact with reagent, specimens, tips, cuvettes. Ensure to use automatic pipettes and laboratory gloves.

### SPECIMEN

Serum or plasma collected on heparin, free from hemolysis, urine.

Do not use anticoagulants: EDTA, citrates and oxalates as they inhibit amylase activity.

Serum / plasma can be stored for 7 days at 15-25°C or for one month at 2-8°C.<sup>7</sup>

Urine can be stored for 2 days at 15-25°C or for 10 days at 2-8°C.<sup>9</sup> Amylase is very unstable in acid urine. Adjust pH to approximately 7.0 before storage.

Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-Reagent is ready to use.  
 Deionized water is recommended as a reagent blank.

### REFERENCE VALUES <sup>8</sup>

serum / plasma	U/l	$\mu$ kat/l
	20 – 104	0.34 – 1.77
urine	U/l	$\mu$ kat/l
	32 – 641	0.54 – 10.90

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use the following controls with each batch of samples:

CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) for determination in serum; CORMAY URINE CONTROL LEVEL 1 (Cat. No 5-161) and LEVEL 2 (Cat. No 5-162) for determination in urine. For the calibration of automatic analyzers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and/or LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

The calibration curve should be prepared every 4 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analyzer ACCENT-300. Results may vary if a different instrument or a manual procedure is used.

▪ **Sensitivity:** 9 U/l (0.15  $\mu$ kat/l).

▪ **Linearity:** up to 1500 U/l (25  $\mu$ kat/l).

▪ **Specificity / Interferences**

Haemoglobin up to 2.5 g/dl, ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl and triglycerides up to 1000 mg/dl do not interfere with the test.

▪ **Precision**

Repeatability (run to run) n = 20	Mean [U/l]	SD [U/l]	CV [%]
level 1	120.65	0.96	0.80
level 2	367.60	1.50	0.41
Reproducibility (day to day) n = 80	Mean [U/l]	SD [U/l]	CV [%]
level 1	71.65	2.50	3.50
level 2	315.93	9.68	3.06

### Method comparison

A comparison between amylase values determined at ACCENT-200 (y) and at ACCENT-300 (x) using 52 samples gave following results:

$$y = 1.0241x - 0.3934 \text{ U/l};$$

$$R = 0.9999 \quad (R - \text{correlation coefficient})$$

### WASTE MANAGEMENT

Please refer to local legal requirements.

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