## ACCENT-300 CHOLINESTERASE

Nr kat. **7-356** (PL)

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Cat. No **7-356** (EN)

### INTENDED USE

Diagnostic kit for determination of cholinesterase activity, intended to use in automatic analyzer: ACCENT-300. The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

There are two cholinesterases (CHE and ACHE) differing in substrate specificity, tissue of origin and biological role. The term cholinesterase (ACHE), also known as acetylcholine acetylhydrolase, is found in erythrocytes, in the lungs, spleen and in grey matter of the brain. The pseudocholinesterase (CHE), also referred to as acylcholine acylhydrolase, is found in serum, the liver, pancreas, heart and in the white matter of brain. The assay of serum cholinesterase (CHE) is useful to diagnose: liver disorders, hepatitis, cirrhosis, carcinoma with metastasis, sensitivity to succinylcholine administration and pesticide poisoning. Levels decrease in all of the diseases above.

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The method uses butyrylthiocholine as the specific substrate for cholinesterase (CHE). Cholinesterase catalyses the hydrolysis of butyrylthiocholine substrate forming butyrate and thiocholine, in presence of potassium hexacyanoferrate (III). Thiocholine reduces potassium hexacyanoferrate (III) (yellow colour) to potassium hexacyanoferrate (II) (colourless). The decrease in absorbance is directly proportional to the CHE activity in the sample.

### REAGENTS

#### Package

1-Reagent 1 x 37 ml  
2-Reagent 1 x 10 ml

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. The reagents are stable for 12 weeks on board the analyser at 2-10°C.

### Concentrations in the test

#### 1-Reagent

pyrophosphate buffer pH 7.7 65 mmol/l  
hexacyanoferrate (III) 2 mmol/l

#### 2-Reagent

Good's Buffer pH 4.0 20 mmol/l  
butyrylthiocholine iodide 65 mmol/l

### Warnings and notes

- Protect from direct sunlight, evaporation and avoid contamination!

- Do not freeze reagents.
- Immediately after use, recap the bottles and store at 2-8°C.
- Do not use after expiry date.
- Do not interchange caps.
- Reagents should be mixed before use by gentle inverting the bottle several times.
- The appearance of turbidity or control sera values outside the manufacturer's acceptable range may indicate of reagent instability.

### SPECIMEN

Fresh serum free from haemolysis, plasma (EDTA, heparin) not hemolyzed.

Serum / plasma should be separated from red blood cells as soon as possible after blood collection.

Do not use citrate, borate, oxalate and fluoride as an anticoagulant because it inhibits cholinesterase activity.

It is recommended to follow CLSI procedures regarding specimen collecting and handling.

Sample may be stored for up to 15 days at 2-8°C or 12 months at -20°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-Reagent and 2-Reagent are ready to use.

For reagent blank deionized water is recommended.

### REFERENCE VALUES <sup>45</sup>

serum / plasma	37°C	
female	4000 – 12600 U/l	67 – 210 µkat/l
male	5100 – 11700 U/l	85 – 195 µkat/l

In infants up to 6 months of age, cholinesterase activity is 40% to 50% higher than in adults. In young adult (< 35 years) women, the enzyme activity is approximately 64% to 74% of that in adult males. The activity decreases during pregnancy.

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of sample the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173).

For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) or CORMAY MULTICALIBRATOR LEVEL 2 (Cat. No 5-175; 5-177) **is recommended depending on the calibrator lot number.**

The calibration curve should be prepared every 12 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

These metrological characteristics have been obtained using the automatic analyser ACCENT-300. Results may vary if a different instrument or a manual procedure is used.

- Sensitivity:** 199 U/l (3.32 µkat/l).
- Linearity:** up to 20000 U/l (333 µkat/l).

For higher activity dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

- Specificity / Interferences**

Haemoglobin up to 5 g/dl, ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl and triglycerides up to 1000 mg/dl do not interfere with the test.

- Precision**

Repeatability (run to run) n = 10	Mean	SD	CV
	[U/l]	[U/l]	[%]
level 1	5592	51.28	0.92
level 2	7634	81.66	1.07
Reproducibility (day to day) n = 20	Mean	SD	CV
	[U/l]	[U/l]	[%]
level 1	2272	30.69	1.35
level 2	7336	139.01	1.89

- Method comparison**

A comparison between CORMAY reagent (y) and another commercially available assay (x) using 46 samples gave following results:

$$y = 0.9728x + 277.09 \text{ U/l};$$

$$R = 0.9976 \quad (R - \text{correlation coefficient})$$

### WASTE MANAGEMENT

Please refer to local legal requirements.

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