ACCENT-300 FERRUM

Nr kat. 7-358 (PL)

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Cat. No **7-358** (EN)

INTENDED USE

Diagnostic kit for determination of iron concentration, used in automatic analyser ACCENT-300.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Iron is the most abundant trace element in the organism. Most of the iron in humans is located within heme molecule which is incorporated into hemoglobin, myoglobin, catalase, peroxidase and cytochromes. Iron is stored bound to ferritin or hemosiderin and is transported by transferrin. Measurement of iron level is valuable especially in diagnosis of different types of anemia.

METHOD PRINCIPLE

Colorimetric method with ferrozine, without deproteinization.

Iron ions (Fe^{3+}), bounded in blood to transferrin are released in acid solution and in the presence of detergents and reduced to Fe^{2+} by ascorbate. Fe^{2+} forms with 3-(2-pyridyl)-5,6-bis(2-[4-phenyl sulfonic acid])-1,2,4-triazine sodium salt (ferrozine) coloured complex, Cu^{2+} ions are bound by thiourea. The colour intensity is directly related to the iron concentration.

REAGENTS

Package

1-Reagent 3 x 39.5 ml
2-Reagent 2 x 12.5 ml

The reagents are stable up to expiry date printed on the package, if stored at 2-8°C. The reagents are stable for 11 weeks on board the analyser at 2-10°C.

Concentrations in the reagent

1-Reagent

citric acid (pH 1.9) ≤ 240 mmol/l
thiourea ≤ 108 mmol/l
detergent ≤ 7 %

2-Reagent

sodium ascorbate ≤ 150 mmol/l
3-(2-pyridyl)-5,6-bis(2-[5-furyl sulfonic acid])-1,2,4-triazine sodium salt (ferrozine) ≤ 6 mmol/l
preservative
stabilizer

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- 1-Reagent meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Ingredients:

1-Reagent contains polyethylene glycol monoalkyl ether.

Danger

H314 Causes severe skin burns and eye damage.

EUH208 Contains thiourea. May produce an allergic reaction (1 Reagent).

EUH208 Contains 1-[1,3-Bis(hydroxymethyl)-2,5-dioximidazolidin-4-yl]-1,3-bis(hydroxymethyl) urea. May produce an allergic reaction (2-Reagent).

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P310 Immediately call a POISON CENTER or doctor.

SPECIMEN ¹³

Serum free from hemolysis, collected in plastic tubes.

Serum should be separated from red blood cells as soon as possible after blood collection.

Serum can be stored up to 7 days at 15-25°C or up to 3 weeks at 2-8°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.

For reagent blank deionized water is recommended.

Actions required:

When performing assays in analyser ACCENT-300, there is a probability of **cross-contamination** affecting the tests results: FERRUM – UA, FERRUM – UA PLUS, FERRUM – UIBC, FERRUM – dTIBC, FERRUM – URINE PROTEINS, FERRUM – FERRITIN, dTIBC – FERRUM, UIBC - FERRUM. To avoid this effect follow the recommendations contained in the instruction 51_03_24_007_ACCENT-300_CARRYOVER.

REFERENCE VALUES ^{8, 10}

serum	µg/dl	µmol/l
newborns	100 – 250	17.9 – 44.8
infants	40 – 100	7.2 – 17.9
children	50 – 120	9.0 – 21.5
adult females	50 – 170	9.0 – 30.4
adult males	65 – 175	11.6 – 31.3

Samples should be taken in the morning from patients in a fasting state, since iron values decrease by 30% during the course of the day.

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use, |with each batch of samples the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173)

For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) or CORMAY MULTICALIBRATOR LEVEL 2 (Cat. No 5-175; 5-177) is recommended depending on the calibrator lot number. As a 0 calibrator deionised water should be used.

The calibration curve should be prepared every 11 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analyser ACCENT-300. Results may vary if a different instrument or a manual procedure is used.

- Sensitivity:** 7.3 µg/dl (1.31 µmol/l).

- Linearity:** up to 1000 µg/dl (179 µmol/l).

For higher concentrations dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

Specificity / Interferences

Ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl, triglycerides up to 1000 mg/dl and copper up to 500 µg/dl do not interfere with the test. Haemoglobin interferes even in small amount with the determination.

Precision

Repeatability (run to run) n = 10	Mean [µg/dl]	SD [µg/dl]	CV [%]
level 1	59.53	1.31	2.21
level 2	237.15	2.04	0.86
Reproducibility (day to day) n = 20	Mean [µg/dl]	SD [µg/dl]	CV [%]
level 1	242.75	7.25	2.98
level 2	66.26	1.82	2.75

Method comparison

A comparison between ferrum values determined at **Hitachi 912** (y) and at **ACCENT-300** (x) using 54 samples gave following results:

$$y = 0.9253x + 1.4445 \mu\text{g/dl};$$

$$R = 0.997 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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