ACCENT-300 UA PLUS

Nr kat. 7-363

(PL)

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Cat. No 7-363 (EN)

INTENDED USE

Diagnostic kit with ascorbate oxidase for determination of uric acid concentration, used in automatic analyser ACCENT-300.



The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Uric acid is a product of purine catabolism. It is produced in the liver and excreted in the urine. Both, the amount of uric acid production and the efficiency of renal excretion, affect serum urate level. Elevated serum uric acid level is caused usually by gout, leukemia, diabetes mellitus, hyperfunction of parathyroid and thyroid, renal failure, renal calculus. Urate concentration in serum and in urine depends on glomerular filtration, thus is useful for renal function monitoring.

METHOD PRINCIPLE

Enzymatic, colorimetric method with uricase and peroxidase.

The colour intensity is proportional to the uric acid concentration.

REAGENTS

Package
1-Reagent 2 x 40 ml
2-Reagent 1 x 20.5 ml

The reagents, stored at 2-8°C are stable up to expiry date printed on the package. The reagents stored on board of the analyser at 2-10°C are stable for 10 weeks.

Concentrations in the test

buffer PIPES (pH 7.0) 100 mmol/l
4-aminoantipyrine 0.78 mmol/l
ADPS 0.67 mmol/l
ferricyanide potassium 3.8 μmol/l
peroxidase (POD) > 38.34 μkat/l
uricase > 1.65 μkat/l
ascorbate oxidase > 66.7 μkat/l
sodium hydroxide < 1 %
preservative

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- 1-Reagent meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Warning

H315 Causes skin irritation.
H319 Causes serious eye irritation.
P280 Wear protective gloves/protective clothing/eye protection/face protection.

P302+P352 IF ON SKIN: Wash with plenty of soap and water.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

SPECIMEN

24- hours urine, serum, heparinized plasma free from hemolysis.

Do not use EDTA, fluoride and oxalate as anticoagulants!

Urine preparation: To prevent precipitation of salts of uric acid, 10 ml of NaOH (500 g/L) should be added to the collection bottle before collection of a 24-hour specimen. Urine should be diluted with distilled water in the ratio of 1 to 4 (multiply the result by 5).

Serum and plasma can be stored 3-5 days at 2-8°C or 6 months at -20°C. 24-hours urine samples can be stored approximately 3 days at room temperature.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.

Deionised water is recommended as a reagent blank.

Actions required:

When performing assays in analyser ACCENT-300, there is a probability of **cross-contamination** affecting the tests results: UA PLUS – URINE PROTEINS, FERRUM – UA PLUS, HDL DIRECT – UA PLUS, LDL DIRECT – UA PLUS, URINE PROTEINS – UA PLUS, ALBUMIN – UA PLUS. To avoid this effect follow the recommendations contained in the instruction 51_03_24_007_ACCENT-300_CARRYOVER.

REFERENCE VALUES ⁵

serum / plasma	mg/dl	μmol/l
females	2.5 – 6.8	149 – 405
males	3.6 – 7.7	214 – 458
24-hours urine	mg/24h	mmol/24h
	250 – 750	1.49 – 4.46

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) for determination in serum or CORMAY URINE CONTROL LEVEL 1 (Cat. No 5-161) and LEVEL 2 (Cat. No 5-162) for determination in urine.

For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

The calibration curve should be prepared every 10 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analyser ACCENT-300. Results may vary if a different instrument or a manual procedure is used.

Sensitivity (serum / plasma):

1.11 mg/dl (66.02 μmol/l).

Sensitivity (urine):

1.27 mg/dl (75.54 μmol/l).

Linearity (serum / plasma):

33 mg/dl (1962.8 μmol/l).

Linearity (urine):

68 mg/dl (4044.6 μmol/l).

For higher concentration, dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

Specificity / Interferences

Haemoglobin up to 1.25 g/dl, ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl and triglycerides up to 1000 mg/dl do not interfere with the test.

Precision

Repeatability (run to run) n = 10	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	4.72	0.04	0.91
level 2	9.89	0.08	0.83
Reproducibility (day to day) n = 40	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	4.84	0.10	2.11
level 2	9.02	0.13	1.49

Method comparison

A comparison between uric acid values determined at ACCENT-300 (y) and at ADVIA 1650 (x) using 50 serum samples gave following results:

$$y = 1.0241 x + 0.0959 \text{ mg/dl};$$

$$R = 0.9973 \quad (R - \text{correlation coefficient})$$

A comparison between uric acid values determined at ACCENT-300 (y) and at BS-400 (x) using 44 urine samples gave following results:

$$y = 1.0475 x + 0.3669 \text{ mg/dl};$$

$$R = 0.9996 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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PROGRAM NA ANALIZATOR / APPLICATION / АДАПТАЦИЯ:

Parameters

No.	50	Prim.Wave.	546
Test	UA PLUS	Sec.Wave.	700
Method	Endpoint	Sample Vol.	5
Direction	Ascend	R1 Vol.	200
Unit	mg/dl	R2 Vol.	50
Decimals	2	Line. Limit	
Incubation	10	Antigen Check	
Reaction	-2 35	Substrat	0
R1 Blank		Mix. R Blank	
Lower	0	Lower	0
Upper	0	Upper	0
Response		Linearity	
Lower	-2.5	Lower	1.11
Upper	2.5	Upper	33
Sample	45	Full Name	Uric_Acid
Vol.			PLUS
Dilution	5	Print No.	50

Calibration

Rule	Two Point Linear
Replicates	3
Interval	70
Sensitivity	0
Correlation	0
Difference	2.5
Blank Response	0 2.5
Coefficient	0
Difference	
Non-linear SD	0

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- Henry R.J.: Clinical Chemistry, Harper & Row Publishers Inc., New York (1974).