## ACCENT-300 TG mono

Nr kat. **7-373** (PL)

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(EN)

### INTENDED USE

Diagnostic kit for determination of triglycerides concentration, used in automatic analyser ACCENT-300.


The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Triglycerides are built of glycerol molecule esterified with three fatty acids molecules. Triglycerides are delivered with food or are synthesized endogenously in liver. Triglycerides stored in adipose tissue constitute a reserve of energy. Elevated triglycerides serum level is a risk factor of atherosclerosis. Triglycerides measurement is useful for hyperlipidemia diagnosis and treatment or for estimation of atherosclerosis progression.

### METHOD PRINCIPLE

Colorimetric, enzymatic method with glycerophosphate oxidase.

The colour intensity is proportional to the triglycerides concentration.

### REAGENTS

#### Package

1-Reagent 6 x 40 ml

The reagents, stored at 2-8°C are stable up to expiry date printed on the package. The reagents stored on board of the analyser at 2-10°C are stable for 12 weeks.

#### Concentrations in the test

buffer TRIS (PH 8.0)	200 mmol/l
4-aminoantipyrine (4-AA)	< 0.4 mmol/l
ATP	< 1.5 mmol/l
Mg <sup>2+</sup>	< 1.6 mmol/l
4-chlorophenol	< 2.5 mmol/l
chlorophenicol	1.6 mmol/l
potassium hexacyanoferrate (II)	< 1 mmol/l
FAD-2Na	< 1 mmol/l
glycerol kinase (GK)	~2500 U/l
glycerol phosphate oxidase (GPO)	~2500 U/l
peroxidase (POD)	~1900 U/l
lipoprotein lipase (LPL)	~2000 U/l
detergents, preservatives	

### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

### SPECIMEN

Serum, EDTA or heparinized plasma (recommended: heparine lithium, sodium or ammonium salt) free from hemolysis. Blood should be collected only if the patient has been fasting for minimum of 12 hours. Before blood collection patient should stay in rest position for about 30 minutes. Venous blood is recommended for triglycerides measurement.

Plasma triglycerides values have been reported to be 2% to 4% lower than serum triglycerides values.

Serum and plasma can be stored up to 3 days at 2-8°C or up to 3 months at -20°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-Reagent and 2-Reagent are ready to use.

Deionised water is recommended as a reagent blank.

#### Actions required:

When performing assays in analyser ACCENT-300, there is a probability of **cross-contamination** affecting the tests results: ALBUMIN – TG mono, HDL DIRECT – TG mono, TG mono – MG. To avoid this effect follow the recommendations contained in the instruction 51\_03\_24\_007\_ACCENT-300\_CARRYOVER.

### REFERENCE VALUES <sup>7</sup>

serum, plasma	< 150 mg/dl < 1.7 mmol/l
---------------	-----------------------------

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, CORMAY LIPID CONTROL 1 (Cat. No 5-179) and CORMAY LIPID CONTROL 2 (Cat. No 5-180) or CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173).

For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) or LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

The calibration curve should be prepared every 12 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analyser ACCENT-300. Results may vary if a different instrument or a manual procedure is used.

- Sensitivity:** 5.3 mg/dl (0.06 mmol/l).

- Linearity:** up to 1300 mg/dl (14.69 mmol/l).

For higher concentration dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

### Specificity / Interferences

Haemoglobin up to 0.31 g/dl, bilirubin up to 8.6 mg/dl and ascorbate up to 31 mg/l do not interfere with the test.

### Precision

Repeatability (run to run) n = 20	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	84.69	0.70	0.82
level 2	167.3	1.56	0.93
Reproducibility (day to day) n = 80	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	83.65	2.13	2.54
level 2	165.3	3.79	2.29

### Method comparison

A comparison between triglycerides values determined at **ACCENT-300** (y) and at **COBAS INTEGRA 400** (x) using 99 samples gave following results:

$$y = 0.9216x + 11.054 \text{ mg/dl;}$$

$$R = 0.9968 \quad (R - \text{correlation coefficient})$$

### WASTE MANAGEMENT

Please refer to local legal requirements.

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