A-400 IgM

Nr kat. 7-405 (PL)

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Cat. No **7-405** (EN)

INTENDED USE

Diagnostic kit for determination of IgM concentration used in automatic analysers BS-400 and BS-480. The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Immunoglobulins (Igs) are the instrumental proteins of immunity. Immunity is a property of the lymphoid system which is made of organs (spleen, thymus, bone marrow) and of cells (lymphocytes). Circulating immunoglobulins are secreted in the blood by B lymphocytes and they thereby export far-away the specific biological functions of humoral immunity. Immunoglobulin M (IgM) is the first Ig to appear in response to an antigenic stimulus such as an infectious agent. In many cases, the antigen-specific IgM level subsequently falls and remains low as the IgG response appears.

METHOD PRINCIPLE

The IgM present in a sample form with the specific antibody an immunological complex. The increase of turbidity after the addition of antiserum measured at $\lambda=340$ nm is proportional to IgM concentration in the sample.

REAGENTS

Package
 1-Reagent 1 x 40 ml
 2-Reagent 1 x 9 ml

Buffer (1-Reagent) stored at 2-25°C and antiserum (2-Reagent) stored at 2-8°C are stable until expiry date printed on the package.

Reagent components

1-Reagent: tricine buffer, polymer, inorganic salt, preservative.

2-Reagent: anti-human IgM antiserum, buffer, inorganic salt, preservative.

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Store closed.
- Do not freeze the reagents.
- Nanoparticle-based reagents can settle over time. It may be necessary to delicately mix by repeated turning.
- Human-origin products have been tested for HBsAg and antibodies to HIV 1, HIV 2, HCV and found to be non-reactive. However this material should be handled as potentially infectious.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN ⁴

Serum. Nonhemolyzed, nonlipaemic, fresh serum is recommended.

Serum can be stored up to 3 days at 2-8°C or up to 6 months at -20°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use. For reagent blank 0.9% NaCl is recommended.

REFERENCE VALUES ⁴

adults	0.50 – 3.00 g/l
children (1 year – 12 years)	0.45 – 2.50 g/l
children (1 month – 12 months)	0.20 – 1.50 g/l

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY IMMUNO-CONTROL III (Cat. No 4-291) with each batch of samples.

For the calibration of automatic analysers systems the CORMAY IMMUNO-MULTICAL (Cat. No 4-287) is recommended. 0.9% NaCl should be used as a calibrator 0. The calibration curve should be prepared with every change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

These metrological characteristics have been obtained using an automatic analysers Hitachi or BS-480. Results may vary if a different instrument is used.

- Measurement range:** 0.038 g/l to 7 g/l.

- Specificity / Interferences:**

Hemoglobin up to 0.32 g/dl, bilirubin up to 22 mg/dl, triglycerides up to 312 mg/dl, heparin up to 0.5 g/l, sodium fluoride up to 4 g/l, EDTA up to 5 g/l, sodium citrate up to 5 g/l do not interfere with the test.

- Precision**

Repeatability (run to run) n = 30	Mean [g/l]	SD [g/l]	CV [%]
level 1	0.60	0.004	0.73
level 2	1.25	0.012	0.87
level 3	18.92	0.014	0.72

Reproducibility (day to day) n = 60	Mean [g/l]	SD [g/l]	CV [%]
level 1	0.56	0.024	3.99
level 2	1.30	0.037	2.81
level 3	19.95	0.050	2.49

- Method comparison**

A comparison between IgM values determined at BS-480 (y) and at Advia 1650 (x) using 49 samples gave following results:

$$y = 0.9925x + 0.0416 \text{ g/l};$$

$$R = 0.990 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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• Calibration		Judgment Criteria		Auto QC	
Calibration		Sensitivity		Blank Abs.	
Rule	Logit-Log 5P	Factor Diff.		Error Limit	
Replicate	1	SD		Corr. Coeff.	
• QC		Cum. Sum Check		Interval	
Rules		1.0 - 2.7			
Westgard Multi-rule		• 1.0 - 3.0			
v	1-2S	0.5 - 5.1			
v	1-3S				
v	2-2S				

• **BS-480**

Chem IgM		No.	032	Sample Type SERUM	
Chemistry Immunoglobulin M		Print name IgM		Reaction Direction Increase	
Reaction Type Endpoint		Sec Wave 700		Decimal 0.01	
Pri Wave 340		Reaction Time 81 82		Reagent Vol	
Unit g/l		Diluent		Diluent	
Blank Time 48 49		R1 200 μL		R2 40 μL	
Standard		R3		R4	
Decreased		R1 2.5 μL		Aspirated 20 μL	
Increased		R2 2.5 μL		Diluent 180 μL	
		R3		R4	
		R4		Sample Blank <input type="checkbox"/>	
Linearity Range (Standard)		Linearity Limit		Substrate Depletion	
Linearity Range (Decreased)		Mixed Blank Abs		Uncapping Time	
Linearity Range (Increased)		R1 Blank Abs		Reagent Alarm Limit	
Blank Response		Blank Response		Enzyme Linear Extension	
Twin Chemistry		Prozone Check		Antigen Addition	
		Rate Check		Q3	
		Q1		Q4	
		PC		ABS	

Calibration Settings		Auto Calibration	
Math Model	Logit-Log 5P	<input type="checkbox"/> Bottle Changed	
Factor		<input type="checkbox"/> Lot Changed	
Replicates	1	<input type="checkbox"/> Cal Time	
Acceptance Limits		SD	
Cal Time	672 Hour	Repeatability	
Slope Diff			
Sensitivity			
Deter Coeff			

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