A-400 UA

Nr kat. **7-408** (PL)

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(EN)

INTENDED USE

Diagnostic kit with ascorbate oxidase for determination of uric acid concentration used in automatic analysers BS-400 and BS-480.



The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Uric acid is a product of purine catabolism. It is produced in the liver and excreted in the urine. Both, the amount of uric acid production and the efficiency of renal excretion, affect serum urate level. Elevated serum uric acid level is caused usually by gout, leukemia, diabetes mellitus, hyperfunction of parathyroid and thyroid, renal failure, renal calculus. Urate concentration in serum and in urine depends on glomerular filtration, thus is useful for renal function monitoring.

METHOD PRINCIPLE

Enzymatic, colorimetric method with uricase and peroxidase.

The colour intensity is proportional to the uric acid concentration.

REAGENTS

Package

1-Reagent 2 x 40 ml
2-Reagent 1 x 21.5 ml

The reagents, stored at 2-8°C are stable up to expiry date printed on the package.

Concentrations in the reagent

1-Reagent

ascorbate oxidase ≤ 104 µkat/l
peroxidase (POD) ≤ 22.4 µkat/l
4-aminoantipyrine ≤ 1.2 mmol/l
sodium hydroxide ≤ 0.8 %
buffer PIPES (pH 7.0) ≤ 120 mmol/l
stabilizers, preservatives, detergent

2-Reagent

buffer PIPES (pH 7.0) ≤ 60 mmol/l
ADPS ≤ 2 mmol/l
uricase ≤ 9.9 µkat/l
ferricyanide potassium ≤ 22.8 µmol/l
sodium hydroxide ≤ 0.4 %
stabilizers, preservatives, detergent

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- 1-Reagent meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Warning

H315 Causes skin irritation.

H319 Causes serious eye irritation.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P302+P352 IF ON SKIN: Wash with plenty of soap and water.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

SPECIMEN

24- hours urine, serum, heparinized plasma free from hemolysis.

Do not use EDTA, fluoride and oxalate as anticoagulants! Urine preparation: To prevent precipitation of salts of uric acid, 10 ml of NaOH (500 g/L) should be added to the collection bottle before collection of a 24-hour specimen. Urine should be diluted with distilled water in the ratio of 1 to 4 (multiply the result by 5).

Serum and plasma can be stored 3-5 days at 2-8°C or 6 months at -20°C. 24-hours urine samples can be stored approximately 3 days at room temperature.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.

Deionised water is recommended as a reagent blank.

REFERENCE VALUES ⁵

serum / plasma	mg/dl	µmol/l
females	2.5 – 6.8	149 – 405
males	3.6 – 7.7	214 – 458
24-hours urine	mg/24h	mmol/24h
	250 – 750	1.49 – 4.46

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) for determination in serum or CORMAY URINE CONTROL LEVEL 1 (Cat. No 5-161) and LEVEL 2 (Cat. No 5-162) for determination in urine.

For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended. Calibrator and deionised water should be used for calibration.

The calibration curve should be prepared with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analysers BS-400 and BS-480. Results may vary if a different instrument or a manual procedure is used.

LoB (Limit of Blank):

0.01 mg/dl (0.59 µmol/l) - BS-400, BS-480

LoD (Limit of Detection):

0.04 mg/dl (2.38 µmol/l) - BS-400, BS-480

LoQ (Limit of Quantitation):

0.16 mg/dl (9.52 µmol/l) - BS-400 (serum/ plasma)

0.09 mg/dl (5.35 µmol/l) - BS-400 (urine)

0.08 mg/dl (4.76 µmol/l) - BS-480 (serum/ plasma)

0.08 mg/dl (4.76 µmol/l) - BS-480 (urine)

Linearity:

up to 36 mg/dl (2141.28 µmol/l) - BS-400, BS-480 (serum/ plasma)

up to 56 mg/dl (3330.88 µmol/l) - BS-400 (urine)

up to 61 mg/dl (3628.28 µmol/l) - BS-480 (urine)

For higher concentration, in serum or plasma, dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

Specificity / Interferences

Haemoglobin up to 5 g/dl, ascorbate up to 30 mg/dl for determinations in serum, ascorbate up to 50 mg/dl for determinations in urine, bilirubin up to 20 mg/dl and triglycerides up to 1000 mg/dl do not interfere with the test.

Precision

Repeatability (run to run)		Mean [mg/dl]	SD [mg/dl]	CV [%]
BS-400 (n = 20)	level 1	4.84	0.02	0.39
	level 2	9.11	0.02	0.22
Reproducibility (day to day)		Mean [mg/dl]	SD [mg/dl]	CV [%]
BS-400 (n = 80)	level 1	4.81	0.04	0.9
	level 2	9.51	0.05	0.5

Method comparison

A comparison between uric acid concentration at **BS-400** (y) and at **BECKMAN COULTER AU680** (x) using 60 serum samples gave following results:

$$y = 0.9919x + 0.2168 \text{ mg/dl};$$

R = 0.998 (R – correlation coefficient)

A comparison between uric acid concentration at **BS-400** (y) and at **BECKMAN COULTER AU680** (x) using 30 urine samples gave following results:

$$y = 1.0429x - 0.9941 \text{ mg/dl};$$

R = 0.996 (R – correlation coefficient)

A comparison between uric acid concentration at **BS-400** (y) and at **BECKMAN COULTER AU680** (x) using 33 plasma samples gave following results:

$$y = 0.9742x + 0.2383 \text{ mg/dl};$$

R = 0.997 (R – correlation coefficient)

A comparison between uric acid concentration at **BS-480** (y) and at **BECKMAN COULTER AU680** (x) using 30 serum samples gave following results:

$$y = 0.9954x + 0.2649 \text{ mg/dl};$$

R = 0.999 (R – correlation coefficient)

A comparison between uric acid concentration at **BS-480** (y) and at **BECKMAN COULTER AU680** (x) using 30 urine samples gave following results:

$$y = 0.997x + 0.7884 \text{ mg/dl};$$

R = 0.997 (R – correlation coefficient)

A comparison between uric acid concentration at **BS-480** (y) and at **BECKMAN COULTER AU680** (x) using 30 plasma samples gave following results:

$$y = 0.966x + 0.3531 \text{ mg/dl};$$

R = 0.993 (R – correlation coefficient)

WASTE MANAGEMENT

Please refer to local legal requirements.

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• Calibration

Calibration

Rule	Multi-point Linear
Replicate	2
K	

Judgment Criteria

Sensitivity		Blank Abs.	
Factor Diff.		Error Limit	
SD		Corr. Coeff.	

• QC

Rules

Westgard Multi-rule			
v	1-2S	v	R-4S
v	1-3S	v	4-1S
v	2-2S	v	10-X

Cum. Sum Check: 1.0 - 2.7, 1.0 - 3.0, 0.5 - 5.1

Auto QC

Interval	
----------	--

• **BS-480**

Chem	UA	No.	009	Sample Type	SERUM/URINE
Chemistry	Uric Acid	Print name	UA		
Reaction Type	Endpoint	Reaction Direction	Increase		
Pri Wave	546	Sec Wave	800		
Unit	mg/dl	Decimal	0.01		
Blank Time	48	49	Reaction Time	71	73
Standard	Sample Vol	Aspirated	Diluent	Reagent Vol	Diluent
Decreased	4 μL	20 μL	180 μL	R1 160 μL	
Increased	800 μL			R2 40 μL	
				R3	
				R4	
Linearity (Standard)	Range 0.08	36	Linearity Limit		
Linearity (Decreased)	Range		Substrate Depletion		
Linearity (Increased)	Range		Mixed Blank Abs	-33000	33000
R1 Blank Abs	-33000	33000	Uncapping Time		Day(s)
Blank Response	-33000	33000	Reagent Alarm Limit		
Twin Chemistry			<input type="checkbox"/> Enzyme Linear Extension		
<input type="checkbox"/> Prozone Check	<input type="checkbox"/> Rate Check	<input type="checkbox"/> Antigen Addition			
Q1 0	Q2 0	Q3 0	Q4 0		
PC 0	ABS 0				

Calibration Settings

Math Model	Multi-point Linear		
Factor		Replicates	2

Auto Calibration

<input type="checkbox"/> Bottle Changed
<input type="checkbox"/> Lot Changed
<input type="checkbox"/> Cal Time

Acceptance Limits

Cal Time		Hour
Slope Diff		SD
Sensitivity		Repeatability
Deter Coeff		

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