## A-400 CK

Nr kat. **7-420** (PL)

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Cat. No **7-420** (EN)

### INTENDED USE




Diagnostic kit for determination of creatine kinase activity, intended to use in automatic analysers: BS-400 and BS-480. The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Creatine kinase (CK) catalyzes the transfer of phosphate group between creatine phosphate and adenosine diphosphate (ADP). The product of this reaction is adenosine triphosphate (ATP) – molecular source of energy. CK is a dimer, composed of two different subunits called M and B. Three different isoenzymes formed from these subunits are found in brain and smooth muscle (BB), skeletal muscle (MM) and cardiac muscle (MM and MB). Increased level of CK is usually the result of muscle injury, myocardial or pulmonary infarction.

### METHOD PRINCIPLE

Optimized kinetic method according to International Federation of Clinical Chemistry (IFCC).

The rate of absorbance changes at  $\lambda=340$  nm is directly proportional to creatine kinase activity.

### REAGENTS

<b>Package</b>	
1-Reagent	2 x 40 ml
2-Reagent	1 x 18 ml

The reagents, stored at 2-8°C are stable up to expiry date printed on the package. The reagents stored on board of the analyser at 2-10°C are stable for 12 weeks.

### Concentrations in the test

<b>1-Reagent</b>	
imidazole buffer	100 mmol/l
glucose	20 mmol/l
N-acetylcysteine	20 mmol/l
magnesium acetate	10 mmol/l
EDTA	2 mmol/l
NADP	2 mmol/l
ADP	2 mmol/l
AMP	5 mmol/l
HK	> 2.5 U/ml
<b>2-Reagent</b>	
diadenosinepentaphosphate	10 $\mu$ mol/l
glucose-6-phosphate-dehydrogenase (G6P-DH)	> 1.5 U/ml
creatine phosphate	30 mmol/l
preservative	

### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Do not use reagents past the expiry date.
- Do not freeze reagents.
- Do not interchange caps among reagents.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- 1-Reagent meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Ingredients:

1-Reagent contains imidazole.

### Danger

H360 May damage fertility or the unborn child.  
P201 Obtain special instructions before use.  
P202 Do not handle until all safety precautions have been read and understood.  
P308+P313 IF exposed or concerned:

Get medical advice/attention.

P405 Store locked up.

P501 Dispose of the contents/containers in accordance with the current legislation on waste treatment.

### SPECIMEN

Serum, free from hemolysis.

CK activity is unstable and is rapidly lost during storage. Probes should be stored tightly closed and protected from light.

Specimens can be stored up to 4-8 hours at 15-25 °C or 1-2 days at 2-8°C or 1 month at -20°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-Reagent and 2-Reagent are ready to use.

Deionized water is recommended as a reagent blank.

### Actions required:

When performing assays at analyser **BS-400**, there is a probability of **cross-contamination** affecting the tests results: CK-MB – CK, CK – CREATININE. To avoid this effect follow the recommendations contained in the instruction 51\_03\_24\_002\_BS-400\_CARRYOVER.

### REFERENCE VALUES <sup>8</sup>

Serum	37°C	
female	< 167 U/l	< 2.78 $\mu$ kat/l
male	< 190 U/l	< 3.17 $\mu$ kat/l

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) with each batch of samples.

**For the calibration of automatic analyser system BS-400** the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

**For the calibration of automatic analyser system BS-480** the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended. Deionised water should be used as a calibrator 0.

The calibration curve should be prepared every 12 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analysers: BS-400 and BS-480. Results may vary if a different instrument or a manual procedure is used.

#### Sensitivity

12.3 U/l (0.21  $\mu$ kat/l) – BS-400  
8.94 U/l (0.15  $\mu$ kat/l) – BS-480

#### Linearity

up to 1800 U/l (30.0  $\mu$ kat/l) – BS-400  
up to 1820 U/l (30.3  $\mu$ kat/l) – BS-480

Samples with higher CK activity dilute 1:1 with 0.9% NaCl and repeat the assay. Multiply the result by 2.

#### Specificity / Interferences

Haemoglobin up to 0,156 g/dl, bilirubin up to 20 mg/dl, ascorbate up to 62 mg/l and triglycerides up to 1000 mg/dl do not interfere with the test.

#### Precision

Repeatability (run to run)		Mean [U/l]	SD [U/l]	CV [%]
<b>BS-400</b> (n = 20)	level 1	141.77	1.93	1.36
	level 2	513.47	2.81	0.55
<b>BS-480</b> (n = 10)	level 1	146.32	1.26	0.86
	level 2	463.27	1.58	0.34
Reproducibility (day to day)		Mean [U/l]	SD [U/l]	CV [%]
<b>BS-400</b> (n = 20)	level 1	143.00	2.37	1.66
	level 2	523.00	7.48	1.43
<b>BS-480</b> (n = 10)	level 1	143.29	2.05	1.43
	level 2	451.61	5.96	1.32

#### Method comparison

A comparison between CK values determined at **BS-400** (y) and at **Cobas Integra 400** (x) using 39 samples gave following results:

$$y = 0.9014x + 0.4543 \text{ U/l};$$

$$R = 0.998 \quad (R - \text{correlation coefficient})$$

A comparison between CK values determined at **BS-480** (y) and at **Advia 1650** (x) using 58 samples gave following results:

$$y = 1.0108x + 1.3718 \text{ U/l};$$

$$R = 1.000 \quad (R - \text{correlation coefficient})$$

### WASTE MANAGEMENT

Please refer to local legal requirements.

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**BS-480**

Chem	CK	No.	021	Sample Type	SERUM
Chemistry	CK	Print name	CK		
Reaction Type	Kinetic	Reaction Direction	Increase		
Pri Wave	340	Sec Wave	450		
Unit	U/L	Decimal	0.01		
Blank Time	0 0	Reaction Time	68 78		
Standard	Sample Vol 6 μL Aspirated 6 μL Diluent μL	Reagent Vol	R1 200 μL Diluent μL		
Decreased	6 μL 20 μL 180 μL	R2	40 μL μL		
Increased	μL μL μL	R3	μL μL		
	μL μL μL	R4	μL μL		
Linearity Range (Standard)	8.94 1820	Linearity Limit	0.2		
Linearity Range (Decreased)		Substrate Depletion			
Linearity Range (Increased)		Mixed Blank Abs	-33000 33000		
R1 Blank Abs	-33000 33000	Uncapping Time	84 Day(s)		
Blank Response	-33000 33000	Reagent Alarm Limit			
Twin Chemistry		<input type="checkbox"/> Enzyme Linear Extension			
<input type="checkbox"/> Prozone Check	<input type="checkbox"/> Rate Check	<input type="checkbox"/> Antigen Addition			
Q1	0	Q2	0	Q3	0 Q4 0
PC	0	ABS	0		
<b>Calibration Settings</b>		<b>Auto Calibration</b>			
Math Model	Multi-point Linear	<input type="checkbox"/> Bottle Changed			
Factor		<input type="checkbox"/> Lot Changed			
Replicates	3	<input type="checkbox"/> Cal Time			
<b>Acceptance Limits</b>					
Cal Time	2016 Hour	SD			
Slope Diff		Repeatability			
Sensitivity					
Deter Coeff					

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