A-400 URINE PROTEINS

Nr kat. 7-442 (PL)

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Cat. No 7-442 (EN)

INTENDED USE

Diagnostic kit for determination of total proteins in urine and cerebrospinal fluid intended to use in automatic analysers BS-400 and BS-480.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

In healthy people with properly functioning kidneys, proteins are actively reabsorbed in the proximal tubules and only small amounts of proteins (several mg per day) are excreted in urine. The measurement of total proteins concentration in urine is used in the diagnosis and treatment of heart and thyroid diseases, which are characterized by proteinuria or albuminuria.

The measurement of total proteins concentration in cerebrospinal fluid (CSF) is especially useful in detecting increased permeability of the blood-brain barrier and in detecting increased intrathecal synthesis of immunoglobulins. Increased concentration of protein in CSF may indicate brain tumors, intracerebral hemorrhage, brain injury, bacterial and viral encephalitis and multiple sclerosis.

METHOD PRINCIPLE

Direct, colorimetric method with pyrogallol red.

At an acidic pH the protein amino acid groups, with the pyrogallol red-molybdate complex, form a coloured compound which shows a maximum absorbance at a wavelength of 600 nm. Colour intensity is proportional to the concentration of proteins in the sample.

REAGENTS

Package

1-Reagent 2 x 32.2 ml

The reagent when stored at 15-25°C is stable up to expiry date printed on the package. The reagents stored on board of the analyser at 2-10°C are stable for 12 weeks.

Concentrations in the reagent

Succinate buffer ≤ 60 mmol/l
Pyrogallol red ≤ 0.07 mmol/l
Sodium molybdate ≤ 0.05 mmol/l
Stabilizers, preservative

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Do not freeze the reagent.
- Reagent should be mixed before use by gentle inverting the bottle several times.
- The appearance of turbidity or control urine values outside the manufacturer's acceptable range may indicate of reagent instability.

SPECIMEN

Urine: Urine used for analysis may come from the first morning sample, random sample or timed collection sample according to the classification by ECLM².

In order to collect and prepare the samples only dedicated tubes and containers should be used.

Do not use preservatives.

Freshly collected urine should be kept at room temperature for about an hour and then cooled to 4°C. Direct reduction of the temperature on freshly collected urine can cause precipitation of minerals.

Samples with visible turbidity should be centrifuged before analysis. Determination of uncentrifuged samples may give increased results.

Urine samples are stable for 2 days at 2-8°C⁶.

Nevertheless it is recommended to perform the assay with freshly collected samples!

Cerebrospinal fluid: Cerebrospinal fluid should be centrifuged before analysis. The presence of blood in samples of cerebrospinal fluid may result in false results of protein determination. For proper interpretation of the results, cerebrospinal fluid must be determined simultaneously with a sample of blood taken from a patient at the same time.

The stability of the cerebrospinal fluid is 3 days at 2-8°C, 6 months at the temperature -20°C⁶.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent is ready to use.

0.9% NaCl is recommended as a reagent blank.

Actions required:

BS-400: When performing assays at analyzer BS-400 there is a probability of **cross-contamination** affecting the tests results: CREA ENZYMATIC – URINE PROTEINS II GEN, GGT – URINE PROTEINS II GEN, LDL DIRECT II GEN – URINE PROTEINS II GEN, HDL DIRECT II GEN – URINE PROTEINS II GEN, FERRUM – URINE PROTEINS II GEN, dTIBC – URINE PROTEINS II GEN, URINE PROTEINS II GEN – UA, CREATININE – URINE PROTEINS II GEN, URINE PROTEINS II GEN – UA PLUS, URINE PROTEINS - CREATININE. To avoid this effect follow the recommendations contained in the instruction 51_03_24_002_BS-400_CARRYOVER.

BS-480: When performing assays at analyzer BS-480 there is a probability of **cross-contamination** affecting the tests results: URINE PROTEINS II GEN – dTIBC.

Calculation

For the calculation of proteins excreted over 24 hours, multiply the concentration (mg/dl) by the volume (dl) of the 24 hours urine.

REFERENCE VALUES^{7,10}

urine (adults)	< 15 mg/dl (0.15 g/l)	
urine 24-h (adults)	< 100 mg (0.10 g)	
cerebrospinal fluid	mg/dl	g/l
children 0 - 4 weeks	20 – 80	0.20 – 0.80
children >4 weeks, adults	15 – 45	0.15 – 0.45

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use with each batch of samples the CORMAY URINE CONTROL LEVEL 1 (Cat. No 5-161) and LEVEL 2 (Cat. No 5-162).

For the calibration of automatic analysers the CORMAY URINE PROTEINS CALIBRATORS (Cat. No 5-181) is recommended. 0.9% NaCl should be used as a calibrator 0.

The calibration curve should be prepared every 12 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analysers BS-400 and BS-480. Results may vary if a different instrument or a manual procedure is used.

Sensitivity:

1.6 mg/dl (0.016 g/l) – BS-400

1.8 mg/dl (0.018 g/l) – BS-480

Linearity:

up to 170 mg/dl (1.70 g/l) – BS-400

up to 158 mg/dl (1.58 g/l) – BS-480

For higher concentration dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

Specificity / Interferences

Hemoglobin up to 0.004 g/dl, ascorbic acid up to 20 mg/dl, creatinine up to 6 g/l, bilirubin up to 5 mg/dl, conjugated bilirubin up to 60 mg/dl, uric acid up to 85 mg/dl, glucose up to 35 g/l, citrates up to 250 mg/dl, oxalates up to 90 mg/dl, calcium ions up to 130 mg/dl, magnesium ions up to 1.8 g/l, phosphate ions up to 1.2 g/l, urea up to 50 g/l do not affect the results of the determination.

High concentration of iron (II) ions in the test sample may cause interferences¹².

Acetaminophen and some antibiotics from penicillins and aminoglycosides can interfere^{4,11-13}.

Precision

Repeatability (run to run)		Mean [mg/dl]	SD [mg/dl]	CV [%]
BS-400 n=10	level 1	20.35	0.15	0.74
	level 2	57.29	0.21	0.36
BS-480 n=10	level 1	19.26	0.09	0.49
	level 2	59.10	0.56	0.94
Reproducibility (day to day)		Mean [mg/dl]	SD [mg/dl]	CV [%]
BS-400 n=20	level 1	19.11	0.59	3.09
	level 2	60.58	0.48	0.80
BS-480 n=20	level 1	21.80	0.26	1.20
	level 2	56.88	0.91	1.60

Method comparison

A comparison between total proteins values determined at **BS-400** (y) and at **ADVIA 1650** (x) using 126 urine samples gave following results:

$$y = 1.1263 x - 4.1552 \text{ mg/dl};$$

$$R = 0.997 \quad (R - \text{correlation coefficient})$$

A comparison between total proteins values determined at **BS-400** (y) and at **ADVIA 1650** (x) using 112 CSF samples gave following results:

$$y = 1.0705 x + 0.597 \text{ mg/dl};$$

$$R = 0.997 \quad (R - \text{correlation coefficient})$$

A comparison between total proteins values determined at **BS-480** (y) and at **ADVIA 1650** (x) using 42 urine samples gave following results:

$$y = 1.0163 x - 1.8527 \text{ mg/dl};$$

$$R = 0.988 \quad (R - \text{correlation coefficient})$$

A comparison between total proteins values determined at **BS-480** (y) and at **ADVIA 1650** (x) using 46 CSF samples gave following results:

$$y = 1.0826 x + 2.5936 \text{ mg/dl};$$

$$R = 0.976 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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• Calibration		Judgment Criteria			
Calibration		Sensitivity		Blank Abs.	
Rule	Two-point Linear	Factor Diff.		Error Limit	
Replicate	3	SD		Corr. Coeff.	
• QC					
Rules		Cum. Sum Check		Auto QC	
Westgard Multi-rule		Interval			
<input type="checkbox"/> 1-2S	<input type="checkbox"/> R-4S	1.0 - 2.7			
<input type="checkbox"/> 1-3S	<input type="checkbox"/> 4-1S	• 1.0 - 3.0			
<input type="checkbox"/> 2-2S	<input type="checkbox"/> 10-X	0.5 - 5.1			

BS-480

Chem UP		No. 056		Sample Type URINE/CSF	
Chemistry Urine Proteins		Print name UP			
Reaction Type Endpoint		Reaction Direction Increase			
Pri Wave 605		Sec Wave			
Unit mg/dl		Decimal 0.01			
Blank Time 10 11		Reaction Time 43 44			
Standard Sample Vol 4 μL Aspirated 20 μL Diluent 180 μL		Reagent Vol R1 180 μL Diluent		R2	
Decreased 4 μL 20 μL 180 μL		R3		R4	
Increased		R3		R4	
<input type="checkbox"/> Sample Blank <input type="checkbox"/> Auto Retun					
Linearity Range (Standard) 1.8 158		Linearity Limit			
Linearity Range (Decreased)		Substrate Depletion			
Linearity Range (Increased)		Mixed Blank Abs -33000 33000			
R1 Blank Abs -33000 33000		Uncapping Time 84 Day(s)			
Blank Response -33000 33000		Reagent Alarm Limit			
Twin Chemistry		<input type="checkbox"/> Enzyme Linear Extension			
<input type="checkbox"/> Prozone Check <input type="checkbox"/> Rate Check		• Antigen Addition			
Q1 0 Q2 0		Q3 0 Q4 0			
PC 0 ABS 0					

Calibration Settings		Auto Calibration	
Math Model	Two-point Linear	<input type="checkbox"/>	Bottle Changed
Factor		<input type="checkbox"/>	Lot Changed
Replicates	3	<input type="checkbox"/>	Cal Time
Acceptance Limits			
Cal Time	2016 Hour		
Slope Diff		SD	
Sensitivity		Repeatability	
Deter Coeff			

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