A-400 TG

Nr kat. **7-453** (PL)

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(EN)

INTENDED USE



Diagnostic kit for determination of triglycerides concentration intended to use in automatic analyzers: BS-400 and BS-480. The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Triglycerides are built of glycerol molecule esterified with three fatty acids molecules. Triglycerides are delivered with food or are synthesized endogenously in liver. Triglycerides stored in adipose tissue constitute a reserve of energy. Elevated triglycerides serum level is a risk factor of atherosclerosis. Triglycerides measurement is useful for hyperlipidemia diagnosis and treatment or for estimation of atherosclerosis progression.

METHOD PRINCIPLE

Colorimetric, enzymatic method with glycerophosphate oxidase.

The colour intensity is proportional to the triglycerides concentration.

REAGENTS

Package	
1-Reagent	3 x 40 ml
2-Reagent	1 x 31 ml

The reagents, stored at 2-8°C are stable up to expiry date printed on the package. The reagents stored on board of the analyser at 2-10°C are stable for 11 weeks.

Concentrations in the test

buffer PIPES (pH 7.0)	40 mmol/l
4-aminoantipyrine (4-AA)	0.4 mmol/l
ATP	1.5 mmol/l
Mg ²⁺	1.6 mmol/l
ADPS	0.6 mmol/l
glycerol kinase (GK)	> 66.67 µkat/l
glycerol-3-phosphate oxidase (GPO)	> 60.00 µkat/l
peroxidase (POD)	> 20.00 µkat/l
lipoprotein lipase (LPL)	> 16.67 µkat/l
preservative	

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN

Serum, EDTA or heparinized plasma (recommended: heparine lithium, sodium or ammonium salt) free from hemolysis.

Blood should be collected only if the patient has been fasting for minimum of 12 hours. Before blood collection patient should stay in rest position for about 30 minutes. Venous blood is recommended for triglycerides measurement.

Plasma triglycerides values have been reported to be 2% to 4% lower than serum triglycerides values.

Serum and plasma can be stored up to 3 days at 2-8°C or 3 months at -20°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.

Deionised water is recommended as a reagent blank.

Actions required:

When performing assays in the analyser **BS-400**, there is a probability of **cross-contamination** affecting the tests results: HDL DIRECT II GEN – TG, TG – LIPASE II GEN, TG – LIPASE. To avoid this effect follow the recommendations contained in the instruction 51_03_24_002_BS-400_CARRYOVER.

When performing assays in the analyser **BS-480**, there is a probability of **cross-contamination** affecting the tests results: TG – LIPASE II GEN.

REFERENCE VALUES²:

serum, plasma	< 150 mg/dl
	< 1.7 mmol/l

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, CORMAY LIPID CONTROL 1 (Cat. No 5-179) and CORMAY LIPID CONTROL 2 (Cat. No 5-180) or CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173).

For the calibration of automatic analyser system BS-400 the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

For the calibration of automatic analyser system BS-480 CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended. Deionised water should be used as a calibrator 0.

The calibration curve should be prepared every 11 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analysers BS-400 and BS-480. Results may vary if a different instrument or a manual procedure is used.

Sensitivity

11 mg/dl (0.124 mmol/l) – BS-400
7.0 mg/dl (0.079 mmol/l) – BS-480

Linearity:

up to 2000 mg/dl (22.6 mmol/l) – BS-400
up to 2400 mg/dl (27.1 mmol/l) – BS-480

For higher concentration dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

Specificity / Interferences

Haemoglobin up to 2.5 g/dl, bilirubin up to 20 mg/dl and ascorbate up to 62 mg/dl do not interfere with the test.

Precision

Repeatability (run to run)		Mean [mg/dl]	SD [mg/dl]	CV [%]
BS-400 n=20	level 1	102.93	0.89	0.86
	level 2	511.31	3.48	0.68
BS-480 n=10	level 1	78.42	0.57	0.73
	level 2	178.30	1.31	0.73
Reproducibility (day to day)		Mean [mg/dl]	SD [mg/dl]	CV [%]
BS-400 n=56	level 1	87.28	1.77	2.03
	level 2	266.91	2.98	1.12
BS-480 n=20	level 1	78.47	0.54	0.68
	level 2	176.61	1.59	0.90

Method comparison

A comparison between triglycerides concentration at **BS-400** (y) and at **OLYMPUS AU400** (x) using 41 samples gave following results:

$$y = 1.0208 x - 14.513 \text{ mg/dl};$$

$$R = 0.999 \quad (R - \text{correlation coefficient})$$

A comparison between triglycerides concentration at **BS-480** (y) and at **COBAS INTEGRA 400 Plus** (x) using 39 samples gave following results:

$$y = 1.0355 x - 8.4544 \text{ mg/dl};$$

$$R = 0.997 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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• BS-480

Chem	TG	No.	007	Sample Type	SERUM
Chemistry	TRIGLYCERIDES	Print name	TG		
Reaction Type	Endpoint	Reaction Direction	Increase		
Pri Wave	505	Sec Wave	800		
Unit	mg/dl	Decimal	0.1		
Blank Time	48 49	Reaction Time	69 71		
Standard	Sample Vol 2.5 μL Aspirated 2.5 μL Diluent 180 μL	Reagent Vol	R1 160 μL R2 40 μL R3 μL R4 μL	Diluent	μL μL μL μL
Decreased					
Increased					
	Sample Blank v Auto Retun				
Linearity Range (Standard)	7.0 2400	Linearity Limit			
Linearity (Decreased) Range		Substrate Depletion			
Linearity (Increased) Range		Mixed Blank Abs	-33000 33000		
R1 Blank Abs	-33000 33000	Uncapping Time	77 Day(s)		
Blank Response	-33000 33000	Reagent Alarm Limit			
Twin Chemistry		<input type="checkbox"/> Enzyme Linear Extension			
	<input type="checkbox"/> Prozone Check <input type="radio"/> Rate Check	• Antigen Addition			
Q1	0	Q2	0	Q3	0
PC	0	ABS	0	Q4	0
Calibration Settings		Auto Calibration			
Math Model	Multi-point Linear	<input type="checkbox"/> Bottle Changed			
Factor		<input type="checkbox"/> Lot Changed			
	Replicates 3	<input type="checkbox"/> Cal Time			
Acceptance Limits					
Cal Time	1848 Hour	SD			
Slope Diff		Repeatability			
Sensitivity					
Deter Coeff					

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