

BS-400

Chemistry Analyzer



A-400 VALPROIC ACID

DIAGNOSTIC KIT FOR DETERMINATION OF VALPROIC ACID CONCENTRATION

INTRODUCTION

Valproic acid (2-propylpentanoic acid) is anticonvulsant drug used in all types of epilepsy in children and adults. It is absorbed rapidly and almost completely, however it reaches steady-state after 2-5 days. Over 90% of valproic acid is eliminated by liver metabolism. It pharmacokinetically interacts with other anti-epileptic drugs. The reason for monitoring of valproic acid concentration is lack of correlation between dose and concentration at good correlation between concentration and therapeutic effect and also significant diversity of pharmacokinetics, dependent on age and accompanying anti-epileptic therapy.

METHOD PRINCIPLE

Immunoturbidimetric method; inhibition of agglutination. Increase of absorbance is inversely proportional to the concentration of valproic acid in the sample. Valproic acid, which is present in the sample forms immune complexes with specific antibodies contained in the first reagent. After adding a second reagent containing polystyrene latex particles coated with valproic acid, agglutination is inhibited proportionally to valproic acid concentration in the sample.

REAGENTS

| | |
|----------------|-----------|
| Package | |
| 1-Reagent | 1 x 31 ml |
| 2-Reagent | 1 x 12 ml |

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. The reagents are stable for 4 weeks on board the analyser at 2-10°C. Protect from light and avoid contamination!

Concentrations in the test

Suspension of latex particles coated with valproic acid, buffer with monoclonal antibodies to valproic acid.

Warnings and notes

- Product for in vitro diagnostic use only.
- The reagents should be used by suitably qualified laboratory personnel only in accordance with intended purpose.
- Do not freeze reagents!
- Before use reagent bottles should be shaken gently by inverting several times to dislodge air bubbles and ensure reagent homogeneity.
- Human source material from which this product has been derived has been tested for HBsAg and antibodies to HIV 1 and HIV 2 and found to be non-reactive. However this material and all patient samples should be handled as though capable of transmitting infectious disease.
- The reagents contain sodium azide (< 0.1%) as a preservative. Avoid contact with skin and mucous membranes.
- This assay should not be run immediately after gentamicin, digoxin, phenytoin assays. It is recommended to assay valproic acid in separated batch of samples.

SPECIMEN

Serum.
Samples may be stored up to 3 days at 2-8°C or longer at -20°C. Any additional clotting or precipitation that occurs due to freeze/thawing should be removed by centrifugation prior to analysing the valproic acid of that sample. Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

These reagents may be used in automatic analyser BS-400. 1-Reagent and 2-Reagent are ready to use. For reagent blank CALIBRATOR 1 is recommended.

APPLICATION

BASIC

| Test information | | Reagent volume | |
|------------------|---------------|----------------|-----|
| No. | 80 | R1 | 240 |
| Test | VPA | R2 | 75 |
| Full Name | VALPROIC ACID | R3 | |
| Std. No. | 80 | R4 | |

Sample volume

| | | | |
|-----------|---|----|----|
| Standard | 3 | 15 | 10 |
| Increased | 6 | 15 | 10 |
| Decreased | | | |

Reaction Parameters

| Reac. Type | End-point | Direction | Increase |
|------------|-----------|------------|----------|
| Pri. Wave | 605 | Rgt. Blank | 49 50 |
| Sec. Wave | | Reac. Time | 68 70 |

Result Setup

| | | | |
|---------|-------|-------|---|
| Decimal | 0.01 | Slope | 1 |
| Unit | µg/ml | Inter | 0 |

Judgment Criteria

| | | | | | |
|-------------|---|---|-------------|------|-----|
| Absorbance | 0 | 0 | Lin. Range | 10.1 | 150 |
| Incre. Test | 0 | | Lin. Limit | | |
| Decre. Test | 0 | | Subs. Limit | | |

Prozone Rate Antigen

| | | | | | | | |
|----|---|----|---|-----|---|----|---|
| Q1 | 0 | Q2 | 0 | Q3 | 0 | Q4 | 0 |
| PC | 0 | | | ABS | 0 | | |

CALIBRATION

Calibration

| | |
|-----------|--------------|
| Rule | Logit Log 5P |
| Replicate | 1 |
| K | |

Judgment Criteria

| | |
|--------------|--------------|
| Sensitivity | Blank Abs. |
| Factor Diff. | Error Limit |
| SD | Corr. Coeff. |

THERAPEUTIC RANGE ⁷

| | |
|---------------------------|----------------|
| therapeutic concentration | 50 – 100 µg/ml |
| toxic concentration | > 100 µg/ml |

Some patients achieve the desired therapeutic response at levels outside this range, so it is recommended to consider the need to establish an individual therapeutic ranges for each patient

QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY TDM CONTROLS (Cat. No 5-107) with each batch of samples. For the calibration of automatic analysers systems the CORMAY VALPROIC ACID CALIBRATORS kit (Cat. No 5-151) is recommended.

The calibration curve should be prepared every 7 days, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analysers BS-400 and/or Rx Daytona. Results may vary if a different instrument is used.

▪ Sensitivity

10.1 µg/ml.

▪ Linearity

150 µg/ml.

▪ Specificity / Interferences

Haemoglobin up to 1 g/dl, bilirubin up to 25 mg/dl, triglycerides up to 1000 mg/dl and intralipid up to 800 mg/dl do not interfere with the test.

The specificity of the assay was evaluated by assaying compounds whose chemical structure or concurrent usage could potentially interfere with the valproic acid assay.

| Substance | Concentration (µg/ml) | Cross-Reactivity |
|-------------------------------------|-----------------------|------------------|
| Carbamazepine | 1000 | -0.17% |
| Clonazepam | 100 | 4.01% |
| Diazepam | 100 | 1.30% |
| Ethosuximide | 1000 | 0.18% |
| Phenobarbital | 750 | 0.07% |
| Phenytoin | 1000 | 0.39% |
| Primidone | 1000 | 0.21% |
| 2-N-Propyl-2-pentenoic Acid | 100 | 3.91% |
| 2-N-Propyl-3-pentenoic Acid | 50 | 10.06% |
| 2-N-Propyl-4-pentenoic Acid | 10 | 38.30% |
| 2-Hydroxy-2-N-propylpentanoic Acid | 100 | 3.69% |
| 3-Hydroxy-2-N-propylpentanoic Acid | 10 | 15.20% |
| 5-Hydroxy-2-N-propylpentanoic Acid | 25 | 8.88% |
| 3-Oxo-2-N-propylpentanoic Acid | 50 | 4.38% |
| 4-Oxo-2-N-propylpentanoic Acid | 15 | 26.93% |
| 2-(1'-Propenyl)-2-pentenoic Acid | 50 | 4.12% |
| 3-Oxo-2-(2'-Propenyl)pentanoic Acid | 100 | -1.65% |
| 2-(2'-Propenyl)-4-pentenoic Acid | 25 | 8.76% |
| 2-N-Propylglutaric Acid | 50 | 3.50% |
| 2-N-Propylmalonic Acid | 500 | 0.16% |
| 2-N-Propylsuccinic Acid | 100 | 4.64% |
| Salicylate | 3000 | -0.01% |
| 2-Ethyl-2-phenylmalonamide | 100 | -0.10% |
| 3-Ketovalproic acid | 150 | 0.40% |

▪ Precision

| Repeatability (run to run) n = 10 | Mean [µg/ml] | SD [µg/ml] | CV [%] |
|--------------------------------------|-----------------|---------------|-----------|
| level 1 | 24.7 | 1.57 | 6.34 |
| level 2 | 81.4 | 1.53 | 1.87 |
| level 3 | 136.0 | 2.97 | 2.18 |

| Reproducibility (day to day) n = 20 | Mean [µg/ml] | SD [µg/ml] | CV [%] |
|--|-----------------|---------------|-----------|
| level 1 | 25.9 | 2.14 | 8.27 |
| level 2 | 78.8 | 1.74 | 2.21 |
| level 3 | 137.0 | 3.80 | 2.77 |

▪ Method comparison

A comparison between CORMAY reagent (y) and commercially available assay (x) using 55 samples gave following results:

$$y = 1.06x + 0.35 \text{ µg/ml};$$

$$R = 0.99$$

(R – correlation coefficient)

WASTE MANAGEMENT

Please refer to local legal requirements.

LITERATURE

1. Kodama Y., Koike Y., Kimoto H., Yasuraga F., Takeyama M., Teraoka I., Fujii I.: Binding parameters of valproic acid to serum protein in healthy adults at steady state. *Therapeutic Drug Monitoring* 1992, 14; 55-60.
2. Tisdale J.E., Tsuyuki R.T., Oles K.S., Penry J.K.: Relationship between serum concentrations and dose of valproic acid during monotherapy in adult outpatients. *Therapeutic Drug Monitoring* 1992, 14; 416-423.
3. Chadwick D.W.: Valproate monotherapy in the Management of generalized and partial seizures. *Epilepsia* 1987, 78 (supplement 2): 512-7.
4. Chadwick D.W.: Concentration effect relationships of Valproic Acid. *Clinical Pharmacokinetics* 1985, 10; 155-163.
5. Newman, D.J., Henneberry H., Price C.P.: Particle Enhanced Light Scattering Immunoassay. *Ann. Clin. Biochem.*, 1992, 29; 22-42.
6. Moyer T.P., Shaw L.M., Therapeutic Drugs and Their Management. In: Burtis CA, Ashwood ER, Bruns DE, editors. *Tietz Textbook of Clinical Chemistry and Molecular Diagnostics*. 4th ed, St. Louis: W.B Saunders Company; 2006. p.1253-4; 2315f.
7. Alan H.B. Wu. erditor. *Tietz Clinical Guide to Laboratory Tests*, 4th ed. St. Louis: W.B Saunders Company; 2006. p.1506-7.

Date of issue: 01. 2017.

MANUFACTURER

PZ CORMAY S.A.
22 Wiosenna Street,
05-092 Łomianki, POLAND
tel.: +48 (0) 22 751 79 10
fax: +48 (0) 22 751 79 14
<http://www.cormay.pl>

01/17/01/17