A-400 CREA ENZYMATYC

Nr kat. 7-477 (PL)

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LITERATURA

- Newman DJ, Pnce CP, Renal function and nitrogen metabolites. In: Burtis CA, Ashwood ER, Bruns DE, editors. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 4th ed. St. Louis: W.B Saunders Company; 2006. p. 797-801.
- Alan H.B. Wu. editor. Tietz Clinical Guide to Laboratory Tests, 4th ed. St. Louis: W.B Saunders Company; 2006. p.316.
- Mazzachi BC, Peake M, Erhardt V, Reference range and method comparison for enzymatic and Jaffe Creatinine assays in plasma and serum and early morning urine. Clin Lab 2000; 46: 53-5.
- Susumu Osawa, Medical Technology 1982, Vol. 10, No 7, 575-579.
- Minoru Konno, Medical Technology 1984, Vol. 12, No 3, 270-276.
- Schlebusch H, Liappis N, Klein G. Ultrasensitive CRP and Creatinine: Reference intervals from infancy to childhood. Clin Chem Lab Med. 2001; 39 Special supplement pp S1-S448; May 2001. PO-T042.

Data wydania: 10. 2020.

A-400 CREA ENZYMATIC

Cat. No **7-477** (EN)

INTENDED USE

Diagnostic kit for determination of creatinine concentration, used in automatic analysers BS-400 and BS-480.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Creatinine is a product of creatine nonenzymatic dehydration in skeletal muscle. The amount of creatinine generated and excreted by kidney is proportional to muscle mass and usually is higher in men than women. Daily creatinine generation remains fairly constant, with the exception of crushing injury or degenerative diseases that cause massive damage to muscle. Creatinine blood and urine level depends on glomerular filtration so creatinine clearance is excellent index of renal function.

METHOD PRINCIPLE

Enzymatic, colorimetric method.

creatinine + H₂O $\xrightarrow{\text{creatininase}}$ creatine

creatinine + H₂O $\xrightarrow{\text{creatininase}}$ sarcosine + urea

sarcosine + O₂ + H₂O $\xrightarrow{\text{sarcosine oxidase}}$ glycine + HCHO + H₂O₂

2H₂O₂ + ESPMT + 4-AA $\xrightarrow{\text{peroxidase}}$ quinoneimine dye + 4H₂O

The colour intensity measured at 546 nm is proportional to the creatinine concentration.

REAGENTS

Package

1-Reagent 2 x 40 ml
2-Reagent 1 x 27.5 ml

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. The reagents are stable for 10 weeks on board the analyser at 2-10°C.

Components and concentrations

1-Reagent

Good's buffer ≤ 5%
creatininase ≤ 5%
N-ethyl-N-(3-sulfopropyl)-3-methylaniline (ESPMT) ≤ 5%
sarcosine oxidase ≤ 0.01%
ascorbate oxidase ≤ 1%
detergents, stabilizers and preservatives

2-Reagent

Good's buffer ≤ 5%
creatininase ≤ 1%
peroxidase ≤ 5%
4-ammoantipyrine (4-AA) ≤ 0.01%
stabilizers and preservatives

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN

Serum and urine.

Urine preparation: before analysis urine sample should be diluted with 0,9% NaCl 2-10 times. Multiply the result by dilution factor.

Serum can be stored up to 1 day at 2-8°C. For longer storage samples should be frozen at -20°C.

Urine can be stored up to 1 day at 20-25°C, 4 days at 2-8°C. For longer storage samples should be frozen at -20°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.

For reagent blank deionized water is recommended.

Actions required:

When performing assays at analyser **BS-400**, there is a probability of **cross-contamination** affecting the tests results: CREA ENZYMATIC – URINE PROTEINS II GEN, GLUCOSE - CREA ENZYMATIC. To avoid this effect follow the recommendations contained in the instruction 51_03_24_002_BS-400_CARRYOVER

REFERENCE VALUES^{2,3}

serum / plasma	mg/dl	μmol/l
newborns	0.3 – 1.0	26.5 – 88.4
infants	0.2 – 0.4	17.7 – 35.4
children	0.2 – 0.8	17.7 – 70.7
female	0.5 – 1.0	44.2 – 88.4
male	0.7 – 1.2	61.9 – 106.1
urine (morning)	mg/dl	mmol/l
female	29 – 226	2.56 – 20.0
male	40 – 278	3.54 – 24.6

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use with each batch of samples the following controls: CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) for determination in serum. CORMAY URINE CONTROL LEVEL 1 (Cat. No 5-161) and LEVEL 2 (Cat. No 5-162) for determination in urine. For the calibration of automatic analyser system **BS-400** the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) or LEVEL 2 (Cat. No 5-175; 5-177) is recommended. Deionised water should be used as a calibrator 0.

For the calibration of automatic analyser system **BS-480** the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is

recommended. Deionised water should be used as a calibrator 0.

The calibration curve should be prepared every 8 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

These metrological characteristics have been obtained using automatic analysers BS-400 and BS-480. Results may vary if a different instrument or a manual procedure is used.

Sensitivity:

0.08 mg/dl (7.07 μmol/l) – BS-400

0.09 mg/dl (7.96 μmol/l) – BS-480

Linearity:

up to 25.5 mg/dl (2254 μmol/l) – BS-400

up to 33 mg/dl (2917 μmol/l) – BS-480

Specificity / Interferences

Haemoglobin up to 5 g/dl, ascorbic acid up to 62 mg/l, bilirubin up to 20 mg/dl, triglycerides up to 1000 mg/dl and creatine up to 20 mg/dl do not interfere with the test both at serum and urine.

Precision

Repeatability (run to run)		Mean [mg/dl]	SD [mg/dl]	CV [%]
BS-400 (n = 20)	level 1	0.46	0.01	1.27
	level 2	5.23	0.02	0.46
BS-480 (n = 10)	level 1	1.85	0.03	1.39
	level 2	5.33	0.08	1.44
Reproducibility (day to day)		Mean [mg/dl]	SD [mg/dl]	CV [%]
BS-400 (n = 20)	level 1	1.34	0.04	3.14
	level 2	5.92	0.13	2.15
BS-480 (n = 10)	level 1	1.83	0.02	1.07
	level 2	5.28	0.03	0.62

Method comparison

A comparison between creatinine values determined at **BS-400** (x) and at **Prestige 24i** (y) using 58 serum samples gave following results:

$y = 1.0324x - 0.0073 \text{ mg/dl}$;

$R = 1.000$ (R – correlation coefficient)

A comparison between creatinine values determined at **BS-480** (x) and at **Cobas Integra 400 Plus** (y) using 43 serum samples gave following results:

$y = 0.9596x + 0.0704 \text{ mg/dl}$;

$R = 0.996$ (R – correlation coefficient)

WASTE MANAGEMENT

Please refer to local legal requirements.

LITERATURE

- Newman DJ. Pnce CP. Renal function and nitrogen metabolites. In: Burtis CA. Ashwood ER. Bruns DE. editors. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. 4th ed. St. Louis: W.B Saunders Company; 2006. p. 797-801.
- Alan H.B. Wu. editor. Tietz Clinical Guide to Laboratory Tests. 4th ed. St. Louis: W.B Saunders Company; 2006. p.316.
- Mazzachi BC. Peake M. Erhardt V. Reference range and method comparison for enzymatic and Jaffe Creatinine assays in plasma and serum and early morning urine. Clin Lab 2000; 46: 53-5.
- Susumu Osawa. Medical Technology 1982. Vol. 10. No 7. 575-579.
- Minoru Konno. Medical Technology 1984. Vol. 12. No 3. 270-276.
- Schlebusch H. Liappis N. Klein G. Ultrasensitive CRP and Creatinine: Reference intervals from infancy to childhood. Clin Chem Lab Med. 2001; 39 Special supplement pp S1-S448; May 2001. PO-T042.

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- Newman DJ, Pnce CP, Renal function and nitrogen metabolites. In: Burtis CA, Ashwood ER, Bruns DE, editors. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. 4th ed, St. Louis: W.B Saunders Company; 2006. p. 797-801.
- Alan H.B. Wu. editor. Tietz Clinical Guide to Laboratory Tests, 4th ed. St. Louis: W.B Saunders Company; 2006, p.316.
- Mazzachi BC, Peake M, Erhard V, Reference range and method comparison for enzymatic and Jaffe Creatinine assays in plasma and serum and early morning urine. Clin Lab 2000; 46: 53-5.
- Susumu Osawa, Medical Technology 1982, Vol. 10, No 7, 575-579.
- Minoru Konno, Medical Technology 1984, Vol. 12, No 3, 270-276.
- Schlebusch H, Liappis N, Klein G. Ultrasensitive CRP and Creatinine: Reference intervals from infancy to childhood. Clin Chem Lab Med. 2001; 39 Special supplement pp S1-S448; May 2001. PO-T042.

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• Calibration

Calibration		Judgment Criteria	
Rule	Two-point Linear	Sensitivity	
Replicate	3	Factor Diff.	
K		SD	

• QC

Rules		Auto QC	
Westgard Multi-rule		Cum. Sum Check	Interval
v 1-2S	v R-4S	1.0 - 2.7	
v 1-3S	v 4-1S	• 1.0 - 3.0	
v 2-2S	v 10-X	0.5 - 5.1	

• **BS-480**

Chem	CREA ENZ	No.	078	Sample Type	SERUM/URINE
Chemistry	CREA ENZYMATIC	Print name	CREA ENZ		
Reaction Type	Endpoint	Reaction Direction	Increase		
Pri Wave	546	Sec Wave	800		
Unit	mg/dl	Decimal	0.01		
Blank Time	48 49	Reaction Time	81 82		
Standard	Sample Vol 6 μL Aspirated 33 μL Diluent μL	Reagent Vol	R1 180 μL Diluent μL		
Decreased	6 μL 20 μL 180 μL	R2	60 μL μL		
Increased	μL μL μL	R3	μL μL		
	Sample Blank v Auto Retun	R4	μL μL		
Linearity Range (Standard)	0.09 33	Linearity Limit			
Linearity Range (Decreased)		Substrate Depletion			
Linearity Range (Increased)		Mixed Blank Abs	-33000 33000		
R1 Blank Abs	-33000 33000	Uncapping Time	56 Day(s)		
Blank Response	-33000 33000	Reagent Alarm Limit			
Twin Chemistry		Enzyme Linear Extension			
	Prozone Check o Rate Check	Antigen Addition			
Q1	0	Q2	0	Q3	0
PC	0	ABS	0	Q4	0

Calibration Settings		Auto Calibration	
Math Model	Multi-point Linear		Bottle Changed
Factor		Replicates	3
			Lot Changed
			Cal Time
Acceptance Limits			
Cal Time	672 Hour	SD	
Slope Diff		Repeatability	
Sensitivity			
Deter Coeff			

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