## A-400 BILE ACIDS

Nr kat. 7-484 (PL)

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Cat. No 7-484 (EN)

### INTENDED USE

Diagnostic kit for determination of total bile acids concentration used in automatic analysers BS-400 and BS-480.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.


### INTRODUCTION

Bile acids are the main product of degradation of endogenous cholesterol formed in the liver. Total bile acids are metabolized in the liver and are a valuable indicator of normal or abnormal liver function. Serum total bile acids are increased in patients with viral hepatitis, liver cirrhosis and liver cancer.

### METHOD PRINCIPLE

Enzymatic method with 3- $\alpha$ -hydroxysteroid dehydrogenase (3- $\alpha$ -HSD).

Bile acids under the influence of 3-hydroxysteroid dehydrogenase (3- $\alpha$ -HSD) in the presence of thio-NAD are converted to 3-ketosteroids and thio-NADH. The reaction is reversible and 3- $\alpha$ -HSD can convert 3-ketosteroids and NADH to bile acids and NAD.

The rate of thio-NADH formation can be monitored at 405 nm and is proportional to the bile acids activity.

### REAGENTS

#### Package

1-Reagent 2 x 19 ml  
 2-Reagent 2 x 8 ml

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. The reagents are stable for 9 weeks on board the analyser at 2-10°C.

### Concentrations in the test

#### 1-Reagent

Thio-NAD > 0.1 mmol  
 Buffer

#### 2-Reagent

3- $\alpha$ -HSD > 2 kU/l  
 NADH > 0.1 mmol  
 Buffer

### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Avoid contact with skin and mucous membranes.
- Yellow or yellow-brown color of the reagent does not affect the reagents performance.
- Reagents from different lots must not be interchanged.
- Samples from patients treated with ursodeoxycholic acid (UDCA) are not suitable for the determination of total bile acid concentrations.

### SPECIMEN

Serum.

Total bile acids concentration is increased after meals, therefore samples should be collected under fasting conditions. Serum samples are stable for a 7 days at 4 °C or for 3 month at -20 °C. Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-Reagent and 2-Reagent are ready to use.  
 For reagent blank deionized water is recommended.

### REFERENCE VALUES <sup>3</sup>

serum	2.5 – 6.8 $\mu$ mol/l (1.25 – 3.4 $\mu$ g/ml)
-------	---

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY BILE ACIDS CONTROLS (Cat. No 5-149) with each batch of samples.

For the calibration of automatic analysers systems the CORMAY BILE ACIDS CALIBRATOR (Cat. No 3-125) is recommended. Deionised water should be used as a calibrator 0.

The calibration curve should be prepared every 9 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using analysers BS-400 and/or BS-480 and/or Hitachi 717. Results may vary if a different instrument or a manual procedure is used.

#### Sensitivity:

1.9  $\mu$ mol/l (0.95  $\mu$ g/ml) – BS-400  
 0.98  $\mu$ mol/l (0.49  $\mu$ g/ml) – BS-480

#### Linearity:

up to 140  $\mu$ mol/l (70  $\mu$ g/ml) – BS-400  
 up to 130  $\mu$ mol/l (65  $\mu$ g/ml) – BS-480

For higher concentration, dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

### Specificity / Interferences

Haemoglobin up to 0.5 g/dl, bilirubin up to 50 mg/dl, ascorbic acid up to 50 mg/dl and triglycerides up to 750 mg/dl do not interfere with the test.

### Precision

Repeatability (run to run) n = 10	Mean [ $\mu$ mol/l]	SD [ $\mu$ mol/l]	CV [%]
level 1	31.42	0.13	0.42
level 2	88.13	0.23	0.26
Reproducibility (day to day) n = 10	Mean [ $\mu$ mol/l]	SD [ $\mu$ mol/l]	CV [%]
level 1	30.78	0.92	2.97
level 2	85.83	1.10	1.28

### Method comparison

A comparison between bile acids values determined at BS-400 (y) and at OLYMPUS AU 400 (x) using 55 samples gave following results:

$$y = 1.0487x - 0.098 \mu\text{mol/l};$$

$$R = 1.000 \quad (R - \text{correlation coefficient})$$

A comparison between bile acids values determined at BS-480 (y) and at BS-800 (x) using 50 samples gave following results:

$$y = 1.0064x + 0.2202 \mu\text{mol/l};$$

$$R = 0.998 \quad (R - \text{correlation coefficient})$$

### WASTE MANAGEMENT

Please refer to local legal requirements.

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<b>• Calibration</b>	
<b>Calibration</b>	<b>Judgment Criteria</b>
Rule <input type="text" value="Two-point Linear"/>	Sensitivity <input type="text"/>
Replicate <input type="text" value="3"/>	Factor Diff. <input type="text"/>
K <input type="text"/>	SD <input type="text"/>
<b>• QC</b>	Blank Abs. <input type="text"/>
<b>Rules</b>	Error Limit <input type="text"/>
Westgard Multi-rule	Corr. Coeff. <input type="text"/>
<input type="checkbox"/> 1-2S	Cum. Sum Check <input type="text"/>
<input type="checkbox"/> 1-3S	<input type="checkbox"/> 1.0 - 2.7
<input type="checkbox"/> 2-2S	<input checked="" type="checkbox"/> 1.0 - 3.0
	<input type="checkbox"/> 0.5 - 5.1
	Interval <input type="text"/>
	<b>Auto QC</b>

### BS-480

Chem <input type="text" value="Bile Acids"/>	No. <input type="text" value="085"/>	Sample Type <input type="text" value="SERUM"/>
Chemistry <input type="text" value="Bile Acids"/>	Print name <input type="text" value="Bile Acids"/>	
Reaction Type <input type="text" value="Fixed-time"/>	Reaction Direction <input type="text" value="Increase"/>	
Pri Wave <input type="text" value="412"/>	Sec Wave <input type="text" value="660"/>	
Unit <input type="text" value="umol/L"/>	Decimal <input type="text" value="0.01"/>	
Blank Time <input type="text" value="0"/> <input type="text" value="0"/>	Reaction Time <input type="text" value="68"/> <input type="text" value="78"/>	
Standard <input type="text" value="5"/> <input type="text" value="5"/> <input type="text" value="5"/>	Reagent Vol <input type="text" value="210"/> <input type="text" value="70"/>	Diluent <input type="text"/>
Decreased <input type="text" value="5"/> <input type="text" value="20"/> <input type="text" value="180"/>	R1 <input type="text" value="210"/> <input type="text" value="70"/>	R2 <input type="text" value="70"/>
Increased <input type="text" value="5"/> <input type="text" value="20"/> <input type="text" value="180"/>	R3 <input type="text" value="70"/>	R4 <input type="text" value="70"/>
<input type="checkbox"/> Sample Blank <input type="checkbox"/> Auto Retun		
Linearity Range (Standard) <input type="text" value="0.98"/> <input type="text" value="130"/>	Linearity Limit <input type="text"/>	
Linearity Range (Decreased) <input type="text"/>	Substrate Depletion <input type="text"/>	
Linearity Range (Increased) <input type="text"/>	Mixed Blank Abs <input type="text" value="-33000"/> <input type="text" value="33000"/>	
R1 Blank Abs <input type="text" value="-33000"/> <input type="text" value="33000"/>	Uncapping Time <input type="text" value="63"/> Day(s)	
Blank Response <input type="text" value="-33000"/> <input type="text" value="33000"/>	Reagent Alarm Limit <input type="text"/>	
Twin Chemistry <input type="text"/>	<input type="checkbox"/> Enzyme Linear Extension	
<input type="checkbox"/> Prozone Check <input type="checkbox"/> Rate Check	<b>• Antigen Addition</b>	
Q1 <input type="text" value="0"/> <input type="text" value="0"/>	Q3 <input type="text" value="0"/> <input type="text" value="0"/>	Q4 <input type="text" value="0"/>
PC <input type="text" value="0"/> <input type="text" value="0"/>	ABS <input type="text" value="0"/>	

<b>Calibration Settings</b>	<b>Auto Calibration</b>
Math Model <input type="text" value="Two-point Linear"/>	<input type="checkbox"/> Bottle Changed
Factor <input type="text"/>	<input type="checkbox"/> Lot Changed
Replicates <input type="text" value="3"/>	<input type="checkbox"/> Cal Time
<b>Acceptance Limits</b>	
Cal Time <input type="text" value="1512"/> Hour	
Slope Diff <input type="text"/>	SD <input type="text"/>
Sensitivity <input type="text"/>	Repeatability <input type="text"/>
Deter Coeff <input type="text"/>	

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