

**AccuDiag™**  
**Echinococcus IgG**  
**ELISA Kit**

Cat # 8202-35



Test	Echinococcus IgG ELISA
Method	Enzyme Linked Immunosorbent Assay
Principle	Sandwich Complex
Detection Range	Qualitative : Positive, Negative Control
Sample	5 µL serum
Total Time	~ 20 min.
Shelf Life	12 Months from the manufacturing date
Specificity	100%
Sensitivity	100%

**INTENDED USE**

The **Echinococcus ELISA Test** is an enzyme linked immunosorbent assay (ELISA) for the qualitative identification of IgG antibodies to *Echinococcus* in serum.

**SUMMARY AND EXPLANATION**

Echinococcosis (hydatidosis) is the infection caused by cestodes of the genus *Echinococcus*. Humans are potential intermediate hosts and can become infected by ingesting eggs passed in the feces of an infected animal. The resulting disease is called hydatidosis, or hydatid disease.

Four species are known pathogens of the disease: *E. granulosus*, *E. multilocularis*, *E. oligarthrus* and *E. vogeli*. The infection caused by *E. granulosus* is referred to as cystic hydatid disease (CHD) and results in cysts in various organs, especially the liver and lungs. These cysts may become quite large and contain hundreds or thousands of scoleces called hydatid sand. The degree of antibody response to these cysts will vary depending on their location and degree of calcification. Liver cysts typically produce a higher antibody response than lung cysts. Infection due to *E. multilocularis* is referred to as alveolar hydatid disease (AHD), and also occurs as cysts that may spread throughout the infected tissue.

Since *Echinococcus* eggs are not shed by infected humans, serological determination has been important in the diagnosis of hydatid disease. A number of tests have been used, including latex agglutination (LA), indirect hemagglutination (IHA), complement fixation (CF), agar gel diffusion (AGD) and enzyme linked immunosorbent assay (ELISA).

**TEST PRINCIPLE**

The principle of the Echinococcus ELISA test involves three incubation steps. Before the first incubation, the microwells are coated with Echinococcus cyst antigen. The patient's sera are added, and if there are any antibodies present, they

will bind to the wells during the first incubation. Next, the wells must be washed of any test sample, and added to the wells at this point is the Enzyme Conjugate. During the second incubation, the Enzyme Conjugate will bind to any antibodies present. Before the third incubation step, more washing are necessary. Then a chromogen (tetramethylbenzidine or TMB) is added. With the presence of the Enzyme Conjugate and the peroxidase causing the consumption of peroxide, the chromogen changes to a blue color. The blue color turns to a bright yellow color after the addition of the stop solution, which ends the reaction. ELISA readers can be used to obtain results, or the results can be read visually.

**SPECIMEN COLLECTION AND PREPARATION**

Serological specimens should be collected under aseptic conditions. Coagulate blood and remove serum. Hemolysis is avoided through prompt separation of serum from the clot. Serum should be stored at 2-8°C if it is to be used within a few days. Serum may be held for 3 to 6 months by storage at -20° C or lower. Lipemic and strongly hemolytic serum should be avoided. Do not heat inactivate serum and avoid repeated freezing and thawing of samples.

Test samples: Make a 1:64 dilution of patients' sera using the dilution buffer (e.g. 5 µl sera and 315 µl dilution buffer).

**MATERIALS AND COMPONENTS**

**Materials provided with the test kits**

1. **Plate:** Microwells containing Echinococcus antigens - 96 test wells in a test strip holder.
2. **Enzyme Conjugate:** One (1) bottle containing 11 ml of Protein A conjugated to peroxidase.
3. **Positive Control:** One (1) vial containing 1 ml of diluted positive rabbit serum.
4. **Negative Control:** One (1) vial containing 1 ml of diluted negative human serum.
5. **TMB Substrate Solution:** One (1) bottle containing 11 ml of the chromogen tetramethylbenzidine (TMB).
6. **Wash Concentrate 20X:** One (1) bottle containing 25 ml of concentrated buffer and surfactant.
7. **Dilution Buffer:** Two (2) bottles containing 30 ml of buffered protein solution.
8. **Stop Solution:** One (1) bottle containing 11 ml of 0.73 M phosphoric acid.

**Materials required but not provided**

1. Pipettes
2. Squeeze bottle for washing strips (narrow tip is recommended)
3. Reagent grade water and graduated cylinder
4. Tubes for sample dilution
5. Sample Dilution Tubes
6. Absorbent paper
7. ELISA plate reader with a 450 nm and a 620-650 nm filter (optional if results are read visually).

**Preparation**

Wash Buffer - Remove cap and add contents of bottle to 475 ml of reagent grade water. Place diluted wash buffer into a squeeze bottle with a narrow tip opening.

Note: Washings consist of filling to the top of each well, shaking out the contents and refilling.

Avoid generating bubbles in the wells during the washing steps.

**ASSAY PROCEDURE**

1. Break off number of wells needed (two for controls plus number of samples) and place in strip holder.
2. Add 100 µl (or two drops) of the negative control to well #1, 100 µl of the positive control to well #2 and 100 µl of the diluted (1:64) test samples to the remaining wells.  
 Note: Negative and positive controls are supplied prediluted. Do not dilute further.
3. Incubate at room temperature (15 to 25 °C) for 10 minutes.

4. Shake out contents and wash 3 times with the diluted wash buffer.
5. Add 100µl of Enzyme Conjugate to each well.
6. Incubate at room temperature for 5 minutes.
7. Shake out contents and wash 3 times with wash buffer.
8. Slap wells vigorously against paper towels to remove excess moisture.
9. Add 100 µl of the Chromogen to every well.
10. Incubate at room temperature for 5 minutes.
11. Add 100 µl of the Stop Solution and mix by tapping strip holder.

## RESULTS

**Visually:** Look at each well against a white background (e.g. paper towel) and record as clear or +, ++ or +++ reaction.

**ELISA Reader:** Zero reader on air. Set for bichromatic readings at 450/620-650 nm.

## Troubleshooting

Negative control has excessive color after development.

**Reason:** inadequate washings

**Correction:** wash more vigorously. Remove excessive liquid from the wells by tapping against an Absorbent towel. Do not allow test wells to dry out.

## Interpretation of Results - ELISA Reader

Zero ELISA reader on air. Read all wells at 450/650-620 nm.

**Positive** - Absorbance reading equal to or greater than 0.3 OD units.

**Negative** - Absorbance reading less than 0.3 OD units.

A negative OD reading indicates that the patient has no detectable level of antibodies. This may be due to lack of infection or poor immune response by the patient.

## Interpretation of Results -Visual

Compare results to the controls. A sample should be interpreted as positive if the degree of color development is significant and obvious.

## QUALITY CONTROL

The use of controls allows validation of kit stability. The kit should not be used if any of the controls are out of range.

Expected values for the controls are:

**Negative** - 0.0 to 0.3 OD units

**Positive** - 0.5 OD units and above

## PERFORMANCE CHARACTERISTICS

		Reference Method *	
		+	-
Diagnostic Automation, Inc.	+	18	0
	-	0	65

Positive Agreement: 100% (18/18)

Negative Agreement: 100% (65/65)

\*Reference Method refers to a commercially available ELISA.

## LIMITATIONS OF PROCEDURE

Serologic results are an aid in diagnosis but cannot be used as the sole method of diagnosis.

## EXPECTED VALUES

The number of individuals showing positive results can vary significantly between populations and geographic regions. If possible, each laboratory should establish an expected range for its patient population.

## PRECAUTIONS



1. Do not use solutions if they precipitate or become cloudy. Wash concentrate may show crystallization upon storage at 2-8 °C. Crystallization will disappear after dilution to working strength.
2. Do not use serum that may have supported microbial growth, or is cloudy due to high lipid content. Samples high in lipids should be clarified before use.
3. Treat all sera as if capable of being infectious. Negative control has been tested and found negative for Hepatitis B surface antigen and for the antibody to HIV by required test methods. This product should be used under appropriate safety conditions that would be used for any potentially infectious agent.
4. Do not add azides to the samples or any of the reagents.

## STORAGE

1. Reagents, strips and bottled components should be stored at 2-8 °C
2. Squeeze bottle containing diluted wash buffer may be stored at room temperature.

## REFERENCES

1. Schantz, P. and B. Gottstein. "Echinococcosis (Hydatidosis)." *Immunodiagnosis of Parasitic Diseases, Vol. 1, Helminthic Diseases*. Ed. Walls and Schantz. Academic Press, 1986. pp. 69-107
2. Krogstad D., G. Visvesvara R. Walls, J. Smith. "Tissue Helminths." *Manual of Clinical Microbiology*, 4th Ed.. American Society for Microbiology, Washington, DC. 1985, pp. 657-659
3. Kagan, I. "Serodiagnosis of Parasitic Diseases." *Manual of Clinical Laboratory Immunology*. 3rd Ed.. American Society for Microbiology, Washington, DC. 1986, p. 472
4. Ash, L. and T. Orihel. "Parasites - A Guide to Laboratory Procedures and Identification." 1987, ASCP Press, pp. 176-177
5. Coltorti, E. "Standardization and Evaluation of an Enzyme Immunoassay as a Screening Test for the Seroprevalence of Human Hydatidosis." *Am J Trop Med Hyg* #35(5), 1986, pp. 1000-1005
6. Coltorti, E. et. al. "Field Evaluation of an Enzyme Immunoassay for Detection of Asymptomatic Patients in a Hydatid Control Program." *Am J Trop Med Hyg*, #38(3), 1988, pp. 603-607

<p><b>ISO 13485</b> <b>ISO 9001</b></p> 			
 <p><b>Diagnostic Automation/Cortez Diagnostics, Inc.</b> 23961 Craftsman Road, Suite E/F, Calabasas, California 91302 USA</p>			
<b>Date Adopted</b>	<b>Cat # 8202-35</b>		
<b>2013-10-09</b>	<b>AccuDiag™- Echinococcus IgG ELISA -2013</b>		
<table border="1" style="display: inline-table;"> <tr> <td>EC</td> <td>REP</td> </tr> </table>	EC	REP	<p><b>CEpartner4U, Esdoornlaan 13, 3951DB Maarn. The Netherlands.</b> <a href="http://www.cepartner4u.eu">www.cepartner4u.eu</a></p>
EC	REP		
Revision C Date: 12-27-2013			