## CRP ULTRA

Nazwa zestawu	(PL) Nr kat.
CORMAY CRP ULTRA	6-300
CORMAY CRP ULTRA 60	3-339
CORMAY CRP ULTRA 120	3-340
OS-CRP ULTRA	9-437
B50-CRP ULTRA	5-553

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### INTENDED USE

Diagnostic kit for determination of C-reactive protein concentration intended to use in several automatic analyzers. The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

CRP (C-reactive protein) is an acute phase protein whose concentration is seen to increase as a result of the inflammatory process, most notably in response to pneumococcal (bacterial) infectious, histolytic disease and a variety of disease states. CRP to be used as a marker or general diagnostic indicator of infections and inflammation, in addition to serving as a monitor of patient response to therapy and surgery. Furthermore, regular measurements of CRP in infants can be a useful aid in the early diagnosis of infectious disease.

### METHOD PRINCIPLE

When an antigen-antibody reaction occurs between CRP in a sample and anti-CRP antibody which has been sensitized to latex particles, agglutination results. This agglutination is detected as an absorbance change (572 nm), with the magnitude of the change being proportional to the quantity of CRP in the sample. The actual concentration is then determined by interpolation from a calibration curve prepared from calibrators of known concentration.

### REAGENTS

Package	CORMAY CRP ULTRA	CORMAY CRP ULTRA 60	CORMAY CRP ULTRA 120
1-REAGENT	1 x 46.5 ml	3 x 50 ml	3 x 100 ml
2-REAGENT	1 x 48.5 ml	3 x 50 ml	3 x 100 ml
	<b>OS-CRP ULTRA</b>	<b>B50-CRP ULTRA</b>	
1-REAGENT	2 x 41 ml	2 x 18.5 ml	
2-REAGENT	2 x 41 ml	2 x 18.5 ml	

The reagents are stable up to the kit expiry date printed on the package when stored at 2-10°C. The reagents are stable for 11 weeks on board the analyser at 2-10°C.

### Concentrations in the test

suspension of latex particles sensitized with anti-CRP antibodies (rabbit) (pH 7.3) 0.20 w/v%  
glycine buffer solution (pH 7.0)  
preservative

### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- After measurements are taken, reagent bottles should be capped and kept at 2-10°C. Care should be taken not to interchange the caps of reagent bottles.
- Reagents with different lot numbers should not be interchanged or mixed.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

### SPECIMEN

Serum or plasma (Na-EDTA, K-EDTA, Na-Heparin, Li-Heparin, citrate).

If the test cannot be done immediately, the sample should be placed in a tightly sealable container and stored at -20°C. Repeated freezing and thawing should be avoided. Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

The reagents are ready to use.  
Applications for analysers are available on request.

### REFERENCE VALUES

serum, plasma	
adults	< 0.5 mg/dl (< 5 mg/l)
children (2 months – 15 years)	0.01 – 0.28 mg/dl (0.1 – 2.8 mg/l)
newborns (0 – 3 weeks)	0.01 – 0.41 mg/dl (0.1 – 4.1 mg/l)

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY IMMUNO-CONTROL I (Cat. No 4-288) with each batch of samples.

For the calibration of automatic analysers systems the CORMAY CRP ULTRA CALBRATORS kit (Cat. No 4-276) is recommended. Calibrators and 0.9% NaCl should be used for calibration.

The calibration curve should be prepared every 3 weeks, with every change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analysers Hitachi 912 and Hitachi 917. Results may vary if a different instrument is used.

- Sensitivity:** 0.1 mg/dl
- Linearity:** up to 32 mg/dl.

For higher concentrations dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

### Specificity / Interferences

Haemoglobin up to 0.5 g/dl, bilirubin up to 30 mg/dl, triglycerides up to 500 mg/dl and RF up to 500 IU/ml do not interfere with the test.

### Precision

Repeatability (run to run) n = 20	Mean	SD	CV
	[mg/dl]	[mg/dl]	[%]
level 1	3.93	0.06	1.45
level 2	10.75	0.22	2.00
Reproducibility (day to day) n = 21	Mean	SD	CV
	[mg/dl]	[mg/dl]	[%]
level 1	0.047	0.003	6.97
level 2	0.218	0.007	3.34
level 3	0.976	0.012	1.23

### Method comparison

A comparison between CRP values determined at Hitachi 912 (y) and at Cobas Integra 400 (x) using 39 samples gave following results:

$$y = 1.0274x + 0.0479 \text{ mg/dl};$$

$$R = 0.9972 \quad (R - \text{correlation coefficient})$$

### WASTE MANAGEMENT

Please refer to local legal requirements.

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