Lyophilised Control Serum

**POTASSIUM CONTROL SET
(Assayed Control Serum)**

2 levels

Human based serum for quality control of quantitative in vitro determination of potassium in human serum on photometric systems.

REF

Cont.

910190 2 x 3 ml control set, 2 levels

Additionally offered:
910110 5 x 10 ml Potassium reagent kit

CONTROL COMPOSITION

The controls consist of potassium in human serum.

CONTROL PREPARATION

Reconstitute each vial with 3.0 mL distilled water; allow standing at room temperature for 20 minutes, and mixing gently by swirling until thoroughly dissolved.

CONTROL STABILITY AND STORAGE

| | | |
|-----------------------|--------------|-----------------------|
| | temperature: | Stability: |
| Lyophilized: | 2 – 8 °C | up to expiration date |
| After reconstitution: | 2 – 8 °C | 2 weeks |

CLOSE IMMEDIATELY AFTER USE
PROTECT FROM LIGHT!

WARNINGS AND PRECAUTIONS

Each donor unit of serum in the preparation of this material was tested by FDA-approved methods and found negative for the Human Immunodeficiency Virus Antibody (HIV I/II Ab), Hepatitis B Surface Antigen (HBsAg), and Hepatitis C Virus Antibody (HCV). Because no method can offer complete assurance as to the absence of infectious agents, however this material and all patient samples should be handled as though capable of transmitting infectious disease and disposed of accordingly.

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