α-1 Acid Glycoprotein

Diagnostic reagent for the quantitative in vitro determination of α-1 Acid Glycoprotein in human serum by turbidimetric assay

REF	Content
A00502	1x 10 mL α-1 Acid Glycoprotein Antibody Reagent 5x 25 mL PEG6 Buffer
Additionally offered:	
A00704	5x 1 mL Protein Calibrator 5 Level Series
A00580	1x 1 mL Protein Calibrator High
A00703	1x 5 mL Protein Calibrator High
A00701	1x 1 mL Protein Calibrator Low
A00702	1x 5 mL Protein Calibrator Low
A00590	1x 1 mL Protein Control
A00800	1x 5 mL Protein Control
A08591	1x 1 mL Protein Control Low
A08823	1x 5 mL Protein Control Low

GENERAL INFORMATION

Method	Immunoturbidimetric
Reaction	Nonlinear, endpoint
Wavelength	340 nm
Assay Temperature	18 – 37 °C
Sample	Serum
Measuring Range	approx. 0 – 300 mg/dL
Sensitivity	4 mg/dL (Cobas Mira)
Hook Effect	without sample dilution: > 1300 mg/dL (Cobas Mira) with sample dilution: > 1700 mg/dL (Cobas Mira)
Manual Test Procedure	Tests/Kit*
without sample dilution	166
with sample dilution	250

Automated Test Procedure

Instrument dependent – please ask for applications

* calculated on amount of antibody reagent; additional buffer on request

REAGENT COMPOSITION

COMPONENTS FINAL CONCENTRATION

α-1 Acid Glycoprotein Antibody Reagent	
Turbidimetric grade antibody raised in the goat, monospecific for α-1 Acid Glycoprotein	variable
Sodium azide	0.095 %

PEG6 Buffer

Phosphate buffered saline	
PEG	6 %
Sodium azide	0.095 %

REAGENT PREPARATION

The reagents are liquid and ready to use.

REAGENT STABILITY AND STORAGE

Conditions:	Protect from light. Close immediately after use.
Stability:	at 2 – 8 °C up to the expiration date at 18 – 25 °C 1 month

Do not freeze!

SAMPLE STABILITY AND STORAGE

Stability:	at 2 – 8 °C 48 hours at – 20 °C 3 months
-------------------	---

Freeze only once!

MANUAL TEST PROCEDURE

Test Procedure without Sample Dilution:

Samples/Controls: ready to use

Calibration curve: Use Protein Calibrator High to generate a calibration curve by making 1:2 serial dilutions of the calibrator with 0.9% saline as diluent or use the 5 level calibrator series. Use 0.9% saline as zero point.

Pipette into test tubes:	Calibrators	Samples/Controls
Buffer	900 µL	900 µL
Cal./Ctrls/Samples	5 µL	5 µL
Mix. Read A1 of calibrators and samples/controls at 340 nm. Then add:		
Antibody Reagent	60 µL	60 µL
Mix. Incubate 5 minutes at assay temperature. Read A2 of calibrators and samples/controls at 340 nm. Calculate: $\Delta A = (A2 - A1)$		

Test Procedure with Sample Dilution:

Sample/ Control: dilute 1:10 in saline 0.9%

Calibration curve: Use Protein Calibrator High to generate a calibration curve by making 1:10, 1:20, 1:40, 1:80, 1:160 dilutions with 0.9% saline as diluent. Use 0.9% saline as zero point.

Pipette into test tubes	Calibrators	Samples/Controls
Buffer	900 µL	900 µL
Cal./Ctrls/Samples	30 µL	30 µL
Mix. Read A1 of calibrators and samples/controls at 340 nm. Then add:		
Antibody Reagent	40 µL	40 µL
Mix. Incubate 5 minutes at assay temperature. Read A2 of calibrators and samples/controls at 340 nm. Calculate: $\Delta A = (A2 - A1)$		

CALCULATION

Calculate and plot $\Delta A = (A2 - A1)$ of the calibrators versus assigned concentration values on a linear-linear graph paper. Calculate ΔA optical densities of samples and control(s) and read values in mg/dL on the reference curve. Samples yielding absorbances above highest calibrator should be retested after further dilution.

REFERENCE RANGE

Men	50 – 130 mg/dL
Women	40 – 120 mg/dL

It is recommended that each laboratory establishes its own normal range.

TEST PRINCIPLE

The assay of α-1 Acid Glycoprotein is based on turbidimetric measurement. Turbidity is caused by the formation of antigen-antibody insoluble immuno complexes.

DIAGNOSTIC IMPLICATIONS

As an early acute phase reactant, α-1 Acid Glycoprotein is especially useful in monitoring tumour recurrence. Levels are also helpful in differentiating acute phase responses (elevated levels) from estrogen effects (normal or depressed levels). In addition it is an excellent protein to assay along with Haptoglobin in assessing in vivo hemolysis. An elevated α-1 Acid Glycoprotein level but normal Haptoglobin suggests an acute phase response with mild to moderate in vivo hemolysis.

PERFORMANCE CHARACTERISTICS

SENSITIVITY

4 mg/dL (Cobas Mira)

ACCURACY

Control	Assigned Value (mg/dL)	Measured Value (mg/dL)
Immunology 1 (Ciba Corning)	57 (46 - 68)	57.2
Immunology 2 (Ciba Corning)	121 (97 - 145)	125.5
Liquicheck 1 (Biorad)	46 (37 - 55)	44.4
Liquicheck 2 (Biorad)	110 (88 - 132)	102.3
Seronorm L (Nycomed)	35 (28 - 42)	46.6
Seronorm N (Nycomed)	94 (75 - 113)	90.3
Seronorm H (Nycomed)	147 (118 - 176)	134.0

PRECISION

Intra-Assay Precision

3 serum samples were consecutively measured 20x on the Cobas Mira.

Expected Value	n	Mean	S.D.	C.V.
Low	20	27.9	1.30	4.66
Medium	20	71.2	0.81	1.14
High	20	149.8	3.67	2.45

Inter-Assay Precision

A control was measured daily on the Cobas Mira after calibration. The variation coefficient was calculated..

Control	n	Mean	S.D.	C.V.
Dialab	27	78.3	1.99	2.55

METHOD COMPARISON

A comparison with Nephelometry gave the following results:
 $y = 1.554x - 9.1048$; $r = 0.9958$

INTERFERING SUBSTANCES

No interference up to:

Sodium-Citrate	1000 mg/dL	Hemoglobin	1000 mg/dL
Triglycerides	2500 mg/dL	Heparin	50 mg/dL
Bilirubin	20 mg/dL		

QUALITY CONTROL

All commercially available Control sera with α-1 Acid Glycoprotein values measured by this method may be used. We recommend the Dialab Protein Control and the Protein Control Low.

CALIBRATION

The assay requires the use of α-1 Acid Glycoprotein serum Calibrators. We recommend the Dialab Protein Calibrator 5 Level Series, the Protein Calibrator High or the Protein Calibrator Low.

AUTOMATION

Applications for automated systems (with and without sample dilution) are available upon request.

WARNINGS AND PRECATIONS

- The α-1 Acid Glycoprotein reagents are intended for in vitro diagnostic use only.
- Sodium azide has been reported to form lead or copper azide in laboratory plumbing which may explode on percussion.
- Each donor unit used in the preparation of the standards and controls was found to be negative for the presence of HIV antibodies, as well as for Hepatitis B surface antigen, using a method approved by the FDA

WASTE MANAGEMENT

Please refer to local requirements.

REFERENCES

- Schmid, K. in FW Putman, Editor, The Plasma Proteins, Vol 1, second edition, Academic Press, New York, 1975, pp184-228
- Johnson, A.M. et al., J. Clin. Invest., 48(1969) 2293
- Dati, F. et al., Lab. Med. 13 (1989) 87

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