# Lipoprotein (a)

Diagnostic reagent for the quantitative in vitro determination of Lp (a) in human serum by turbidimetric assay

## REF

## Content

**A00527** 1x 10 mL Lp(a) Latex Reagent  
 5x 25 mL Lp(a) Buffer

Additionally offered:

A03736 4x 1 mL Lp(a) Calibrator 4 Level Series  
 A00733 1x 1 mL Lp(a) Calibrator High  
 A00734 1x 5 mL Lp(a) Calibrator High  
 A00731 1x 1 mL Lp(a) Calibrator Low  
 A00732 1x 5 mL Lp(a) Calibrator Low  
 A00816 1x 1 mL Lp(a) Control High  
 A00817 1x 5 mL Lp(a) Control High  
 A00814 1x 1 mL Lp(a) Control Low  
 A00815 1x 5 mL Lp(a) Control Low

## GENERAL INFORMATION

**Method** Immunoturbidimetric  
**Reaction** Nonlinear, endpoint  
**Wavelength** 600 nm  
**Assay Temperature** 18 – 37 °C  
**Sample** Serum  
**Measuring Range** approx. 0 – 80 mg/dL  
**Sensitivity** 6 mg/dL (Cobas Mira)  
**Hook Effect** without sample dilution: > 780 mg/dL (Cobas Mira)  
 with sample dilution: > 780 mg/dL (Cobas Mira)  
**Manual Test Procedure Tests/Kit\***  
 without sample dilution 125  
 with sample dilution 166

## Automated Test Procedure

Instrument dependent – please ask for applications

\* calculated on amount of antibody reagent; additional buffer on request

## REAGENT COMPOSITION

### COMPONENTS FINAL CONCENTRATION

**Lp(a) Latex Reagent**  
 Polystyrene latex particles sensitized with  
 Lp(a) antibodies variable  
 Sodium azide 0.095 %  
**Lp(a) Buffer**  
 Sodium Chloride 0.9 %  
 Sodium azide 0.095 %  
 Detergent

## REAGENT PREPARATION

The reagents are liquid and ready to use.

## REAGENT STABILITY AND STORAGE

**Conditions:** Protect from light. Close immediately after use.  
**Stability:** at 2 – 8 °C up to the expiration date  
 at 18 – 25 °C 1 month

Do not freeze!

## SAMPLE STABILITY AND STORAGE

**Stability:** at 2 – 8 °C 48 hours  
 at – 20 °C 3 months

Freeze only once!

## MANUAL TEST PROCEDURE

### Test Procedure without Sample Dilution:

Samples/Controls: ready to use

Calibration curve: Use Lp(a) Calibrator High to generate a calibration curve by making 1:2 serial dilutions of the calibrator with 0.9% saline as diluent or use the 4 level calibrator series. Use 0.9% saline as zero point.

Pipette into test tubes:	Calibrators	Samples/Controls
Buffer	900 µL	900 µL
Cal./Ctrls/Samples	3 µL	3 µL
Mix. Read A1 of calibrators and samples/controls at 600 nm. Then add:		
Latex Reagent	80 µL	80 µL
Mix. Incubate 5 minutes at assay temperature. Read A2 of calibrators and samples/controls at 600 nm. Calculate: $\Delta A = (A2 - A1)$		

### Test Procedure with Sample Dilution:

Sample/ Control: dilute 1:10 in saline 0.9%

Calibration curve: Use Lp(a) Calibrator High to generate a calibration curve by making 1:10, 1:20, 1:40, 1:80 and 1:160 dilutions with 0.9% saline as diluent. Use 0.9% saline as zero point.

Pipette into test tubes	Calibrators	Samples/Controls
Buffer	900 µL	900 µL
Cal./Ctrls/Samples	20 µL	20 µL
Mix. Read A1 of calibrators and samples/controls at 600 nm. Then add:		
Latex Reagent	60 µL	60 µL
Mix. Incubate 5 minutes at assay temperature. Read A2 of calibrators and samples/controls at 600 nm. Calculate: $\Delta A = (A2 - A1)$		

## CALCULATION

Calculate and plot  $\Delta A = (A2 - A1)$  of the calibrators versus assigned concentration values on a linear-linear graph paper. Calculate  $\Delta A$  optical densities of samples and control(s) and read values in mg/dL on the reference curve. Samples yielding absorbances above highest calibrator should be retested after further dilution.

## REFERENCE RANGE

0 – 30 mg/dL

It is recommended that each laboratory establishes its own normal range.

## TEST PRINCIPLE

The assay of Lp(a) is based on turbidimetric measurement. Turbidity is caused by the formation of antigen-antibody insoluble immuno complexes.

## DIAGNOSTIC IMPLICATIONS

Lp(a) is a human serum protein whose structure is close to that of LDL. Its density lies between those of LDL and HDL. The Lp(a) concentration in blood varies from almost undetectable levels to more than 100 mg/dL. The wide differences in Lp(a) levels are largely due to hereditary factors and cannot be controlled by dietary or lifestyle changes. The presence of high LP (a) levels in serum is a significant marker of increased risk for atherosclerosis and coronary heart disease. Epidemiological studies have shown, that people with normal serum cholesterol and a serum Lp(a) level over 30 mg/dL have a double risk of coronary heart disease. The risk is 8 times higher when LDL and Lp(a) levels are simultaneously elevated.

## PERFORMANCE CHARACTERISTICS

### SENSITIVITY

6 mg/dL (Cobas Mira)

### ACCURACY

Controls were assayed in duplicate on a Cobas Mira.

Control	Assigned Value (mg/dL)	Measured Value (mg/dL)
Dialab Control	40 (34 – 46)	39.2

## PRECISION

### Intra-Assay Precision

The test was performed on the Cobas Mira with a low, medium and high Lp(a) positive serum.

Expected Value	n	Mean	S.D.	C.V
Low	9	19.3	0.50	2.6
Medium	9	43.1	1.10	2.5
High	9	73.3	2.29	3.1

### Inter-Assay Precision

Two control sera with a low and high Lp(a) value were measured on the Cobas Mira during two weeks after calibration.

Expected Value	n	Mean	S.D.	C.V
Low	9	20	0.97	4.91
Medium	9	46	2.29	4.98

## METHOD COMPARISON

A comparison with Nephelometry gave the following results:

$$y = 0.910x - 1.973; r = 0.989$$

## INTERFERING SUBSTANCES

No interference up to:

Triglycerides	2500 mg/dL	Hemoglobin	500 mg/dL
Bilirubin	30 mg/dL	Plasminogen	200 mg/dL
APO B	200 mg/dL		

## QUALITY CONTROL

All Control sera with Lp(a) values measured by this method may be used. We recommend the Dialab Lp(a) Control High and the Lp(a) Control Low.

## CALIBRATION

The assay requires the use of Lp(a) serum Calibrators. We recommend the Dialab Lp(a) Calibrator 4 Level Series, the Lp(a) Calibrator High or the Lp(a) Calibrator Low.

## AUTOMATION

Applications for automated systems (with and without sample dilution) are available upon request.

## WARNINGS AND PRECATIONS

- The Lp(a) reagents are intended for in vitro diagnostic use only.
- Sodium azide has been reported to form lead or copper azide in laboratory plumbing which may explode on percussion.
- Each donor unit used in the preparation of the standards and controls was found to be negative for the presence of HIV antibodies, as well as for Hepatitis B surface antigen, using a method approved by the FDA

## WASTE MANAGEMENT

Please refer to local requirements.

## REFERENCES

- Poulik, M. D., and Weiss, M. L., in F. W. Putman, Editor, "The Plasma Proteins", vol. 2 second Edition, Academic Press, New York, pp 52 - 108

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