ASO Calibrator Super High

Liquid human based calibrator for use as a calibrator in immunoturbidimetric assays for ASO

REF	Cont.	
A00729	1x 1 mL	ASO Calibrator Super High
A00730	1x 5 mL	ASO Calibrator Super High

COMPOSITION

Dialab ASO Calibrator Super High is based on human serum, which contains high levels of ASO with phosphate buffered saline containing 1 g % bovine serum albumin.

Sodium azide 0.095 %

PREPARATION

The ASO Calibrator Super High is liquid, ready to use.

STABILITY AND STORAGE

Storage: at 2 – 8 °C
 Stability: up to the expiration date
 Stability after opening: 6 weeks at 2 – 8 °C
 DO NOT FREEZE!

DIRECTIONS FOR USE

- Use this Calibrator in the same manner as specified with the reagent and/or instrument being used.
- Allow the Calibrator to equilibrate to room temperature (22 – 28 °C) prior to use.
- The Calibrator should not be allowed to stand for repeated long periods of time (up to 6 hours) at room temperature. Keep vials tightly capped at all times when not in use to avoid microbial contamination. Keep vials refrigerated when not in use.

WARNINGS AND PRECAUTIONS

- For in Vitro Diagnostic use.
- Each individual donation intended for use in manufacture of protein calibrator was tested for hepatitis B surface antigen (HBsAg), anti-hepatitis C (anti-HCV) and anti-HIV 1 and HIV 2 by FDA required tests. Since no test method can assure that products derived from human blood do not contain HIV-1/2 and Hepatitis B and Hepatitis C virus, this material and all patient samples should be handled as though capable of transmitting infectious diseases.
- Reagents containing sodium azide must be handled with due caution: Do not ingest or allow to contact skin or mucous membranes! Sodium azide can form explosive azides when contacting heavy metals such as copper or lead.

WASTE MANAGEMENT

Please refer to local legal requirements.

LOT SPECIFIC ASSAY DATA

Values and expiry date are lot specific.

LOT: S2086		EXP: 2017/02	
Parameter	Assigned Value		
ASO *	395.0 IU/mL *		
ASO **	652.0 IU/mL **		
* traceable to WHO standardisation			
** traceable to Dade Behring standardisation			

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