## RF Control

Liquid human based control serum for use as a control in immunoturbidimetric assays for RF (Rheumatoid Factor)

REF	Cont.	
A00810	1x 1 mL	RF Control
A00811	1x 5 mL	RF Control

### COMPOSITION

Dialab RF Control is a liquid dilution of human plasma with high values of RF with PBS.

Sodium azide 0.095 %

### PREPARATION

The RF Control is liquid, ready to use.

### STABILITY AND STORAGE

Storage: at 2 – 8 °C  
 Stability: up to the expiration date  
 Stability after opening: 6 weeks at 2 – 8 °C  
**DO NOT FREEZE!**

### DIRECTIONS FOR USE

1. Use this control in the same manner as specified with the reagent and/or instrument being used.
2. Allow the control to equilibrate to room temperature (22 – 28 °C) prior to use.
3. The control should not be allowed to stand for repeated long periods of time (up to 6 hours) at room temperature. Keep vials tightly capped at all times when not in use to avoid microbial contamination. Keep vials refrigerated when not in use

### WARNINGS AND PRECAUTIONS

1. For in Vitro Diagnostic use.
2. Each individual donation intended for use in manufacture of protein calibrator was tested for hepatitis B surface antigen (HBsAg), anti-hepatitis C (anti-HCV) and anti-HIV 1 and HIV 2 by FDA required tests. Since no test method can assure that products derived from human blood do not contain HIV-1/2 and Hepatitis B and Hepatitis C virus, this material and all patient samples should be handled as though capable of transmitting infectious diseases.
3. Reagents containing sodium azide must be handled with due caution: Do not ingest or allow to contact skin or mucous membranes! Sodium azide can form explosive azides when contacting heavy metals such as copper or lead.

### WASTE MANAGEMENT

Please refer to local legal requirements.

### LOT SPECIFIC ASSAY DATA

Values and expiry date are lot specific.

LOT: C1043		EXP: 2017/04
Parameter	Target Value	Range
RF	105 IU/mL	89.0 – 121 IU/mL
Aggregated human IgG Method + Latex Method, based on WHO Standardisation		

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