## Fibronectin Calibrator

Lyophilised human based calibrator for use as a calibrator in immunoturbidimetric assays for Fibronectin

**REF**

**Cont.**

**A02712** 1x 0.5 mL Fibronectin Calibrator (Iyo.)

### COMPOSITION

Dialab Fibronectin Calibrator is a lyophilised preparation of stabilised, buffered human plasma containing less than 5 % bovine plasma.  
 Sodium azide 0.095 %

### PREPARATION

1. Reconstitute the Fibronectin Calibrator with 0.5 mL of distilled or deionized water.
2. Close vial and allow to stand until fully reconstituted (a minimum of 5 minutes).
3. Gently mix by rotation or inversion until dissolution is complete. DO NOT SHAKE.
4. Allow to equilibrate for a total of 30 minutes before use.

### STABILITY AND STORAGE

Storage: at 2 – 8 °C  
 Stability: up to the expiration date  
 Stability after reconstitution: 48 hours at 2 – 8 °C  
 DO NOT FREEZE!

### DIRECTIONS FOR USE

1. Use this Calibrator in the same manner as specified with the reagent and/or instrument being used.
2. Allow the Calibrator to equilibrate to room temperature (22 – 28 °C) prior to use.
3. The Calibrator should not be allowed to stand for repeated long periods of time (up to 6 hours) at room temperature. Keep vials tightly capped at all times when not in use to avoid microbial contamination. Keep vials refrigerated when not in use.

### WARNINGS AND PRECAUTIONS

1. For in Vitro Diagnostic use.
2. Each individual donation intended for use in manufacture of protein calibrator was tested for hepatitis B surface antigen (HBsAg), anti-hepatitis C (anti-HCV) and anti-HIV 1 and HIV 2 by FDA required tests. Since no test method can assure that products derived from human blood do not contain HIV-1/2 and Hepatitis B and Hepatitis C virus, this material and all patient samples should be handled as though capable of transmitting infectious diseases.
3. Reagents containing sodium azide must be handled with due caution: Do not ingest or allow to contact skin or mucous membranes! Sodium azide can form explosive azides when contacting heavy metals such as copper or lead.

### WASTE MANAGEMENT

Please refer to local legal requirements.

### LOT SPECIFIC ASSAY DATA

Values and expiry date are lot specific.

<b>LOT: 26715</b>	<b>EXP: 2016/07</b>
<b>Parameter</b>	<b>Assigned Value</b>
Fibronectin	62.00 mg/dL

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