Pediatric Control

Liquid human based control serum for use as a control in immunoturbidimetric assays

REF	Cont.	
A03808	1x 1 mL	Pediatric Control
A03809	1x 5 mL	Pediatric Control

COMPOSITION

Dilution of delipitated and defibrinated normal human plasma with PBS.
 Sodium azide 0.095 %

PREPARATION

The Control is liquid, ready to use.

STABILITY AND STORAGE

Storage: at 2 – 8 °C
 Stability: up to the expiration date
 Stability after opening: 6 weeks at 2 – 8 °C
DO NOT FREEZE!

DIRECTIONS FOR USE

- Use this Calibrator in the same manner as specified with the reagent and/or instrument being used.
- Allow the Calibrator to equilibrate to room temperature (22 – 28 °C) prior to use.
- The Calibrator should not be allowed to stand for repeated long periods of time (up to 6 hours) at room temperature. Keep vials tightly capped at all times when not in use to avoid microbial contamination. Keep vials refrigerated when not in use

WARNINGS AND PRECATIONS

- For in Vitro Diagnostic use.
- Each individual donation intended for use in manufacture of protein calibrator was tested for hepatitis B surface antigen (HBsAg), anti-hepatitis C (anti-HCV) and anti-HIV 1 and HIV 2 by FDA required tests. Since no test method can assure that products derived from human blood do not contain HIV-1/2 and Hepatitis B and Hepatitis C virus, this material and all patient samples should be handled as though capable of transmitting infectious diseases.
- Reagents containing sodium azide must be handled with due caution: Do not ingest or allow to contact skin or mucous membranes! Sodium azide can form explosive azides when contacting heavy metals such as copper or lead.

WASTE MANAGEMENT

Please refer to local legal requirements.

LOT SPECIFIC ASSAY DATA

Values and expiry date are lot specific.

LOT: PED009		EXP: 2016/01
Parameter	Target Value	Range
α-2 Macroglobulin	10.3 mg/dL	8.8 – 11.8 mg/dL
Transferrin	13.0 mg/dL	11.1 – 15.0 mg/dL
IgA	9.9 mg/dL	8.4 – 11.4 mg/dL
IgG	44.9 mg/dL	38.2 – 51.6 mg/dL
IgM	6.0 mg/dL	5.1 – 6.9 mg/dL
Kappa Light Chain **	12.0 mg/dL	10.2 – 13.8 mg/dL
Lambda Light Chain **	5.8 mg/dL	4.9 – 6.7 mg/dL

Values traceable to ERM-DA470k/IFCC.
 ** Values based on Dade Behring Material.

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