## RF Calibrator 5 Level Series

Liquid human based calibrators for use as calibrators in immunoturbidimetric assays for RF

**REF**

A11720

**Cont.**

5x 1 mL RF Calibrator 5 Level Series

### COMPOSITION

Dilution of human plasma containing a high level of RF with saline. The dilutions are liquid stabilized.

Sodium azide 0.095 %

### PREPARATION

The RF Calibrator Series is liquid, ready to use.

### STABILITY AND STORAGE

Storage: at 2 – 8 °C  
 Stability: up to the expiration date  
 Stability after opening: 6 weeks at 2 – 8 °C  
 DO NOT FREEZE!

### DIRECTIONS FOR USE

- Use these Calibrators in the same manner as specified with the reagent and/or instrument being used.
- Allow the Calibrators to equilibrate to room temperature (22 – 28 °C) prior to use.
- The Calibrators should not be allowed to stand for repeated long periods of time (up to 6 hours) at room temperature. Keep vials tightly capped at all times when not in use to avoid microbial contamination. Keep vials refrigerated when not in use

### WARNINGS AND PRECAUTIONS

- For in Vitro Diagnostic use.
- Each individual donation intended for use in manufacture of protein calibrator was tested for hepatitis B surface antigen (HBsAg), anti-hepatitis C (anti-HCV) and anti-HIV 1 and HIV 2 by FDA required tests. Since no test method can assure that products derived from human blood do not contain HIV-1/2 and Hepatitis B and Hepatitis C virus, this material and all patient samples should be handled as though capable of transmitting infectious diseases.
- Reagents containing sodium azide must be handled with due caution: Do not ingest or allow to contact skin or mucous membranes! Sodium azide can form explosive azides when contacting heavy metals such as copper or lead.

### WASTE MANAGEMENT

Please refer to local legal requirements.

### LOT SPECIFIC ASSAY DATA

Values and expiry date are lot specific.

LOT: S1085		EXP: 2015/05
Parameter	Level	Assigned Value
RF	Level 1	35.00 IU/mL
RF	Level 2	67.00 IU/mL
RF	Level 3	128.0 IU/mL
RF	Level 4	253.0 IU/mL
RF	Level 5	502.0 IU/mL
Values are traceable to WHO standardisation.		

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