

A rapid test for the qualitative detection of Benzodiazepines in human urine. For medical and other professional *in vitro* diagnostic use only.

【INTENDED USE】

The DIA Rapid Test Dipstick (Urine) is a rapid chromatographic immunoassay for the detection of Diazepam in urine at a cut-off concentration of 300ng/ml. This test will detect other related compounds, please refer to the Analytical Specificity table in this package insert.

This assay provides only a qualitative, preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are used.

【SUMMARY】

Benzodiazepines are medications that are frequently prescribed for the symptomatic treatment of anxiety and sleep disorders. They produce their effects via specific receptors involving a neurochemical called gamma aminobutyric acid (GABA). Because they are safer and more effective, Benzodiazepines have replaced Barbiturates in the treatment of both anxiety and insomnia. Benzodiazepines are also used as sedatives before some surgical and medical procedures, and for the treatment of seizure disorders and alcohol withdrawal.

Risk of physical dependence increases if Benzodiazepines are taken regularly (e.g., daily) for more than a few months, especially at higher than normal doses. Stopping abruptly can bring on such symptoms as trouble sleeping, gastrointestinal upset, feeling unwell, loss of appetite, sweating, trembling, weakness, anxiety and changes in perception. Only trace amounts (less than 1%) of most Benzodiazepines are excreted unaltered in the urine; most of the concentration in urine is conjugated drug. The detection period for the Benzodiazepines in the urine is 3-7 days.

The DIA Rapid Test Dipstick (Urine) is a rapid urine-screening test that can be performed without the use of an instrument. The test utilizes the antibody to selectively detect elevated levels of Benzodiazepines in urine. The DIA Test Dipstick (Urine) yields a positive result when the Benzodiazepines in urine exceeds the cut-off level.

【PRINCIPLE】

The DIA Rapid Test Dipstick (Urine) is an immunoassay based on the principle of competitive binding. Drugs which may be present in the urine specimen compete against the drug conjugate for binding sites on the antibody.

During testing, a urine specimen migrates upward by capillary action. Benzodiazepines, if present in the urine specimen below the cut-off level, will not saturate the binding sites of the antibody in the test. The antibody coated particles will then be captured by immobilized Benzodiazepines-protein conjugate and a visible colored line will show up in the test line region. The colored line will not form in the test line region if the Benzodiazepines level exceeds the cut-off level, because it will saturate all the binding sites of anti-Benzodiazepines antibody.

A drug-positive urine specimen will not generate a colored line in the test line region because of drug competition, while a drug-negative urine specimen or a specimen containing a drug concentration less than the cut-off will generate a line in the test line region. To serve as a procedural control, a colored line will always appear at the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred.

【REAGENTS】

The test contains mouse monoclonal anti-Benzodiazepines antibody coupled particles and Benzodiazepines-protein conjugate. A goat antibody is employed in the control line system.

【PRECAUTIONS】

- For medical and other professional *in vitro* diagnostic use only. Do not use after the expiration date.
- The test should remain in the sealed pouch until use.
- All specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The used test should be discarded according to local regulations.

【STORAGE AND STABILITY】

Store as packaged at room temperature or refrigerated (2-30°C). The test is stable through the expiration date printed on the sealed pouch or label of the closed canister. The test must remain in the sealed pouch or closed canister until use. **DO NOT FREEZE.** Do not use beyond the expiration date. NOTE: Once the canister has been opened, the remaining test(s) are stable for 50 days only.

【SPECIMEN COLLECTION AND PREPARATION】

Urine Assay

The urine specimen must be collected in a clean and dry container. Urine collected at any time of the day may be used. Urine specimens exhibiting visible particles should be centrifuged, filtered, or allowed to settle to obtain a clear specimen for testing.

Specimen Storage

Urine specimens may be stored at 2-8°C for up to 48 hours prior to testing. For prolonged storage, specimens may be frozen and stored below -20°C. Frozen specimens should be thawed and mixed before testing.


【MATERIALS】

- Test Dipsticks
 - Package insert
- Specimen collection container
 - Materials Required **But Not Provided**
 - Timer

【DIRECTIONS FOR USE】

Allow the test, urine specimen, and/or controls to reach room temperature (15-30°C) prior to testing.

- Bring the pouch to room temperature before opening it. Remove the test dipstick from the sealed pouch and use it within one hour.
- With arrows pointing toward the urine specimen, **immerse the test dipstick vertically in the urine specimen for at least 10-15 seconds.** Do not pass the maximum line (MAX) on the Test Dipstick when immersing the strip. See the illustration below.
- Place the test dipstick on a non-absorbent flat surface, start the timer and wait for the colored line(s) to appear. **Read results at 5 minutes.** Do not interpret the result after 10 minutes



【INTERPRETATION OF RESULTS】

(Please refer to the illustration above)

NEGATIVE: * Two lines appear. One colored line should be in the control line region (C), and another apparent colored line should be in the test line region (T). This negative result indicates that the Benzodiazepine concentration is below the detectable cut-off level.

***NOTE:** The shade of color in the test line region (T) may vary, but it should be considered negative whenever there is even a faint colored line.

POSITIVE: One colored line appears in the control line region (C). No line appears in the test line region (T). This positive result indicates that the Benzodiazepine concentration exceeds the detectable cut-off level.

INVALID: Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

【QUALITY CONTROL】

A procedural control is included in the test. A colored line appearing in the control line region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

Control standards are not supplied with this kit; however, it is recommended that positive and negative controls be tested as good laboratory testing practice to confirm the test procedure and to verify proper test performance.

【LIMITATIONS】

- The DIA Rapid Test Dipstick (Urine) provides only a qualitative, preliminary analytical result. A secondary analytical method must be used to obtain a confirmed result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method.^{1,2}
- It is possible that technical or procedural errors, as well as other interfering substances in the urine specimen may cause erroneous results.
- Adulterants, such as bleach and/or alum, in urine specimens may produce erroneous results regardless of the analytical method used. If adulteration is suspected, the test should be repeated with another urine specimen.
- A positive result indicates presence of the drug or its metabolites but does not indicate level of intoxication, administration route or concentration in urine.
- A negative result may not necessarily indicate drug-free urine. Negative results can be obtained when drug is present but below the cut-off level of the test.
- Test does not distinguish between drugs of abuse and certain medications.

【EXPECTED VALUES】

This negative result indicates that the Benzodiazepines concentration is below the detectable level of 300ng/ml. Positive result means the concentration of Benzodiazepines is above the level of 300ng/ml. The DIA Rapid Test Dipstick has a sensitivity of 300ng/ml

【PERFORMANCE CHARACTERISTICS】

Accuracy

A side-by-side comparison was conducted using The DIA Rapid Test Dipstick (Urine) and GC/MS at the cut-off of 300ng/ml. Testing was performed on 250 clinical specimens previously collected from subjects present for Drug Screen Testing. The following results were tabulated:

Method	GC/MS		Total Results
	Positive	Negative	
DIA Rapid Test Dipstick	121	1	122
	2	126	128
Total Results	123	127	250
% Agreement	98.4%	99.2%	98.8%

Analytical Sensitivity

A drug-free urine pool was spiked with Diazepam at the following concentrations: 0ng/ml, 150ng/ml, 225 ng/ml, 300ng/ml, 375ng/ml, 450ng/ml and 900 ng/ml. The result demonstrates >99% accuracy at 50% above and 50% below the cut-off concentration. The data are summarized below:

Diazepam Concentration (ng/mL)	Percent of Cut-off	n	Visual Result	
			Negative	Positive
0	0	30	30	0
150	50%	30	30	0
225	75%	30	27	3
300	Cut-off	30	15	15
375	+25%	30	3	27
450	+50%	30	0	30
900	3X	30	0	30

Analytical Specificity

The following table lists compounds that are positively detected in urine by the DIA Rapid Test Dipstick (Urine) at 5 minutes.

Compound	Concentration (ng/mL)	Compound	Concentration (ng/mL)
Alprazolam	100	Flunitrazepam	200
a-hydroxyalprazolam	1,500	(±) Lorazepam	3,000
Bromazepam	900	RS-Lorazepam glucuronide	200
Chlordiazepoxide	900	Midazolam	6,000
Clobazam	200	Nitrazepam	200
Clonazepam	500	Norchlordiazepoxide	100
Clorazepate	900	Nordiazepam	900
Delorazepam	900	Oxazepam	300
Desalkylflurazepam	200	Temazepam	100
Diazepam	300	Triazolam	3,000
Estazolam	6,000		

Precision

A study was conducted at three hospitals by laypersons using three different lots of product to demonstrate the within run, between run and between operator precision. An identical panel of coded specimens containing, according to GC/MS, no Diazepam, 25% Diazepam above and below the cut-off and 50% Diazepam above and below the 300ng/ml cut-off was provided to each site. The following results were tabulated:

Diazepam Concentration (ng/mL)	n per Site	Site A		Site B		Site C	
		-	+	-	+	-	+
0	10	10	0	10	0	10	0
150	10	10	0	10	0	10	0
225	10	9	1	9	1	9	1
375	10	1	9	1	9	1	9
450	10	0	10	0	10	0	10

Effect of Urinary Specific Gravity

Fifteen urine specimens of normal, high, and low specific gravity ranges were spiked with 150ng/ml and 450ng/ml of Diazepam. The DIA Test Dipstick (Urine) was tested in duplicate using the fifteen neat and spiked urine specimens. The results demonstrate that varying ranges of urinary specific gravity do not affect the test results.

Effect of Urinary pH

The pH of an aliquoted negative urine pool was adjusted to a pH range of 5 to 9 in 1 pH unit increments and spiked with Diazepam to 150ng/ml and 450ng/ml. The spiked, pH-adjusted urine was tested with the DIA Rapid Test Dipstick (Urine) in duplicate. The results demonstrate that varying ranges of pH do not interfere with the performance of the test.

Cross-Reactivity

A study was conducted to determine the cross-reactivity of the test with compounds in either drug-free urine or Diazepam positive urine. The following compounds show no cross-reactivity when tested with the DIA Rapid Test Dipstick (Urine) at a concentration of 100 µg/ml.

Non Cross-Reacting Compounds

Acetaminophen	Deoxycorticosterone	MDE	β-Phenylethylamine
Acetophenetidin	Dextromethorphan	Mepiperidine	Phenylpropanolamine
N-Acetylprocainamide	Diclofenac	Meprobamate	Prednisolone
Acetylsalicylic acid	Diffunisal	Methadone	Prednisone
Aminopyrine	Digoxin	L-Methamphetamine	Procaine
Amitypyline	Diphenhydramine	Morphoxyphenamine	Promazine
Amobarbital	Doxylamine	(±) - 3,4-Methylenedioxy-amphetamine	Promethazine
Amoxicillin	Ecgonine	(±) - 3,4-Methylenedioxy-methamphetamine	D,L-Proparalolol
Ampicillin	Ecgonine methyl ester	Morphine-3-β-D glucuronide	D-Propoxyphene
L-Ascorbic acid	(-)-ψ-Ephedrine	Morphine Sulfate	D-Pseudoephedrine
D,L-Amphetamine sulfate	{1R,2S} (-) Ephedrine	Nalidixic acid	Quinacrine
Apomorphine	(L) - Epinephrine	Naloxone	Quinidine
Aspartame	Erythromycin	Naltrexone	Quinine
Atropine	β-Estradiol	Naproxen	Ranitidine
Benzilic acid	Estrone-3-sulfate	Niacinamide	Salicylic acid
Benzoinic acid	Ethyl-p-aminobenzoate	Nifedipine	Secobarbital
Benzoylcegonine	Fenoprofen	Nirocodeine	Serotonin
Benzphetamine	Furosemide	Norethindrone	Sulfamethazine
Bilirubin	Genistic acid	Hydrochlorothiazide	Sulindac
(±) - Brompheniramine	Hemoglobin	Hydrocodone	Tetracycline
Caffeine	Hydralazine	Hydrocortisone	Tetrahydrocortisone, 9-Acetate
Cannabinoid	Hydrochlorothiazide	Hydrocodone	Tetrahydrocortisone
Cannabidiol	Hydrocodone	Hydrocortisone	3(β-D-glucuronide)
Carlorhydrate	Hydrocortisone	O-Hydroxyhippuric acid	Tetrahydrozoline
Chloramphenicol	O-Hydroxyhippuric acid	Oxycodone	Thiamine
Chlorothiazide	p-Hydroxyamphetamine	p-Hydroxy-metamphetamine	Thioridazine
(±) - Chlorpheniramine	p-Hydroxy-metamphetamine	Chlorquine	D,L-Tyrosine
Chlorpromazine	Chlorquine	3-Hydroxytyramine	Tolbutamide
Cholesterol	lbutrofen	Clomipramine	Triamterene
Clomipramine	Imipramine	Clonidine	Trifluoperazine
Clonidine	Iproniazid	Cocacethylene	Perphenazine
Cocacethylene	(±) - Isoproterenol	Cocaine	Phencyclidine
Codine	Isoxsuprine	Codine	Phenelzine
Cortisone	Ketoprofen	Cortisone	Phenobarbital
(-) Cotinine	Labetalol	Creatinine	Phentermine
Creatinine	Loperamide		Trans-2-phenylcyclopropylamine hydrochloride
	Maprotiline		Verapamil
			Zomepirac

【BIBLIOGRAPHY】

- Babel RC. *Disposition of Toxic Drugs and Chemicals in Man.* 2nd Ed. Biomedical Publ., Davis, CA. 1992; 488
- Hawks RL, CN Chiang. *Urine Testing for Drugs of Abuse.* National Institute for Drug Abuse (NIDA), Research Monograph 73, 1986

Attention, see instructions for use	For in vitro diagnostic use only	Store between 2-30°C	Do not use if package is damaged	Index of Symbols	EC REP	Authorized Representative
					Use by	Do not reuse
		Lot Number			REF	Catalog #