Instructions for Use

TNF- α ELISA

IVD

REF EIA-4641

Σ 96

DRG

DRG Instruments GmbH, Germany
Frauenbergstraße. 18, D-35039 Marburg
Phone: +49 (0)6421-1700 0, Fax: +49 (0)6421-1700 50
Website: www.drg-diagnostics.de
E-mail: drg@drg-diagnostics.de

Distributed by:

DRG

DRG International, Inc., USA
841 Mountain Ave., Springfield, NJ 07081
Phone: (973) 564-7555, Fax: (973) 564-7556
Website: www.drg-international.com
E-mail: corp@drg-international.com

**Please use only the valid version of the Instructions for Use provided with the kit.
Verwenden Sie nur die jeweils gültige, im Testkit enthaltene, Gebrauchsanweisung.
Si prega di usare la versione valida delle istruzioni per l'uso a disposizione con il kit.
Por favor, se usa solo la version valida de la metodico técnico incluido aqui en el kit.
Utilisez seulement la version valide des Instructions d'utilisation fournies avec le kit.**

Table of Contents / Inhaltsverzeichnis / Tabella die Contenuti / Tabla de Contenidos

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1 INTENDED USE

Immunoenzymetric assay for the in vitro quantitative measurement of human Tumor Necrosis Factor α (TNF- α) in serum.

2 CLINICAL BACKGROUND

2.1 Biological activities

Human Tumor Necrosis Factor Alpha (TNF- α) also named cachectin, is a 157 A.A. unglycosylated polypeptide cytokine mainly produced by activated macrophages (monocytes). Lipopolysaccharide (LPS), the cell-wall component of gram-negative bacteria (endotoxin), is a potent stimulus for TNF- α production by macrophages and TNF- α is an important mediator of the well-known in vivo effects of LPS such as tumour hemorrhagic necrosis, fever, shock and activation of neutrophils. The various biological activities of TNF- α may be classified as:

- *Antitumoral and growth regulatory activities*: TNF- α displays a selective toxicity for tumor and virus-infected cells. Conversely, it is angiogenic and stimulates the growth of cultured fibroblasts.
- *Immunomodulatory and proinflammatory activities*: TNF- α activates macrophages, neutrophils and eosinophils, as well as endothelial cells (which display procoagulant activity). It regulates the production of antibodies by B cells and stimulates cytotoxic T cells. It induces the production of several other inflammatory mediators such as IL-1, IL-6, colony stimulating factors, prostaglandins, platelet-activating factor (PAF), collagenases, etc.
- *Metabolic activities* : TNF- α strongly inhibits lipoprotein lipase and adipocyte gene expression.

2.2 Clinical application

TNF- α has a major pathogenic role: in cachexia associated with chronic infectious or cancerous diseases ; in septic shock where the neutralization of TNF- α protects against the associated acute lethality ; in graft rejection and graft-versus-host disease ; and in parasitic infections where TNF- α may provide some protection but also favours more severe forms of the disease (e.g. the cerebral form of malaria). TNF- α often in combination with other cytokines, has also been involved in several autoimmune diseases and even in the pathogenesis of arteriosclerosis. Abnormal high levels of serum TNF- α have been described in septic shock, graft rejection, parasitic infections, cancer, post hemofiltrations, during in vivo cytokine (IL-2) therapy, etc. Besides an insight into pathogenesis, these determinations might provide an aid in diagnosis (e.g. in graft rejection) and have prognostic value (e.g. in systemic infections).

3 PRINCIPLES OF THE METHOD

The DRG TNF- α -ELISA is a solid phase Enzyme Amplified Sensitivity Immunoassay performed on microtiter plate. The assay uses monoclonal antibodies (MAbs) directed against distinct epitopes of TNF- α . Standards and samples react with the capture monoclonal antibody (MAb 1) coated on microtiter well and with a monoclonal antibody (MAb 2) labelled with horseradish peroxidase (HRP). After an incubation period allowing the formation of a sandwich: coated MAb 1 – human TNF- α – MAb 2 – HRP, the microtiter plate is washed to remove unbound enzyme labelled antibody. Bound enzyme-labelled antibody is measured through a chromogenic reaction. Chromogenic solution (TMB) is added and incubated. The reaction is stopped with the addition of Stop Solution and the microtiter plate is then read at the appropriate wavelength. The amount of substrate turnover is determined colourimetrically by measuring the absorbance, which is proportional to the TNF- α concentration.

A calibration curve is plotted and TNF- α concentration in samples is determined by interpolation from the calibration curve. The use of the ELISA reader (linearity up to 3 OD units) and a sophisticated data reduction method (polychromatic data reduction) result in a high sensitivity in the low range and in an extended calibration range.

4 REAGENTS PROVIDED

	Reagents	96 tests Kit	Reconstitution
MICROTITERPLATE	Microtiter plate with 96 anti-TNF- α (monoclonal antibodies) coated wells	96 wells	Ready for use
Ab HRP CONC	Conjugate: HRP labelled anti-TNF- α (monoclonal antibodies) in TRIS-Maleate buffer with bovine serum albumin and thymol	1 vial 0.75 mL	Add conjugate buffer (see section 6)
CAL 0	Zero calibrator in human plasma, benzamidin and thymol	2 vials lyophil.	Add distilled water (see on the label for the exact volume)
CAL N	Calibrator N = 1 to 5 (see exact values on vial labels) in human plasma, benzamidine and thymol	5 vials lyophil.	Add 2 mL distilled water
CON BUF	Conjugate buffer: TRIS-Maleate buffer with bovine serum albumin, EDTA and thymol	1 vial 6 mL	Ready for use
INC BUF	Incubation buffer: TRIS-Maleate buffer with bovine serum albumin, EDTA and thymol	1 vial 6 mL	Ready for use
WASH SOLN CONC	Wash Solution Conc. (Tris-HCl)	1 vial 10 mL	Dilute 200 x with distilled water (use a magnetic stirrer).
CONTROL N	Control - N = 1 or 2 in human plasma and thymol	2 vials lyophil.	Add 2 mL distilled water
CHROM TMB	Chromogen TMB (Tetramethylbenzidine)	1 vial 12 mL	Ready for use
STOP SOLN	Stopping solution: HCl 1.0N	1 vial 12 mL	Ready for use

- Note:**
1. Use the zero calibrator for sample dilutions.
 2. 1 pg of the calibrator preparation is equivalent to 40 mIU of the NIBSC IS 87/650.

5 SUPPLIES NOT PROVIDED

The following material is required but not provided in the kit:

1. High quality distilled water
2. Pipettes for delivery of: 50 μ L, 200 μ L, 1 mL and 10 mL (the use of accurate pipettes with disposable plastic tips is recommended)
3. Vortex mixer
4. Magnetic stirrer
5. Horizontal microtiter plate shaker capable of 700 rpm \pm 100 rpm
6. Washer for microtiter plates
7. Microtiter plate reader capable of reading at 450 nm, 490 nm and 650 nm (in case of polychromatic reading) or capable of reading at 450 nm and 650 nm (bichromatic reading)

6 REAGENT PREPARATION

A. Calibrators:

Reconstitute the zero calibrator to the volume specified on the vial label with distilled water and the other calibrators with 2 mL distilled water.

B. Controls:

Reconstitute the controls with 2 mL distilled water.

C. Conjugate Solution:

following the number of wells to be used, dilute the concentrated conjugate with the conjugate buffer in a clean glass vial: see below table for the volumes to pipette. Extemporaneous preparation is recommended. Diluted conjugate is stable for max. 1 week at 2 °C - 8 °C.

TABLE CONJUGATE DILUTION

Number of wells	Concentrated conjugate	Conjugate buffer	Working volume
8	50 μ L	500 μ L	550 μ L
16	100 μ L	1000 μ L	1100 μ L
24	150 μ L	1500 μ L	1650 μ L
32	200 μ L	2000 μ L	2200 μ L
48	300 μ L	3000 μ L	3300 μ L
96	600 μ L	6000 μ L	6600 μ L

D. Working Wash solution:

Prepare an adequate volume of Working Wash solution by adding 199 volumes of distilled water to 1 volume of Wash Solution (200x). Use a magnetic stirrer to homogenize. Discard unused Working Wash solution at the end of the day.

7 STORAGE AND EXPIRATION DATING OF REAGENTS

- Before opening or reconstitution, all kits components are stable until the expiry date, indicated on the vial label, if kept at 2 °C - 8°C.
- Unused strips must be stored, at 2 °C - 8 °C, in a sealed bag containing a desiccant until expiration date.
- After reconstitution, standards and controls are stable for 4 days at 2 °C - 8 °C. For longer storage periods, aliquots should be made and kept at -20 °C for maximum 2 months. Avoid successive freeze thaw cycles.
- The concentrated Wash Solution is stable at 18 °C - 25 °C until expiration date.
- Freshly prepared Working Wash solution should be used on the same day.
- After its first use, the conjugate is stable until expiry date, if kept in the original well-closed vial at 2 °C - 8 °C.
- Alterations in physical appearance of kit reagents may indicate instability or deterioration.

8 SPECIMEN COLLECTION AND PREPARATION

- Serum must be removed as soon as possible from the clot of red cells after clotting and centrifugation, and kept at 4 °C. If the samples are not used immediately, they must be kept at -20 °C for maximum 2 months, and at -70 °C for longer storage (maximum one year).
- Avoid subsequent freeze thaw cycles.
- Prior to use, all samples should be at 18 °C - 25 °C. It is recommended to vortex the samples before use.
- Sampling conditions can affect values, therefore, strict precautions have to be taken during sampling to avoid impurities contained in sampling materials that would stimulate TNF- α production by blood cells and thus falsely increase serum TNF- α values.
- Collection tubes must be pyrogen-free.

9 PROCEDURE

9.1 Handling notes

Do not use the kit or components beyond expiry date.

Do not mix materials from different kit lots.

Bring all the reagents 18 °C - 25 °C prior to use.

Thoroughly mix all reagents and samples by gentle agitation or swirling.

Perform calibrators, controls and samples in duplicate. Vertical alignment is recommended.

Use a clean plastic container to prepare the Wash Solution.

In order to avoid cross-contamination, use a clean disposable pipette tip for the addition of each reagent and sample.

For the dispensing of the Revelation Solution and the Stop Solution avoid pipettes with metal parts.

High precision pipettes or automated pipetting equipment will improve the precision.

Respect the incubation times.

To avoid drift, the time between pipetting of the first standard and the last sample must be limited to the time mentioned in section *Time delay*.

Prepare a calibration curve for each run, do not use data from previous runs.

The Chromogen TMB Solution should be colourless. If a blue colour develops within a few minutes after preparation, this indicates that the reagent is unusable, and must be discarded.

Dispense the Chromogen TMB Solution within 15 minutes following the washing of the microtiter plate.

During incubation with Revelation Solution, avoid direct sunlight on the microtiter plate.

9.2 Procedure

1. Select the required number of strips for the run. The unused strips should be resealed in the bag with a desiccant and stored at 2 °C - 8 °C.
2. Secure the strips into the holding frame.
3. Pipette 50 μ L of incubation buffer into all the wells
4. Pipette 200 μ L of each Calibrator, Control and Sample into the appropriate wells.
5. Incubate for 2 hours at 18 °C - 25 °C on a horizontal shaker set at 700 rpm \pm 100 rpm.
6. Aspirate the liquid from each well.
7. Wash the plate 3 times by:
Dispensing 0.4 mL of Wash Solution into each well
Aspirating the content of each well
8. Pipette 100 μ L of zero calibrator into all the wells
9. Pipette 50 μ L of anti- TNF- α -HRP conjugate into all the wells.
10. Incubate for 2 hours at 18 °C - 25 °C on a horizontal shaker set at 700 rpm \pm 100 rpm.
11. Aspirate the liquid from each well.
12. Wash the plate 3 times by:
Dispensing 0.4 mL of Wash Solution into each well
Aspirating the content of each well
13. Pipette 100 μ L of Chromogen TMB into each well within 15 minutes following the washing step.
14. Incubate the microtiter plate for 15 minutes at 18 °C - 25 °C on a horizontal shaker set at 700 rpm \pm 100 rpm, avoid direct sunlight.
15. Pipette 100 μ L of Stop solution into each well.
16. Read the absorbencies at 450 nm and 490 nm (reference filter 630 nm or 650 nm) within 30 minutes 3 hours and calculate the results as described in section 10.

10 CALCULATION OF RESULTS

10.1 Polychromatic Reading

1. In this case, the software will do the data processing.
2. The plate is first read at 450 nm against a reference filter set at 650 nm (or 630 nm).
3. A second reading is performed at 490 nm against the same reference filter.
4. The Software will drive the reader automatically and will integrate both readings into a polychromatic model. This technique can generate OD's up to 10.
5. The principle of polychromatic data processing is as follows:

X_i = OD at 450 nm

Y_i = OD at 490 nm

Using a standard unweighted linear regression, the parameters A & B are calculated: $Y = A \times X + B$

If $X_i < 3$ OD units, then X calculated = X_i

If $X_i > 3$ OD units, then X calculated = $(Y_i - B) / A$

A 4 parameter logistic curve fitting is used to build up the calibration curve.

The TNF- α concentration in samples is determined by interpolation on the calibration curve.

10.2 Bichromatic Reading

1. Read the plate at 450 nm against a reference filter set at 650 nm (or 630 nm).
2. Calculate the mean of duplicate determinations.
3. On semi-logarithmic or linear graph paper plot the OD values (ordinate) for each standard against the corresponding concentration of TNF- α (abscissa) and draw a calibration curve through the standard points by connecting the plotted points with straight lines.
4. Read the concentration for each control and sample by interpolation on the calibration curve.
5. Computer assisted data reduction will simplify these calculations. If automatic result processing is used, a 4 parameter logistic function curve fitting is recommended.

11 TYPICAL DATA

The following data are for illustration only and should never be used instead of the real time calibration curve.

TNF- α -ELISA		OD units Polychromatic model
Standard	0 pg/mL	0.045
	6.8 pg/mL	0.120
	18 pg/mL	0.259
	52 pg/mL	0.619
	176 pg/mL	1.435
	518 pg/mL	3.237

12 PERFORMANCE AND LIMITATIONS

12.1 Detection Limit

Twenty zero standards were assayed along with a set of other standards. The detection limit, defined as the apparent concentration two standard deviations above the average OD at zero binding, was 0.7 pg/mL.

12.2 Specificity

No significant cross-reaction was observed in presence of 50 ng of IL-1 α , IL-1 β , IL-2, IL-3, IL-4, IL-6, IL-7, IL-8, IL-10, TNF- β , IFN- α , IFN- β , IFN- γ , TGF- β , GM-CSF, OSM, MIP-1 α , MIP-1 β , LIF, MCP-1, G-CSF and RANTES.

This TNF- α assay is specific for human natural and recombinant TNF- α .

12.3 Precision

INTRA ASSAY				INTER ASSAY			
Serum	N	<X> \pm SD (pg/mL)	CV (%)	Serum	N	<X> \pm SD (pg/mL)	CV (%)
A	20	91 \pm 6	6.6	A	24	122 \pm 5	4.5
B	20	526 \pm 33	6.3	B	24	431 \pm 14	3.3

SD : Standard Deviation; CV: Coefficient of variation

12.4 Accuracy**RECOVERY TEST**

Sample	Added TNF- α (pg/mL)	Recovered TNF- α (pg/mL)	Recovery (%)
Serum 1	0	6.2	-
	38.4	43.3	97
	83.9	90.0	100
	188.3	192.5	99
	408.2	376.2	91
Serum 2	0	3.8	-
	38.4	45.5	108
	83.9	91.2	104
	188.3	162.2	84
	408.2	379.2	92

DILUTION TEST

Sample	Dilution	Theoretical Concent. (pg/mL)	Measured Concent. (pg/mL)
Serum 1	1	-	436.5
	2	218.3	212.4
	4	109.1	104.8
	8	54.6	59.5
	16	27.3	31.7
Serum 2	1	-	420.2
	2	210.1	211.2
	4	105.0	98
	8	52.5	58.3
	16	26.3	30.7

Samples were diluted with zero standard.

12.5 Time delay between last standard and sample dispensing

As shown hereafter, assay results remain accurate even when a sample is dispensed 30 minutes after the standards have been added to the coated wells.

	T0	30 min	45 min
SC1	202	183	222
SC2	506	520	565

13 INTERNAL QUALITY CONTROL

If the results obtained for Control 1 and/or Control 2 are not within the range specified on the vial label, the results cannot be used unless a satisfactory explanation for the discrepancy has been given.

If desirable, each laboratory can make its own pools of control samples, which should be kept frozen in aliquots. Controls that contain azide will interfere with the enzymatic reaction and cannot be used.

Acceptance criteria for the difference between the duplicate results of the samples should rely on Good Laboratory Practises

It is recommended that controls be routinely assayed as unknown samples to measure assay variability. The performance of the assay should be monitored with quality control charts of the controls.

It is good practise to check visually the curve fit selected by the computer.

14 REFERENCE INTERVALS

These values are given only for guidance; each laboratory should establish its own normal range of values.

For guidance, the results of 30 serum samples from apparently healthy persons with low CRP levels, ranged between 4.6 and 12.4 pg/mL.

15 PRECAUTIONS AND WARNINGS**Safety**

For *in vitro* diagnostic use only.

The human blood components included in this kit have been tested by European approved and/or FDA approved methods and found negative for HBsAg, anti-HCV, anti-HIV-1 and 2. No known method can offer complete assurance that human blood derivatives will not transmit hepatitis, AIDS or other infections. Therefore, handling of reagents, serum or plasma specimens should be in accordance with local safety procedures.

All animal products and derivatives have been collected from healthy animals. Bovine components originate from countries where BSE has not been reported. Nevertheless, components containing animal substances should be treated as potentially infectious.

Avoid any skin contact with all reagents, Stop Solution contains HCl. In case of contact, wash thoroughly with water.

Do not smoke, drink, eat or apply cosmetics in the working area. Do not pipette by mouth. Use protective clothing and disposable gloves.

16 SUMMARY OF THE PROTOCOL

	CALIBRATORS (μL)	SAMPLE(S) CONTROLS (μL)
Incubation buffer	50	50
Calibrators (0-5)	200	-
Samples, Controls	-	200
Incubate for 2 hours at 18 °C - 25 °C with continuous shaking at 700 rpm. Aspirate the contents of each well. Wash 3 times with 400 μ L of Wash Solution and aspirate.		
Zero Calibrator	100	100
Anti-TNF- α -HRP conjugate	50	50
Incubate for 2 hours at 18 °C - 25 °C with continuous shaking at 700 rpm. Aspirate the contents of each well. Wash 3 times with 400 μ L of Wash Solution and aspirate.		
Chromogenic Solution	100	100
Incubate for 15 min at 18 °C - 25 °C with continuous shaking at 700 rpm.		
Stop Solution	100	100
Read on a microtiter plate reader and record the absorbance of each well at 450 nm (and 490 nm) versus 630 (or 650 nm)		

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17 BLIOGRAPHY / LITERATUR / RIFERIMENTI BIBLIOGRAFICI / BIBLIOGRAFÍA

1. BEUTLER B., CERAMI A. (1987) **Cachectin : more than a tumor necrosis factor.**
N. Engl. J. Med., 316 ; 379-385.
2. TRACEY K.J., FONG Y., HESSE D.G., MANOGUE K.R., LEE A.T., KUO G.C., LOWRY S.F. and CERAMI A. (1987) **Anti cachectin/TNF monoclonal antibodies prevent septic shock during lethal bacteraemia.**
Nature, 330 : 662-664.
3. PIGUET P.F., GRAU G.E., ALLET B. and VASSALLI P. (1987) **Tumor necrosis factor/cachectin is an effector of skin and gut lesions of the acute phase of graft-versus-host disease.**
J. Exp. Med., 166 ; 1280-1289.
4. AUKRUST P., LIABAKK N-B., MÜLLER F., LIEN E., ESPEVIK T. and FROLAND S.S. (1994) **Serum Levels of Tumor Necrosis Factor- α (TNF α) and Soluble TNF Receptors in Human Immunodeficiency Virus Type 1 Infection - Correlations to Clinical Immunologic, and Virologic Parameters.**
J. Inf. Dis., 169:420-424.
5. WAAGE A., HALSTENSEN A. and ESPEVIK T. (1987) **Association between tumor necrosis factor in serum and fatal outcome in patients with meningococcal disease.**
Lancet, 1 ; 355-357.
6. LEROUX-ROELS G., OFFNER F., PHILIPPE J. and VERMEULEN A. (1988) **Influence of Blood-Collecting Systems on Concentrations of Tumor Necrosis Factor in Serum and Plasma.**
Clin. Chem., 34 ; 2373-2374.

SYMBOLS USED

Symbol	English	Deutsch	Italiano	Español	Français
	European Conformity	CE-Konformitäts-kennzeichnung	Conformità europea	Conformidad europea	Conformité normes européennes
	Consult instructions for use *	Gebrauchsanweisung beachten *	Consultare le istruzioni per l'uso	Consulte las instrucciones de uso	Consulter les instructions d'utilisation
	<i>In vitro</i> diagnostic medical device *	<i>In-vitro</i> -Diagnostikum *	Dispositivo medico-diagnostico in vitro	Producto sanitario para diagnóstico In vitro	Dispositif médical de diagnostic in vitro
	Catalogue number *	Artikelnummer *	Numero di Catalogo	Número de catálogo	Référence de catalogue
	Batch code *	Chargencode *	Codice del lotto	Codigo de lote	Numéro de lot
	Contains sufficient for <n> tests *	Ausreichend für <n> Prüfungen	Contenuto sufficiente per "n" saggi	Contenido suficiente para <n> ensayos	Contenu suffisant pour "n" tests
	Temperature limit *	Temperaturbegrenzung *	Temperatura di conservazione	Temperatura de conservación	Température de conservation
	Use-by date *	Verwendbar bis *	Utilizzare prima del	Establa hasta	Utiliser jusque
	Manufacturer *	Hersteller *	Fabbricante	Fabricante	Fabricant
	Caution *	Achtung *			
	For research use only	Nur für Forschungszwecke	Solo a scopo di ricerca	Sólo para uso en investigación	Seulement dans le cadre de recherches
<i>Distributed by</i>	Distributed by	Vertreiber	Distributore	Distribuidor	Distributeur
<i>Content</i>	Content	Inhalt	Contenuto	Contenido	Contenu
<i>Volume/No.</i>	Volume / No.	Volumen/Anzahl	Volume/Quantità	Volumen/Número	Volume/Quantité