# "DIAQUICK" hCG Dipstick (Human Chorionic Gonadotropin)

for urine samples

REF	Content
J01402-1CE	- 30 tests individually packed (30 x Ref. No: J01402-1B), 5 mm - 1 package insert
J01402-2CE	- 30 tests individually packed (30 x Ref. No: J01402-2B), 3 mm - 1 package insert
J01402-1B	- 1 test individually packed 5 mm - 1 package insert
J01402-2B	- 1 test individually packed 3 mm - 1 package insert

### For in vitro diagnostic use only

#### GENERAL INFORMATION

<b>Method</b>	sandwich type immunochromatographic assay
<b>Shelf life</b>	24 months from date of production
<b>Storage</b>	2-30°C
<b>Sample</b>	urine
<b>Results</b>	within 3-5 minutes at room temperature
<b>Sensitivity</b>	25 mIU/ml hCG

#### INTENDED USE

The "DIAQUICK" hCG Dipstick (urine) is a rapid chromatographic immunoassay for the qualitative detection of human chorionic gonadotropin (hCG) in urine, as an aid for the early detection of pregnancy.

#### SUMMARY

Human chorionic gonadotropin (hCG) is a glycoprotein hormone produced by the developing placenta shortly after fertilization. In normal pregnancy, hCG can be detected in both serum and urine as early as 7 to 10 days after conception. hCG levels continue to rise very rapidly, frequently exceeding 100 mIU/ml by the first missed menstrual period, and peaking in the 100 -200 IU/ml range about 10-12 weeks into pregnancy. The appearance of hCG in both the urine and serum soon after conception, and its subsequent rapid rise in concentration during early gestational growth, make it an excellent marker for the early detection of pregnancy. The "DIAQUICK" hCG Dipstick (urine) is a rapid urine test to qualitatively detect the presence of hCG in urine specimens at the sensitivity of 25 mIU/ml. The test utilizes a combination of monoclonal and polyclonal antibodies to selectively detect elevated levels of hCG in urine. At the level of claimed sensitivity the "DIAQUICK" hCG Dipstick (urine) shows no cross-reactivity interference from the structurally related glycoprotein hormones hFSH, hLH and hTSH at physiological levels.

#### TEST PRINCIPLE

The "DIAQUICK" hCG Dipstick (urine) is a rapid chromatographic immunoassay for the qualitative detection of human chorionic gonadotropin in urine to aid in the early detection of pregnancy. The test utilizes a combination of antibodies including a monoclonal hCG antibody to selectively detect elevated levels of hCG. The assay is conducted by immersing the test strip in a urine specimen and observing the formation of colored lines. The specimen migrates via capillary action along the membrane to react with the colored conjugate. Positive specimens react with the specific antibody-hCG-colored conjugate to form a colored line at the test line region of the membrane. Absence of this colored line suggests a negative result. To serve as a procedural control, a colored line will always appear in the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred. The presence of this red band serves as verification for sufficient sample volume and proper flow and as a control for the reagents.

#### REAGENTS

The dipstick contains anti-hCG particles and anti-hCG coated on the membrane.

#### WARNINGS AND PRECAUTIONS

- Do not use test strips beyond the expiration date.
- Test strips should remain in the sealed pouch until ready to use.
- All specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The test strips should be discarded in a proper biohazard container after testing.
- The used test should be discarded according to local regulations.

#### STORAGE

Store the test at room temperature or refrigerated (2-30°C). The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. DO NOT FREEZE. Do not use beyond the expiration date.

#### SAMPLE COLLECTION AND PREPARATION

**Urine Assay:** The urine specimen must be collected in a clean and dry container. The first morning urine is preferred since it generally contains the highest concentration of hCG. However, urine collected at any time

of the day may be used. Urine samples exhibiting visible precipitates should be centrifuged, filtered, or allowed to settle to obtain clear supernatant for testing.

**Specimen Storage:** Urine specimens may be stored at 2-8°C for up to 48 hours prior to assay. For prolonged storage, specimens may be frozen and stored below -20°C for up to 6 months with no adverse effect on test performance. Frozen specimens should be thawed and mixed before testing.

#### MATERIALS PROVIDED

- Dipsticks
- Instructions for use

#### MATERIALS REQUIRED BUT NOT PROVIDED

- Specimen collection container
- Timer

#### ASSAY PROCEDURE

**Allow test strips, urine or serum samples and/or controls to equilibrate to room temperature (15-25°C) prior to testing.**

- Remove the test strip from sealed pouch.
- Immerse the test strip in the urine samples with arrows pointing toward the urine. Be sure the sample level is below the "MAX line" on the test strip. See illustration below.
- Leave the test strip immersed for approx. 15 seconds. (the liquid front moving upward can be seen) Remove and place on a non-absorbent flat surface.
- Wait for red lines to appear. Read the result at 3-4 minutes. Do not interpret results after the appropriate read time. It is important that the background is clear before the result is read.


#### INTERPRETATION OF RESULTS

**POSITIVE:** Two distinct red bands will appear, one in the test region and another in the control region.

**NOTE:** A sample hCG concentration below the cut-off level of this test might result in a weak line appearing in the test region (T) after an extended period of time. A line in the test region (T) seen after the read time could be indicative of a low hCG level in the sample. If such results are seen, it is recommended that the test be repeated with a new sample in 48-72 hours or that an alternate confirmation method is used.

**NEGATIVE:** Only a single band appears in the control region. No apparent line appears in the test region.

**INVALID:** Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

**NOTES:** The "DIAQUICK" hCG Dipstick (urine) is a very sensitive test which detects hCG concentrations above 20-25 IU. Using this test it is possible to detect hCG concentrations up to 500,000 mIU/mL. At very high concentrations a high dose hook effect may lead to a faint or even non-detectable test line. Any sample dilution between 1:10 and 1:100 would result in a very distinct test line. Therefore any visible test line should be considered as a positive result.

#### QUALITY CONTROL

A procedural control is included in the test. Correct procedural technique and test strip performance is confirmed when a red band appears at the control region of the membrane. It confirms sufficient specimen volume and correct procedural technique. A clear background is an internal negative procedural control. If a background color appears in the result window and interferes with the ability to read the test result, the result may be invalid.

It is recommended that a positive hCG control (containing 25-250 mIU/ml hCG) and a negative hCG control (containing "0" mIU/ml hCG) should be tested to verify proper test performance when a new shipment of product is received.

### LIMITATIONS

1. The "DIAQUICK" hCG Dipstick (urine) is a preliminary qualitative test, therefore, neither the quantitative value nor the rate of increase in hCG can be determined by this test.
2. Very diluted urine specimens as indicated by low specific gravity, may not contain representative levels of hCG. If pregnancy is still suspected, a first morning urine sample should be obtained 48-72 hours later and tested.
3. False negative results may occur when levels of hCG are below the sensitivity level of the test.
4. Very low levels of hCG (less than 50 mIU/ml) are present in urine or serum shortly after implantation. However, because a significant number of first trimester pregnancies terminate for natural reasons, a test result that is weakly positive should be confirmed by retesting with a first morning urine 48 hours later.
5. A number of conditions other than pregnancy, including trophoblastic disease and certain non-trophoblastic neoplasms including testicular tumors, prostate cancer, breast cancer, and lung cancer, cause elevated levels of hCG. Therefore, the presence of hCG in urine should not be used to diagnose pregnancy unless these conditions have been ruled out.
6. Certain medications alter the hormone level in urine.
7. This test provides a presumptive diagnosis for pregnancy. A confirmed pregnancy diagnosis should only be made by a physician after all clinical and laboratory findings have been evaluated.
8. This test reliably detects intact hCG up to 500,000 mIU/mL. It does not reliably detect hCG degradation products, including free-beta hCG and beta core fragments. Quantitative assays used to detect hCG may detect hCG degradation products and therefore may disagree with the results of this rapid test.

### EXPECTED VALUES

Negative results are expected in healthy non-pregnant women and healthy men. Healthy pregnant women have hCG present in both urine and serum. The amount will vary greatly with gestational age and between individuals. Urine and serum hCG concentration of pregnant women rise very rapidly after implantation, reaching a peak concentration in excess of 200 mIU/ml about 2-3 months after the last menstrual period. The "DIAQUICK" hCG Dipstick (urine) has a sensitivity of 25 mIU/ml for urine and is capable of detecting pregnancy as early as 1 day after the first missed menses.

### PERFORMANCE CHARACTERISTICS

#### Accuracy

A multi-center clinical evaluation was conducted comparing the results obtained using the "DIAQUICK" hCG Dipstick (urine) to another commercially available urine membrane hCG test. The study included 150 urine specimens: both assays identified 72 negative and 78 positive results. The results demonstrated >99.0% overall accuracy of the "DIAQUICK" hCG Dipstick (urine) when compared to the other urine membrane hCG test.

### hCG Reference Method

Method	EIA		Total Results	
	Results	Positive		Negative
DIAQUICK hCG Dipstick	Positive	78	0	78
	Negative	0	72	72
Total Results		78	72	150

Relative Sensitivity: 100,0% (95%-100%)\*  
 Relative Specificity: 100,0% (95%-100%)\*  
 Accuracy: 100,0% (98%-100%)\*  
 \* 95% Confidence Intervals

#### Sensitivity and Specificity

The "DIAQUICK" hCG Dipstick (urine) detects hCG at a concentration of 25 mIU/mL or greater. The test has been standardized to the W.H.O. Third International Standard. The addition of LH (300 mIU/mL), FSH (1,000 mIU/mL), and TSH (1,000 mIU/mL) to negative (0 mIU/mL hCG) and positive (25 mIU/mL hCG) specimens showed no cross-reactivity.

#### Interfering Substances

The following potentially interfering substances were added to hCG negative and positive specimens.

Acetaminophen	20 mg/mL	Caffeine	20 mg/mL
Acetylsalicylic Acid	20 mg/mL	Gentisic Acid	20 mg/mL
Ascorbic Acid	20 mg/mL	Glucose	2 g/dL
Atropine	20 mg/mL	Hemoglobin	1 mg/dL
Bilirubin	2 mg/dL		

None of the substances at the concentration tested interfered in the assay.

### REFERENCES

1. Batzer FR. Hormonal evaluation of early pregnancy, Fertil. Steril. 1980; 34(1): 1-13
2. Catt KJ, ML Dufau, JL Vaitukaitis Appearance of hCG in pregnancy plasma following the initiation of implantation of the blastocyste, J. Clin. Endocrinol. Metab. 1975; 40(3): 537-540
3. Braunstein GD, J Razor, H. Danzer, D Adler, ME Wade Serum human chorionic gonadotropin levels throughout normal pregnancy, Am. J. Obstet. Gynecol. 1976; 126(6): 678-681
4. Lenton EA, LM Neal, R Sulaiman Plasma concentration of human chorionic gonadotropin from the time of implantation until the second week of pregnancy, Fertil. Steril. 1982; 37(6): 773-778
5. Steier JA, P Bergsjö, OL Myking Human chorionic gonadotropin in maternal plasma after induced abortion, spontaneous abortion and removed ectopic pregnancy, Obstet. Gynecol. 1984; 64(3): 391-394
6. Dawood MY, BB Saxena, R Landesman Human chorionic gonadotropin and its subunits in hydatidiform mole and choriocarcinoma, Obstet. Gynecol. 1977; 50(2): 172-181
7. Braunstein GD, JL Vaitukaitis, PP Carbone, GT Ross Ectopic production of human chorionic gonadotropin by neoplasms, Ann. Intern Med. 1973; 78(1): 39-45

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