DIAQUICK FOB-TRF Combi Cassette

An immunochromatographic lateral flow device for the qualitative detection of human haemoglobin (hHb) and human transferrin (hTf) in human faecal samples

REF

Content

- J13101CE**
- 20 Tests, individually packed in foil pouches (20x REF J13101BN)
 - 20 Collection Tubes, filled with 2 mL Buffer each
 - 1 Package Insert

For professional in-vitro diagnostic use only

GENERAL INFORMATION

Method	immunochromatographic antigen detection
Shelf life	24 months from date of production
Storage	2-30 °C
Sample	fresh human faecal samples
Results	after 10 minutes at room temperature

TEST PRINCIPLE

The DIAQUICK FOB-TRF Combi Cassette is a qualitative, membrane based immunochromatographic lateral flow device for the detection of human haemoglobin (hHb) and human transferrin (hTf) in human faecal samples. The membrane is pre-coated with monoclonal antibodies against hHb and hTf in the test line region. After sample application, the sample reacts with latex particles coated with anti-hHb and anti-hTf antibodies stored in the conjugate pad. The mixture migrates upward on the membrane by capillary action. In case of a positive result the specific antibodies present on the membrane will react with the conjugate mixture and generate visibly coloured lines. To serve as a procedural control, a visibly coloured line will always appear in the control line region, indicating that the proper specimen volume was added and correct membrane wicking has occurred.

REAGENT COMPOSITION

FOB-TRF Combi Cassette

Strip, housed in a plastic cassette, with the following components:

- Nitro-cellulose membrane strip with coated monoclonal antibodies against human haemoglobin and human transferrin in the test line region and polyclonal antibodies in the control line region
- Conjugate pad with anti-human haemoglobin antibodies conjugated to blue coloured latex particles, anti-human transferrin antibodies conjugated to red coloured latex particles and a specific binding protein conjugated to green latex particles
- sample and absorbant pad

FOB-TRF Combi Buffer

- Biological buffer
- Sodium azide < 0.1 %

ADDITIONALLY REQUIRED MATERIAL

- Specimen collection equipment
- Pipettes with disposable tips for 15 µL (for liquid samples)
- General laboratory equipment

REAGENT STABILITY AND STORAGE

Conditions: The cassette must remain in the sealed aluminium pouch until use.

Storage: at 2-30 °C; Do not freeze!

Stability: up to the expiration date, when stored in the sealed pouch
 Do not use the test, if the pouch is damaged.


SAMPLE COLLECTION AND PREPARATION

Sample Collection

- Collect faeces in a clean, dry specimen collection container. Use a separate collection tube for each sample. Frozen samples must be totally thawed and brought to room temperature before testing.
- For a better detection of sporadic bleedings it is recommended to collect faecal samples from 3 days in a row and from 3 consecutive defecations.
- Specimen should not be collected during or within 3 days of a menstrual period or if the patient suffers from bleeding haemorrhoids or blood in the urine.
- Alcohol, aspirin and other medications taken in excess may cause GI irritation and occult bleeding and should be discontinued at least 48 hours prior to testing.
- Dietary restrictions are not necessary.

Sample Preparation

- Unscrew the cap of the specimen collection tube, then introduce the collection stick into the faecal specimen on at least 3 different sites to pick up a sufficient quantity (approx. 15 mg). Do not scoop the faecal sample.
- For liquid stool samples, aspirate the faecal specimen with a pipette and add approx. 15 µL into the faeces collection tube.
- Screw on and tighten the cap on the faeces collection tube.
- Vigorously shake the collection tube to mix the sample with the dilution buffer.

SAMPLE STABILITY AND STORAGE

The DIAQUICK FOB-TRF Combi Cassette must be used with human faecal samples.


Stability:	at 2-8 °C	undiluted	diluted in buffer
	at -20 °C	3 days	3 days
		2 months	--

The sample storage at room temperature cannot be recommended for this test. Avoid multiple freeze-thaw cycles, as this can cause erroneous results.

TEST PROCEDURE

1. Bring reagents and samples to room temperature (18-25 °C). The sensitivity of the test may be reduced at low temperatures.
2. Remove the DIAQUICK FOB-TRF Combi Cassette from the sealed pouch.
3. Vigorously shake the Collection Tube to assure good sample dispersion. Break the tip of the tube.
4. Add 4 drops of solution to the specimen well (S). Avoid adding solid particles with the liquid. Start the timer.
5. Read the result at **10 minutes** after dispensing the sample. False positive results can appear if the test is read after that time.

Note: if the test does not start to run due to an excess of solid particles in the solution, try to stir the sample in the specimen well (S) with the collection stick.

QUALITY CONTROL

A procedural control is included in the test. A green line appearing in the control line region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and a correct procedural technique.

Control standards are not supplied with this kit. However, it is recommended that positive and negative controls be tested as a good laboratory practice to confirm the test procedure and to verify proper test performance.

READING AND INTERPRETATION

POSITIVE RESULTS

hTf positive: A red line is visible in the test line region (T) and a green line in the control line region (C). This indicates the presence of hTf in the tested sample above 0.5 µg Tf/g faeces.

hHb positive: A blue line is visible in the test line region (T) and a green line in the control line region (C). This indicates the presence of hHb in the tested sample above 6 µg Hb/g faeces.

hTf & hHb positive: Both, a red and a blue line appear in the test line region (T) and a green line in the control line region (C). This indicates the presence of both hTf and hHb in the tested sample.

NEGATIVE RESULTS

Only the green band in the control line region (C) appears. This indicates that neither hHb nor hTf could be detected in the tested sample or are present below the test cut-off of 6 µg Hb/g faeces or 0.5 µg Tf/g faeces.

INVALID RESULTS

If the green coloured control line does not appear, the result is invalid. Review the package insert and repeat the test with a new device. If the problem persists, discontinue using the test kit and contact your local distributor.

Note: The colour intensity of the test lines will vary depending on the concentration of hHb and hTf in the sample. However, neither the quantitative value nor the rate of increase in hHb or hTf can be determined by this qualitative test.

DIAGNOSTIC IMPLICATIONS

Colorectal cancer (CRC) is an important health problem and ranges among the top leading causes of cancer deaths throughout the world. The disease is believed to develop mainly from non-malignant precursor lesions called adenomas. That CRC and adenomas bleed has provided faecal blood haemoglobin as the biomarker of choice for screening programmes. A significant technical enhancement to the simple guaiac test for blood is achieved by using immunological means to detect Hb. Whilst the haem component of blood, which is detected by the guaiac method, is common to all species, globin is species specific, so the immunological methods are not subject to interference from dietary blood.¹

Tf has also been identified as a potential biomarker for CRC. Bleeding into the gastrointestinal (GI) tract releases both Tf and Hb in a ratio of approx. 1:100 into faeces. Tf has shown greater stability in the faecal matrix and against degradation by bacteria and digestive enzymes, with a recovery ratio in faeces of approx. 1:5 compared to Hb. This makes Hb a more specific marker to bleeding from the distal GI tract, while Tf is a better indicator for a disease in the upper GI tract. The combination of Tf with Hb as faecal CRC biomarkers has shown to produce higher positivity rates in patients with CRC.²

WARNINGS AND PRECAUTIONS

- For professional in vitro diagnostic use only. Do not use after expiration date.
- Do not eat, drink or smoke in the area where the specimens or tests are handled.
- Handle all specimens as if they contained infectious agents. Observe established precautions against microbiological hazards throughout all procedures and follow the standard procedures for the proper disposal of specimens.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are assayed.
- Used tests should be discarded according to local regulations in a proper biohazard container.

PERFORMANCE CHARACTERISTICS

Analytical Sensitivity

Human Haemoglobin: 50 ng/mL buffer or 6 µg/g faeces
 Human Transferrin: 4 ng/mL buffer or 0.5 µg/g faeces

Note: Results in ng/mL cannot be directly compared across different faecal immunochemical tests (FIT). Results in ng/mL depend on the volume of buffer and the mass of the sample, which are both variable. For a direct data comparison, the results must be expressed in µg/g. The conversion can be achieved with the following formula:³

$$\mu\text{g/g faeces} = (\text{ng/mL buffer} \times \text{mL buffer}) / (\text{mg faeces collected})$$

Diagnostic Sensitivity and Specificity

The DIAQUICK FOB-TRF Combi Cassette was compared to other commercially available rapid tests, as well as to a guaiac assay. The diagnostic accuracy was > 99 %.

Cross-Reactivity and Interferences

The cross reactivity and interferences of the DIAQUICK FOB-TRF Combi Cassette was evaluated. No cross reactivity or interference could be detected with common GI pathogens, other organisms or substances occasionally present in faeces:

- | | |
|------------------------------|--------------------------------|
| - Rotavirus | - <i>Campylobacter</i> |
| - Astrovirus | - <i>Giardia lamblia</i> |
| - Adenovirus | - Human lactoferrin |
| - <i>Escherichia coli</i> | - Human transferrin (hHb line) |
| - Bovine and pig haemoglobin | - Human haemoglobin (hTf line) |
| - Bovine and pig transferrin | - Human calprotectin |

No special diet is required prior to testing. The test shows no interference with any food or supplements.

EXPECTED VALUES

Studies show that blood losses of > 1 mL/day may be seen following vigorous brushing of teeth and gums, and in irritation and inflammation of the intestinal tract. Intake of most non-steroidal anti-inflammatory drugs (NSAIDs) and aspirin in low doses usually produces an increased blood loss of 1-2 mL/day. Large aspirin doses of ≥ 1800 mg/day cause daily blood losses of 5-10 mL. Other chronic inflammatory conditions of the GI tract, including inflammatory bowel disease, colitis, Crohns' disease and perianal lesions, also increase blood loss.

It still remains a matter of conjecture whether all early-stage cancers bleed and whether they bleed intermittently, dependant perhaps upon the mechanics of the GI tract and the passage of digested foodstuffs. Intermittent or variable blood loss partially explains why the less-sensitive guaiac tests do not show consistently positive test results in patients, who are later diagnosed with CRC and why, even with highly sensitive tests, 100 % clinical sensitivity is not achieved.¹

WASTE MANAGEMENT


Please refer to local legal regulations.

LIMITATIONS

- The test will only qualitatively indicate the presence of human haemoglobin and human transferrin in human faecal samples. Neither the quantitative value nor the rate of increase in hHb or hTf concentration can be determined by this test.
- The DIAQUICK FOB-TRF Combi Cassette must be used with human faecal samples only.
- This test delivers preliminary results only. All results must be confirmed with alternative clinical tests and information available to the attending physician. A definitive clinical diagnosis should only be made by the physician after all clinical and laboratory finding have been evaluated.
- Various conditions can cause a positive result due to faecal blood loss. Therefore, all positive results should be followed up with additional diagnostic tests by a physician to determine the exact cause and source of the blood in the faeces. Here, endoscopy is the method of choice in diagnosing the cause of upper and lower GI bleeding.
- Negative results do not exclude bleeding since some adenomas and CRCs may bleed intermittently or not at all during certain stages of the disease. Additionally, blood may not be uniformly distributed in faecal samples.
- An excess of sample can cause wrong results (brown bands appear). In this case, dilute less faeces in the buffer and repeat the test.
- Urine, detergents, deodorants and water from the toilet bowl can cause erroneous test results. The use of a stool catcher for the toilet bowl is recommended.
- The test devices are for one-time use only. Use a separate cassette for each sample.

REFERENCES

- European Commission. European guidelines for quality assurance in colorectal cancer screening and diagnosis. First Edition, 2010
- Damian WLL *et al.* Evaluation of the analytical performance of the novel NS-Prime system and examination of temperature stability of fecal transferrin compared with fecal haemoglobin as biomarkers in a colon cancer screening program. *Practical Laboratory Medicine* 2 (2015) 29-36
- Fraser CG *et al.* A Proposal to Standardize Reporting Units for Fecal Immunochemical Tests for Hemoglobin. *J Natl Cancer Inst* 2012;104:810-814
- Chien-Hua Chiang *et al.* A Comparative Study of Three Fecal Occult Blood Tests in Upper Gastrointestinal Bleeding. *Kaohsiung J Med Sci* May 2006, Vol 22, No 5: 223-228
- Ji-Gui Chen *et al.* Colorectal cancer screening: Comparison of transferrin and immuo fecal occult blood test. *World J Gastroenterol* 2012 June 7; 18(21): 2682-2688
- Walker CW. Fecal Occult Blood Tests reduce Colorectal Cancer Mortality. *Am Fam Physician.* 2007 Jun 1;75(11):1652-3.
- Hirofumi Miyoshi *et al.* Accuracy of Detection of Colorectal Neoplasia using an Immunochemical Occult Blood Test in Symptomatic Referred Patients: Comparison of Retrospective and Prospective Studies. *Internal Medicine* Sept. 2000 Vol. 39, No. 9: 701-706.
- Towler BP *et al.* Screening for Colorectal Cancer using the Faecal Occult Blood Test, Hemocult. *Cochrane Database Syst Rev.* 2000(2): CD001216

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
BIBLIOGRAPHIE

1. European Commission. European guidelines for quality assurance in colorectal cancer screening and diagnosis. First Edition, 2010
2. Damian WLL *et al.* Evaluation of the analytical performance of the novel NS-Prime system and examination of temperature stability of fecal transferrin compared with fecal haemoglobin as biomarkers in a colon cancer screening program. *Practical Laboratory Medicine* 2 (2015) 29-36
3. Fraser CG *et al.* A Proposal to Standardize Reporting Units for Fecal Immunochemical Tests for Hemoglobin. *J Natl Cancer Inst* 2012;104:810-814
4. Chien-Hua Chiang *et al.* A Comparative Study of Three Fecal Occult Blood Tests in Upper Gastrointestinal Bleeding. *Kaohsiung J Med Sci* May 2006, Vol 22, No 5: 223-228
5. Ji-Gui Chen *et al.* Colorectal cancer screening: Comparison of transferrin and immuo fecal occult blood test. *World J Gastroenterol* 2012 June 7; 18(21): 2682-2688
6. Walker CW. Fecal Occult Blood Tests reduce Colorectal Cancer Mortality. *Am Fam Physician.* 2007 Jun 1;75(11):1652-3.
7. Hirofumi Miyoshi *et al.* Accuracy of Detection of Colorectal Neoplasia using an Immunochemical Occult Blood Test in Symptomatic Referred Patients: Comparison of Retrospective and Prospective Studies. *Internal Medicine* Sept. 2000 Vol. 39, No. 9: 701-706.
8. Towler BP *et al.* Screening for Colorectal Cancer using the Faecal Occult Blood Test, Hemocult. *Cochrane Database Syst Rev.* 2000(2): CD001216

1.