“DIAQUICK” Multi-Drug Panels

for human urine samples

Multi-3 Drug Panel	BZO,COC,MOP
- REF Z06576CE	Cont.: 30 panels, individually packed (30x REF Z06576B)
- REF Z06576B	Cont.: 1 panel, individually packed
Multi-3/1 Drug Panel	BUP, MOP, MTD
- REF Z09577CE	Cont.: 30 panels, individually packed (30x REF Z06577B)
- REF Z09577B	Cont.: 1 panel, individually packed
Multi-4 Drug Panel	AMP,COC,MOP,THC
- REF Z02575CE	Cont.: 30 panels, individually packed (30x REF Z02575B)
- REF Z02575B	Cont.: 1 panel, individually packed
Multi-5 Drug Panel	BZO,COC,MET,MOP,THC
- REF Z05236CE	Cont.: 30 panels, individually packed (30x REF Z05236B)
- REF Z05236B	Cont.: 1 panel, individually packed
Multi-5/3 Drug Panel	AMP,COC,MET,MOP,THC
- REF Z06502CE	Cont.: 30 panels, individually packed (30x REF Z06502B)
- REF Z06502B	Cont.: 1 panel, individually packed
Multi-5/4 Drug Panel	AMP,COC,MDMA,MOP,THC
- REF Z11504CE	Cont.: 30 panels, individually packed (30x REF Z11504B)
- REF Z11504B	Cont.: 1 panel, individually packed
Multi-5/6 Drug Panel	AMP,BZO,COC,MOP,THC
- REF Z06506CE	Cont.: 30 panels, individually packed (30x REF Z06506B)
- REF Z06506B	Cont.: 1 panel, individually packed
Multi-6 Drug Panel	BZO,COC,MET,MOP,MTD,THC
- REF Z98907CE	Cont.: 30 panels, individually packed (30x REF Z98907B)
- REF Z98907B	Cont.: 1 panel, individually packed
Multi-6/1 Drug Panel	AMP,BZO,COC,MET,MOP,THC
- REF Z03220CE	Cont.: 30 panels, individually packed (30x REF Z03220B)
- REF Z03220B	Cont.: 1 panel, individually packed
Multi-6/3 Drug Panel	BUP, BZO, COC, MTD, OPI, THC
- REF Z08930CE	Cont.: 30 panels, individually packed (30x REF Z08930B)
- REF Z08930B	Cont.: 1 panel, individually packed
Multi-6/4 Drug Panel	AMP,BUP,BZO,MET,MOP,THC
- REF Z08940CE	Cont.: 30 panels, individually packed (30x REF Z08940B)
- REF Z08940B	Cont.: 1 panel, individually packed
Multi-6/6 Drug Panel	BUP,COC,MET,MOP,MTD,THC
- REF Z13960CE	Cont.: 30 panels, individually packed (30x REF Z13960B)
- REF Z13960B	Cont.: 1 panel, individually packed
Multi-6/7 Drug Panel	BUP,BZO,COC,MOP,MTD,THC
- REF Z09970CE	Cont.: 30 panels, individually packed (30x REF Z09970B)
- REF Z09970B	Cont.: 1 panel, individually packed
Multi-6/10 Drug Panel	AMP,BZO,COC,MOP,MTD,THC
- REF Z11911CE	Cont.: 30 panels, individually packed (30x REF Z11911B)
- REF Z11911B	Cont.: 1 panel, individually packed
Multi-7 Drug Panel	AMP,BUP,BZO,COC,MTD,MOP,THC
- REF Z12730CE	Cont.: 30 panels, individually packed (30x REF Z12730B)
- REF Z12730B	Cont.: 1 panel, individually packed
Multi-10 Drug Panel	AMP,BAR,BZO,COC,MDMA,MET,MOP,MTD,TCA,THC
- REF Z04230CE	Cont.: 30 panels, individually packed (30x REF Z04230B)
- REF Z04231CE	Cont.: 10 panels, individually packed (10x REF Z04230B)
- REF Z04230B	Cont.: 1 panel, individually packed
Multi-10/1 Drug Panel	AMP,BAR,BZO,BUP,COC,MDMA,MET,MOP,MTD,THC
- REF Z06235CE	Cont.: 30 panels, individually packed (30x REF Z05235B)
- REF Z06236CE	Cont.: 10 panels, individually packed (10x REF Z05235B)
- REF Z05235B	Cont.: 1 panel, individually packed
Multi-10/2 Drug Panel	AMP,BAR,BZO,COC,MDMA,MOP,MTD,OPI,TCA,THC
- REF Z06102CE	Cont.: 30 panels, individually packed (30x REF Z06102B)
- REF Z06102B	Cont.: 1 panel, individually packed
Multi-10/3 Drug Panel	AMP, BZO, COC, MDMA, MOP, MTD, OPI, PCP, TCA, THC
- REF Z06103CE	Cont.: 30 panels, individually packed (30x REF Z06103B)
- REF Z06103B	Cont.: 1 panel, individually packed
Multi-10/4 Drug Panel	AMP,BAR,BUP,BZO,COC,MDMA,MET,MTD,OPI,THC
- REF Z06104CE	Cont.: 30 panels, individually packed (30x REF Z06104B)
- REF Z06104B	Cont.: 1 panel, individually packed
Multi-10/5 Drug Panel	AMP,BAR,BZO,BUP,COC,MET,MOP,MTD,TCA,THC
- REF Z06105CE	Cont.: 30 panels, individually packed (30x REF Z06105B)
- REF Z06105B	Cont.: 1 panel, individually packed
Multi-10/6 Drug Panel	AMP,BAR,BZO,COC,MET,MOP,MTD,PCP,TCA,THC
- REF Z06106CE	Cont.: 30 panels, individually packed (30x REF Z06106B)
- REF Z06106B	Cont.: 1 panel, individually packed
Multi-10/7 Drug Panel	AMP,BAR,BZO,COC,MET,MOP,MTD,PCP,TCA,THC
- REF Z06107CE	Cont.: 30 panels, individually packed (30x REF Z06107B)
- REF Z06107B	Cont.: 1 panel, individually packed

All products contain a package insert!

For in vitro diagnostic use only. For use by medical professionals only.
 For diagnosis and therapeutic monitoring only.

INTENDED USE

The “DIAQUICK” Multi-Drug Panels (urine) are rapid, lateral flow chromatographic immunoassays for the simultaneous, qualitative detection of the following drugs and their metabolites:

Parameter	Short	Calibrator Substance	Cut-off
Amphetamine	AMP	d-Amphetamine	1,000 ng/mL
Barbiturates	BAR	Secobarbital	300 ng/mL
Buprenorphine	BUP	Buprenorphine	10 ng/mL
Benzodiazepines	BZO	Oxazepam	300 ng/mL
Cocaine	COC	Benzoylcegonine	300 ng/mL
Ecstasy	MDMA	d,l-Methylenedioxymethamphetamine	500 ng/mL
Methamphetamine	MET	d-Methamphetamine	1,000 ng/mL
Methodone	MTD	Methodone	300 ng/mL
Opiate, Morphine, Heroine	MOP	Morphine	300 ng/mL
Opiate, Morphine, Heroine	OPI	Morphine	2,000 ng/mL
Phencyclidine	PCP	Phencyclidine	25 ng/mL
Tricyclic Antidepressants	TCA	Nortriptyline	1,000 ng/mL
Marihuana/Cannabis	THC	11-nor- Δ^9 -THC-9-COOH	50 ng/mL

This test will detect other related compounds, please refer to the Analytical Specificity table in this insert. This assay provides only a preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are obtained. For in vitro diagnostic use only

TEST PRINCIPLE

The “DIAQUICK” Multi-Drug Panels (urine) are immunoassays based on the principle of competitive binding. Drugs which may be present in the urine specimen compete against their respective drug conjugate for binding sites on their specific antibody. During testing, a urine specimen migrates upward by capillary action. A drug, if present in the urine specimen below its cut-off concentration, will not saturate the binding sites of its specific antibody coated on the particles. The antibody coated particles will then be captured by the immobilized drug conjugate and a visible colored line will show up in the test line region of the specific drug strip. The colored line will not form in the test line region if the drug level is above its cut-off concentration because it will saturate all the binding sites of the antibody coated on the particles. A drug-positive urine specimen will not generate a colored line in the specific test line region of the strip because of drug competition, while a drug-negative urine specimen or a specimen containing a drug concentration less than the cut-off will generate a line in the test line region. To serve as a procedural control, a colored line will always appear at the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred.

WARNINGS AND PRECAUTIONS

- For medical and other in vitro diagnostic use only. Do not use after the expiration date.
- The test panel should remain in the sealed pouch until use.
- All specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The used test panels should be discarded according to federal, state and local regulations.

REAGENTS

Each test line contains anti-drug mouse monoclonal antibody and corresponding drug-protein conjugates. The control line contains goat anti-rabbit IgG polyclonal antibodies and rabbit IgG.

STORAGE

The “DIAQUICK” Multi-Drug Panels can be stored refrigerated or at room temperature (2-30°C). The test panel is stable through the expiration date printed on the sealed pouch. The test panel must remain in the sealed pouch until use. DO NOT FREEZE. Do not use beyond the expiration date.

SAMPLE COLLECTION AND PREPARATION


The urine must be collected in a clean and dry container. Urine collected at any time of the day may be used. Urine specimens exhibiting visible precipitations should be centrifuged, filtered or allowed to settle to obtain a clear specimen for testing. Urine specimens may be stored at 2-8°C for up to 48 h prior to testing. For prolonged storage, specimens may be frozen and stored below -20°C. Frozen specimens should be thawed and mixed before testing.

ASSAY PROCEDURE

Allow the test panel, urine specimen, and/or controls to equilibrate to room temperature (15-30°C) prior to testing

1. Remove the test panel from the sealed pouch and use it as soon as possible.

2. Take off the protective cap plugged on the test panel. With arrows pointing towards the urine specimen, immerse the test panel vertically into the urine specimen for 10-15 seconds. Do not allow the urine sample to touch the plastic cassette when immersing the test device into the urine sample. Avoid immersion of the cassette deeper than the mark indicated with the arrows on the device and avoid any direct contact of the sample with the test region.

3. Put the protective cap back onto the test panel. Place the test panel on a non-absorbent flat surface, start the timer and wait for the red line(s) to appear. Read the results at 5 minutes. Do not interpret results after 10 minutes.

INTERPRETATION OF RESULTS

NEGATIVE: A colored line in the control region (C) and a colored line in the test line region (T) for a specific drug indicate a negative results. This indicates that the drug concentration in the urine specimen is below the designated cut-off level for that specific drug.

*NOTE: The shade of color in the test region (T) may vary, but it should be considered negative whenever there is even a faint pink line.

POSITIVE: A colored line in the control line region (C) but no line in the test line region (T) for a specific drug indicates a positive results. This indicates that the drug concentration in the urine specimen exceeds the designated cut-off level.

INVALID: Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test using a new test panel. If the problem persists, discontinue using the lot immediately and contact your local distributor

QUALITY CONTROL

A procedural control is included in the test. A red line appearing in the control region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique. Control standards are not supplied with this kit; however, it is recommended that positive and negative controls be tested as good laboratory practice to confirm the test procedure and to verify proper test performance.

LIMITATIONS

1. The “DIAQUICK” Multi-Drug Panels (urine) provide only a preliminary analytical result. A more specific chemical method must be used to obtain a confirmed result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method.
2. It is possible that technical or procedural errors, as well as other interfering substances in the urine specimen may cause erroneous results.
3. Adulterants, such as bleaching agents in urine specimens may produce erroneous results regardless of the analytical method used. If adulteration is suspected, the test should be repeated with another urine specimen.
4. A positive result indicates presence of the drug or its metabolites but does not indicate

- the level of intoxication, administration route or concentration in urine.
- A negative result may not necessarily indicate drug-free urine. Negative results can be obtained if a drug is present but below the cut-off level of the test.
 - The "DIAQUICK" Multi-Drug Panels (urine) do not distinguish between drugs of abuse and certain medications.
 - A positive result might be obtained from certain foods or food supplements.

PERFORMANCE CHARACTERISTICS

ACCURACY

A side-by-side comparison was conducted using the "DIAQUICK" Multi-Drug Panels (urine) and commercially available drug rapid tests. Testing was performed on approximately 250 specimens. Presumptive positive results were confirmed by GC/MS. The following results were tabulated.

% Agreement with Commercial Kit

	Positive Agreement	Negative Agreement	Total Results
AMP	> 99.9 %	> 99.9 %	> 99.9 %
BAR	> 99.9 %	> 99.9 %	> 99.9 %
BUP	> 99.9 %	> 99.9 %	> 99.9 %
BZO	> 99.9 %	> 99.9 %	> 99.9 %
COC	> 99.9 %	> 99.9 %	> 99.9 %
THC	> 99.9 %	> 99.9 %	> 99.9 %
MTD	> 99.9 %	> 99.9 %	> 99.9 %
MET	> 99.9 %	> 99.9 %	> 99.9 %
MDMA	> 99.9 %	> 99.9 %	> 99.9 %
MOP	> 99.9 %	> 99.9 %	> 99.9 %
OPI	*	*	*
PCP	> 99.9 %	> 99.9 %	> 99.9 %
TCA	*	*	*

* NOTE: Based on GC/MS instead of Commercial Kit.

% Agreement with GC/MS

	Positive Agreement	Negative Agreement	Total Results
AMP	98.4 %	98.4 %	98.4 %
BAR	97.3 %	99.3 %	98.4 %
BUP	99.2 %	> 99.9 %	99.6 %
BZO	98.4 %	99.2 %	98.8 %
COC	98.0 %	99.3 %	98.8 %
THC	> 99.9 %	99.1 %	99.6 %
MTD	99.1 %	98.6 %	98.8 %
MET	98.9 %	93.6 %	95.6 %
MDMA	96.2 %	97.5 %	96.8 %
MOP	96.3 %	95.7 %	96.0 %
OPI	94.3 %	91.3 %	92.8 %
PCP	95.1 %	94.6 %	94.8 %
TCA	94.8 %	89.6 %	92.0 %

ANALYTICAL SENSITIVITY

A drug-free urine pool was spiked with drugs to the concentrations at ± 50% cut-off and ± 25% cut-off. The results are summarized below:

Drug Conc.	n	AMP		BAR		BUP		BZO		COC	
		-	+	-	+	-	+	-	+	-	+
-50% Cut-off	30	30	0	30	0	30	0	30	0	30	0
-25% Cut-off	30	27	3	27	3	26	4	27	3	26	4
Cut-off	30	14	16	16	14	14	16	15	15	15	15
+25% Cut-off	30	3	27	4	26	3	27	4	26	3	27
+50% Cut-off	30	0	30	0	30	0	30	0	30	0	30

Drug Conc.	n	MTD		MET		MDMA		MOP		OPI	
		-	+	-	+	-	+	-	+	-	+
-50% Cut-off	30	30	0	30	0	30	0	30	0	30	0
-25% Cut-off	30	27	3	26	4	25	5	26	4	27	3
Cut-off	30	15	15	14	16	14	16	15	15	15	15
+25% Cut-off	30	3	27	3	27	4	26	3	27	4	26
+50% Cut-off	30	0	30	0	30	0	30	0	30	0	30

Drug Conc.	n	PCP		TCA		THC	
		-	+	-	+	-	+
-50% Cut-off	30	30	0	30	0	30	0
-25% Cut-off	30	26	4	25	5	26	4
Cut-off	30	14	16	15	15	14	16
+25% Cut-off	30	3	27	3	27	3	27
+50% Cut-off	30	0	30	0	30	0	30

CROSS-REACTIVITY

A study was conducted to determine the cross-reactivity of the test with compounds in either drug-free urine or drug positive urine. The following compounds did not show a cross-reactivity when tested with the "DIAQUICK" Multi-Drug Panels (urine) at a concentration of 100 µg/mL.

Non Cross-Reacting Compounds:

Acetophenetidin	Cortisone	Ketamine	d-Pseudoephedrine
N-Acetylprocainamide	Creatinine	Ketoprofen	Quinidine
Acetylsalicylic acid	Deoxycoflicosterone	Labetalol	Quinine
Aminopyrine	Dextromethorphan	Loperamide	Salicylic acid
Amoxicillin	Diclofenac	Meprobamate	Serotonin
Ampicillin	Diffunisal	Methoxyphenamine	Sulfamethazine
l-Ascorbic acid	Digoxin	Methylphenidate	Sulindac
Apomorphine	Diphenhydramine	Nalidixic acid	Tetracycline
Aspartame	Ethyl-p-aminobenzoate	Naproxen	Tetrahydrocortisone,
Atropine	β-Estradiol	Niacinamide	3-acetate
Benzilic acid	Estrone-3-sulfate	Nifedipine	Tetrahydrocortisone
Benzoic acid	Erythromycin	Norethindrone	Tetrahydrozoline
Bilirubin	Fenoprofen	Noscapine	Thiamine
d,l-Brompheniramine	Furosemide	d,l-Octopamine	Thioridazine
Caffeine	Genticic acid	Oxalic acid	d,l-Tyrosine
Cannabidiol	Hemoglobin	Oxolinic acid	Tolbutamide
Chloralhydrate	Hydralazine	Oxymetazoline	Triamterene
Chloramphenicol	Hydrochlorothiazide	Papaverine	Trifluoperazine
Chlorothiazide	Hydrocortisone	Penicillin-G	Trimethoprim
d,l-Chlorpheniramine	o-Hydroxyhippuric acid	Perphenazine	d,l-Tryptophan
Chlorpromazine	3-Hydroxytyramine	Phenazine	Uric acid
Cholesterol	d,l-Isoproterenol	Prednisone	Verapamil
Clonidine	Isosuprine	d,l-Propranolol	Zomepirac

ANALYTICAL SPECIFICITY

The following tables lists the concentration of compounds (ng/mL) that are detected positive in urine by the "DIAQUICK" Multi-Drug Panels (urine) at 5 minutes.

AMPHETAMINES (AMP)	AMP	ECSTASY	MDMA
d,l-Amphetamine sulfate	200	(±) 3,4-Methylenedioxyamphetamine HCl	500
l-Amphetamine	25,000	(±) 3,4-Methylenedioxyamphetamine HCl (MDA)	3,000
d-Amphetamine	1,000	3,4-Methylenedioxyethylamphetamine (MDE)	300
(±) 3,4-Methylenedioxyamphetamine	400	OPIATE, MORPHINE, HEROINE	MOP
Maprotiline	50,000	Codeine	200
Methoxyphenamine	6,000	Ethylmorphine	6,000
Phentermine	800	Hydrocodone	50,000
BARBITURATES	BAR	Hydromorphone	3,000
Allobarbitol	450	Lorphanol	1,500
Amobarbitol	3,000	6-Monoacetylmorphine	400
Alphenol	300	Morphine 3-β-d-glucuronide	800
Aprobarbitol	450	Morphine	300
Barbitol	6,000	Norcodeine	6,000
Butabarbitol	150	Normorphone	50,000
Butalbital	6,000	Oxycodone	30,000
Butethal	450	Oxymorphone	50,000
Cyclopentobarbitol	25,000	Procaine	15,000
5,5-Diphenylhydantoin	6,000	Thebaine	6,000
Pentobarbitol	6,000	MARIJUANA/CANNABIS	THC
Phenobarbitol	300	Cannabinol	20,000
Secobarbitol	300	11-nor-Δ ⁸ -THC-9 COOH	30
Talbutal	30	11-nor-Δ ⁹ -THC-9 COOH	50
BENZODIAZEPINES	BZO	Δ ⁸ -THC	15,000
Alprazolam	100	Δ ⁹ -THC	15,000
a-hydroxyalprazolam	1,500	OPIATE, MORPHINE, HEROINE	OPI
Bromazepam	780	Codeine	2,000
Chlordiazepoxide	780	Ethylmorphine	3,000
Clobazam	200	Hydrocodone	50,000
Clonazepam	390	Hydromorphone	12,500
Clorazepate dipotassium	390	Lorphanol	25,000
Delorazepam	780	6-Monoacetylmorphine	3,000
Desalkylflurazepam	200	Morphine 3-β-d-glucuronide	2,000
Diazepam	1,500	Morphine	2,000
Estazolam	6,250	Norcodeine	25,000
Flunitrazepam	200	Normorphone	50,000
(±) Lorazepam	3,100	Oxycodone	25,000
RS-Lorazepam glucuronide	200	Oxymorphone	25,000
Midazolam	6,250	Procaine	50,000
Nitrazepam	200	Thebaine	25,000
Norchlordiazepoxide	100	PHENCYCLIDINE	PCP
Nordiazepam	780	Phencyclidine	25
Oxazepam	300	4-Hydroxyphencyclidine	6,250
Temazepam	100	TRICYCLIC ANTIDEPRESSANTS	TCA
Triazolam	3,100	Amitriptyline	1,500
BUPRENORPHINE	BUP	Clomipramine	50,000
Buprenorphine	10	Cyclobenzaprine	1,500
Buprenorphine 3-d-Glucuronide	50	Desipramine	200
Norbuprenorphine	50	Doxepine	1,500
3-d-Glucuronide	100	Imipramine	400
COCAINE	COC	Maprotiline	1,500
Benzoylcegonine	300	Nordoxepine	400
Cocain HCl	200	Nortriptyline	1,000
Cocaethylene	12,500	Perphenazine	25,000
Ecgonine HCl	30,000	Promazine	3,000
METHADONE	MTD	Promethazine	25,000
Methadone	300	Trimipramine	3,000
Doxylamine	100,000		
METHAMPHETAMINE	MET		
p-Hydroxymethamphetamine	25,000		
d-Methamphetamine	1,000		
l-Methamphetamine	12,500		
(±) 3,4-Methylenedioxyamphetamine	6,250		
Mephentermine	50,000		

REFERENCES

- Hawks RL, CN Chiang. Urine Testing for Drugs of Abuse. National Institute for Drug Abuse (NIDA), Research Monograph 73, 1986.
- Tietz NW. Textbook of Clinical Chemistry. W.B. Saunders Company. 1986; 1735.
- Stewart DJ, Inaba T, Lucassen M, Kalow W. Clin. Pharmacol. Ther. April 1979; 25 ed: 464, 264-8.
- Ambre J. J. Anal. Toxicol. 1985; 9:241.
- Winger, Gail, A Handbook of Drug and Alcohol Abuse, Third Edition, Oxford Press, 1992, page 146.
- Robert DeCresce. Drug Testing in the workplace, 1989 page 114.
- Glass, IB. The International Handbook of Addiction Behavior. Routledge Publishing, New York, NY. 1991; 216
- B. Cody, J.T., "Specimen Adulteration in drug urinalysis. Forensic Sci. Rev., 1990, 2:63.
- C. Tsai, S.C. et.al., J. Anal. Toxicol. 1998; 22 (6): 474
- Baselt RC. Disposition of Toxic Drugs and Chemicals in Man. 6th Ed. Biomedical Publ., Foster City, CA 2002.
- Hardman JG, Limbird LE. Goodman and Gilman's: The Pharmacological Basis for Therapeutics. 10th Edition. McGraw Hill Medical Publishing, 2001; 208-209.

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BIBLIOGRAPHIE

- Hawks RL, CN Chiang. Urine Testing for Drugs of Abuse. National Institute for Drug Abuse (NIDA), Research Monograph 73, 1986.
- Tietz NW. Textbook of Clinical Chemistry. W.B. Saunders Company. 1986; 1735.
- Stewart DJ, Inaba T, Lucassen M, Kalow W. Clin. Pharmacol. Ther. April 1979; 25: eda: 464, 264-8.
- Ambre J. J. Anal. Toxicol.1985; 9:241.
- Winger, Gail, A Handbook of Drug and Alcohol Abuse, Third Edition, Oxford Press, 1992, page 146.
- Robert DeCresce. Drug Testing in the workplace, 1989 page 114.
- Glass, IB. The International Handbook of Addiction Behavior. Routledge Publishing, New York, NY. 1991; 216
- B. Cody, J.T., "Specimen Adulteration in drug urinalysis. Forensic Sci. Rev., 1990, 2:63.
- C. Tsai, S.C. et.al., J. Anal. Toxicol. 1998; 22 (6): 474
- Baselt RC. Disposition of Toxic Drugs and Chemicals in Man. 6th Ed. Biomedical Publ., Foster City, CA 2002.
- Hardman JG, Limbird LE. Goodman and Gilman's: The Pharmacological Basis for Therapeutics. 10th Edition. McGraw Hill Medical Publishing, 2001; 208-209.