

# D-Dimer for DIAcheck & Coalyser

Diagnostic Reagent for the quantitative in vitro determination of D-Dimer in human plasma on DIAcheck devices and Coalyser.

REF	Cont.	
T099702	1 x 5 mL	D-Dimer Latex Reagent
	1 x 7 mL	D-Dimer Reaction Buffer
	1 x 7 mL	Diluent
	1 x 1 mL	D-Dimer Calibrator (Iyo.)
	1 x 1 mL	D-Dimer Control Low (Iyo.)
T099703	1 x 1 mL	D-Dimer Control High (Iyo.)
	2 x 5 mL	D-Dimer Latex Reagent
	2 x 7 mL	D-Dimer Reaction Buffer
	1 x 7 mL	Diluent
	1 x 1 mL	D-Dimer Calibrator (Iyo.)
	1 x 1 mL	D-Dimer Control Low (Iyo.)
	1 x 1 mL	D-Dimer Control High (Iyo.)

For professional in vitro diagnostic use only.

## GENERAL INFORMATION

Method: Immunospectrometric Assay  
 Wavelength: 400-600 nm  
 Temperature: 37°C  
 Sample: Sodium Citrate Plasma

## Number of Tests

T099702 62 Tests on DIAcheck  
 T099703 125 Tests on DIAcheck

## DIAGNOSTIC IMPLICATION AND TEST PRINCIPLE

Fibrin fragments containing D-Dimer antigen is always present in plasma as a result of plasmin degradation of cross-linked fibrin. After an injury, or when suffering from conditions associated with increased haemostatic activity, there is an increase in plasma D-Dimer concentration. The determination of D-Dimer has become a prevalent aid in the diagnosis of thrombosis. Elevated levels of D-Dimer are found in clinical conditions such as deep vein thrombosis (DVT), pulmonary embolism (PE) and disseminated intravascular coagulation (DIC)<sup>1-4</sup>. A negative D-Dimer test result from a patient with a suspected thrombotic disorder has a high negative predictive value. Dialab D-Dimer for DIAcheck & Coalyser consists of sub-micron sizes polystyrene particles coupled to monoclonal antibodies specific for D-Dimer. When the reagent is exposed to a D-Dimer containing plasma sample, the particles will agglutinate, giving rise to increased light-scattering. When exposed to the appropriate wavelength of light, the increase in measured turbidity, or light-scattering, is proportional to the amount of D-Dimer in the sample.

## REAGENT COMPOSITION AND PREPARATION

**D-Dimer Latex Reagent:** liquid, ready to use, polystyrene particles, coated with monoclonal antibodies, suspended in buffer with stabilisers and sodium azide (<0.1%) as a preservative. As the micro-particles will settle during storage, swirl the vial gently a few times each day before it is used, to ensure a homogenous suspension. Do not shake.

**D-Dimer Reaction Buffer:** liquid, ready to use, containing buffer and sodium azide (<0.1%) as preservative.

**Diluent:** 0.9% Saline solution with sodium azide (<0.1%), ready to use

**D-Dimer Calibrator, D-Dimer Control Low/High:** contain lyophilized citrated plasma of human origin enriched with D-Dimer. Reconstitute with 1mL distilled water. Keep the reconstituted calibrator/control at 15-25°C for 15-30 minutes and verify that the lyophilised cake is completely dissolved before use.

## REAGENT STABILITY AND STORAGE

Conditions: protect from light  
 Storage: at 2 - 8°C  
 Stability: up to the expiration date shown on label

## Opened/reconstituted reagents:

Storage:	at 2 - 8°C	8 – 25°C
Stability of Latex & Buffer	4 weeks	2 weeks
Stability of Calibrator	10 hours	10 hours
Stability of Controls	54 hours	54 hours

## SAMPLE COLLECTION AND STORAGE

Venous blood is collected in 0.11 or 0.13 M trisodium citrate at a ratio of 9 parts blood to 1 part anticoagulant (1:10 ratio). The ratio is critical. If using commercial vacuum tubes, a full draw must be assured. Trauma or stasis during blood sampling should be avoided. The presence of a clot in a specimen is a cause for rejection. Refer to CLSI guideline H21-A5 for further instructions on specimen collection, handling and storage 9.

Plasma samples can be stored at room temperature (18-25°C) for up to 4 hours; refrigerated (2-8°C) for up to 4 hours; frozen at -20°C for up to 2 weeks or at -70°C for up to 6 months. Frozen samples should be thawed rapidly and tested immediately. If testing cannot be performed immediately, the sample may be kept refrigerated (2-8°C) for maximally 2 hours prior to testing. No contact with glass should occur.

## TEST PROCEDURE

	DIAcheck C1, C2, C4	Coalyser
Pipette into an empty cuvette	20 µL plasma	60 µL Plasma
Pipette into cuvette	40 µL Buffer	40 µL Buffer
Incubate at 37°C	2 minutes, activate channel	2 minutes
Pipette into cuvette	80 µL pre-warmed Latex Reagent and mix at least 5 x with pipette. Reading will begin automatically and mOD values will be displayed. Results are displayed in dOD (E)/min.	120 µL pre-warmed Latex

## CALIBRATION

For the calibration curve prepare a dilution series of the calibrator with Diluent (saline solution, 0.9%). Use at least 3 calibration points. A new calibration curve must be run for each new lot of reagent and if the control values are outside the assigned range.

## EXPECTED RESULTS

The D-Dimer results should be used together with other clinical and diagnostic information for forming a diagnosis. The normal level of D-Dimer in the population is typically below 200 ng/mL<sup>1,5</sup>. However, as there is no internationally established standard for D-Dimer, the concentration of D-Dimer in any given specimen may differ when determined using D-Dimer assays from different manufacturers. Thus, each laboratory should establish its own reference range and cut-off values. Elevated levels of D-Dimer are found in patients with deep venous thrombosis (DVT), pulmonary embolism, disseminated intravascular coagulation and trauma<sup>6</sup>. D-Dimer levels increase during pregnancy<sup>7</sup> and with age<sup>8</sup>. D-Dimer results can be reported in units of D-Dimer (ng/mL) or in Fibrinogen Equivalent Units (FEU). 1 ng/mL D-Dimer is approximately 2 FEU, although a stoichiometric calibration would suggest a different theoretical conversion factor.

## QUALITY CONTROL

D-Dimer Control Low and D-Dimer Control High or other commercially DD control plasma should be used for reliable quality control of performance at a frequency in accordance with good laboratory practice (GLP).

## LIMITATIONS AND INTERFERENCE

D-Dimer for DIAcheck & Coalyser is not affected by UF and LMW Heparin up to 100 U/mL, by Bilirubin up to 0,1 g/L, by Triglycerides up to 2,5 g/L and by Haemoglobin up to 4 g/L.

Specimens from patients who have received preparations of mouse monoclonal antibodies for diagnosis or therapy may contain anti-mouse antibodies (HAMA). Such antibodies may cause over-estimation of D-Dimer levels. The presence of rheumatoid arthritis factor may result in falsely elevated D-Dimer values. Turbid or opalescent plasma may cause erratic results and should be interpreted with caution; dilute the sample and re-assay.

The monoclonal antibody in D-Dimer for DIAcheck & Coalyser has been screened for its specificity against cross-linked fibrin degradation products. D-Dimer for DIAcheck & Coalyser has more than 100-fold specificity for D-Dimer (Fibrin or purified D-Dimer), over Fibrinogen, Fibrinogen D or Fragment E.

## WARNINGS AND PRECAUTIONS

- Calibrator and controls contain material of human origin.
- The plasma used in the production is tested free of antibodies to HIV I and II, Hepatitis B and Hepatitis C. No test can however completely exclude the presence of infected material and the product should be treated as potentially infectious.
- Waste is disposed of according to local regulations.
- Wear appropriate clothing.
- Avoid contact with skin and eyes.

## PERFORMANCE CHARACTERISTICS

### Precision:

Within run precision was assessed over multiple runs using specific lots of reagents and control.

The coefficient of variation obtained in this study was ≤ 3%.

## WASTE MANAGEMENT

Please refer to local legal requirements.

## REFERENCES

1. Heit, J.A et al. Determinants of plasma fibrin D-dimer sensitivity for acute pulmonary embolism as defined by pulmonary angiography. Arch pathol Lab Med, 123:235-239, 1999
2. Bounameaux H et al. Plasma measurement of D-dimer as diagnostic aid in suspected venous thromboembolism: an overview. Thromb Haemostas, 71:1-6,1994
3. Pfitzner S.A. et al. Fibrin detected in plasma of patients with disseminated intravascular coagulation by fibrin-specific antibodies consists primarily of high molecular weight factor XIII-crosslinked and plasminmodified complexes partially containing fibrinopeptide A Thromb Haemost 78:1069-1078, 1997
4. Lindahl T.L. et al. Clinical evaluation of diagnostic strategy for deep venous thrombosis with exclusion by low plasma levels of fibrin degradation product D-dimer. Scand J Lab Invest, 58:307-316, 1998
5. Gardiner, C. et al. An evaluation of rapid D-dimer assays for the exclusion of deep vein thrombosis. British Journal of Haematology, 128:842-848, 2005
6. Meissner, M.H., Venous thromboembolism in trauma: a local manifestation of systemic hypercoagulability? J. Trauma, 54(2):224-231,2003.
7. Ballegeer, V. et al. Fibrinolytic response to venous occlusion and fibrin fragment D-dimer levels in normal and complicated pregnancy. Thromb Haemostasis 58: 1030-1032, 1987
8. Kario, K. et al Which factors affect high D-dimer levels in the elderly? Thromb Res, 65(5):501- 508, 1991
9. CLSI. Collection, Transport and Processing of Blood Specimens for testing Plasma-Based Coagulation Assays, 5th Ed, CLSI document H21-A5, Vol. 28 No. 5

## LOT-SPECIFIC DATA

	Lot	Unit	Value (DDU)
D-Dimer Calibrator	17388	ng/mL	2972
D-Dimer Control Low	16021	ng/mL	300
D-Dimer Control High	16004	ng/mL	830

