

## Additive for wash solution

For use instruments with automatic cuvette washer.

Product code:	<b>VA0002SL</b>	2x250ml
	VA0000SL-3	250 ml
	VA0000SL-3	250 ml

For *in vitro* diagnostic use only

## INTENDED USE

Wash solution additive contains concentrated detergents, for regular cleaning and rinsing of clinical chemistry analyzers' probes, cuvettes, and tubings removing reagents, samples and protein deposits.

Additive for wash solution is only compatible with Pictus clinical chemistry analyzers.

Please refer to the instrument's User Manual for further information.

## INGREDIENTS

### 2x 250ml Additive for wash solution / Solution 3

Surfactants .....	< 2.5%
Stabilizers .....	< 0.3%
Preservatives .....	< 0.3%

## HANDLING PRECAUTIONS

- All body fluid samples should be considered potentially infectious materials. Treat all blood and other potentially infectious materials with appropriate precautions. Use gloves, masks and gowns if blood exposure is anticipated.
- Please refer to the MSDS associated with the reagent.
- Use Good Laboratory Practices (GLP) when handling the reagents.

## REAGENT PREPARATION

- Use 2 ml per liter of distilled water.
- Do not shake. Wait 5 minutes for complete mixing.
- Leave the reagent at room temperature for at least 12 hours.

## REAGENT REPLACEMENT

- Person installing the reagent must be a trained laboratory professional.
- Avoid any dust or microbial contamination of the tubing and reagents. Do not pour and mix the remains of a reagent from a container into other one.

## REAGENT STORAGE, STABILITY AND DISPOSAL

- Store wash solution additive between +4 and 32°C.
- The shelf life of the wash solution additive is 24 months from the date of manufacture, if stored at the prescribed temperature range.
- Do not use reagent beyond the expiration date printed on the container label.
- Discard opened container after 90 days.
- Do not use reagent once frozen.
- Dispose of waste product, unused product and contaminated packaging in compliance with local regulations.

In case of any issue contact your local official Diatron Distributor for troubleshooting.

## MANUFACTURER

Diatron MI Plc.  
Táblás utca 39.  
H-1097 Budapest, Hungary  
Tel: +36 12794400  
[www.diatron.com](http://www.diatron.com)

Printing date 18.05.2022  
version 3.0  
Reviewed on 18.05.2022