

Cuvette Cleaner kit

Cuvette Cleaner solution for the automatic washing of cuvettes for specific Pictus® series clinical chemistry analyzers

Product code: **VA0008CCL** 2x100ml
 VA0008CC 100 ml

For *in vitro* diagnostic use only

INTENDED USE

Cuvette Cleaner solution is used for the automatic decontamination and cleaning of cuvettes after each use and for the removal of any reagent, blood component, urine or CSF and lipoprotein deposits.

Cuvette Cleaner solution is only compatible with Pictus clinical chemistry analyzers.

Please refer to the instrument's User Manual for further information.

INGREDIENTS

Mixed amphoteric detergents

Sodium hydroxide < 2.0%

HANDLING PRECAUTIONS

- All body fluid samples should be considered potentially infectious materials. Treat all blood and other potentially infectious materials with appropriate precautions. Use gloves, masks and gowns if blood exposure is anticipated.
- Please refer to the MSDS associated with the reagent.
- Use Good Laboratory Practices (GLP) when handling the reagents.

REAGENT PREPARATION

- In order to prepare working solution; mix carefully 2 mL of Cuvette Cleaner with 1 L of Type 2 laboratory water (preferably deionized, below $\leq 8 \mu\text{S/cm}$).
- Leave the reagent at room temperature for at least 12 hours.

REAGENT REPLACEMENT

- Person installing the reagent must be a trained laboratory professional.
- Avoid any dust or microbial contamination of the tubing and reagents. Do not pour and mix the remains of a reagent from a container into other one.

REAGENT STORAGE, STABILITY AND DISPOSAL

- Store Cleaning reagent between +4 and 32°C.
- The shelf life of the Cuvette Cleaner solutions are 24 months from the date of manufacture, if stored at the prescribed temperature range.
- In use stability is valid until end of shelf life.
- Do not use reagent beyond the expiration date printed on the container label.
- Once opened, recap bottles as soon as possible to avoid solution/reagent deterioration.
- Do not use reagent once frozen.
- The solutions should not be used after prolonged exposure to sunlight or high temperature, or when they appear turbid.
- Dispose of waste product, unused product and contaminated packaging in compliance with local regulations.

In case of any issue contact your local official Diatron Distributor for troubleshooting.

Symbols, signal word, hazard and precautionary statements on VA0008CC (Cuvette Cleaner solution):



Warning

H290, H315, H319,
 P262, P273, P280,
 P302+P352,
 P305+P351+P338

MANUFACTURER

Diatron MI Plc.
 Táblás utca 39.
 H-1097 Budapest, Hungary
 Tel: +36 12794400
 www.diatron.com

Printing date 18.05.2022
 version 3.0
 Reviewed on 18.05.2022