HbA1c Direct

(en) English

REF

REF	Content
Y04602	4x 22.5 mL R1 + 1x 30 mL R2
Y04606	4x 7.5 mL R1 + 1x 10 mL R2
Y04601	4x 45 mL R1 + 1x 60 mL R2
YA1004	1x 30 mL R1 + 1x 10 mL R2
YT1204	1x 30 mL R1 + 1x 10 mL R2
Y06911	1x 30 mL R1 + 1x 10 mL R2
Y06917	1x 30 mL R1 + 1x 10 mL R2
YE1804	2x 30 mL R1 + 1x 20 mL R2
AB20420	1x 60 mL R1 + 1x 20 mL R2
AB20409	4x 60 mL R1 + 4x 20 mL R2

For professional in vitro diagnostic use only.

INTENDED USE

HbA1c Direct is a diagnostic reagent for quantitative in vitro determination of HbA1c in human blood by turbidimetric assay. The determination of HbA1c is most commonly performed for the evaluation of glycaemic control in diabetes mellitus.

DIAGNOSTIC SIGNIFICANCE

Throughout the circulatory life of the red cell, HbA1c is formed continuously by the adduction of glucose to the N-terminal of the hemoglobin beta chain. The process, which is non-enzymatic, reflects the average exposure of hemoglobin to glucose over an extended period. In a classical study, Trivelli et al¹ showed HbA1c in diabetic subjects to be elevated 2-3 folds over the levels found in normal individuals. HbA1c values reflect the average glucose value over the preceding 120 days.⁸ As a consequence a higher HbA1c value indicates poorer glycaemic control. Several investigators have recommended that HbA1c serves as an indicator of metabolic control of the diabetic, since HbA1c levels approach normal values for diabetics in metabolic control.^{2,3,4} HbA1c has been defined operationally as the "fast fraction" hemoglobins (HbA1a, A1b, A1c) that elute first during column chromatography with cation-exchange resins. The non-glycosylated hemoglobin, which consists of the bulk of the hemoglobin has been designated HbA₀.

TEST PRINCIPLE

This method utilizes the interaction of antigen and antibody to determine the HbA1c in whole EDTA blood. HbA1c in test samples is absorbed onto the surface of latex particles, which react with Anti-HbA1c (antigen-antibody reaction) and gives agglutination. The amount of agglutination is measured (increasing reaction) as absorbance at 660nm. The HbA1c value is obtained from a calibration curve.

REAGENT COMPOSITION

Reagent 1	
Latex	
Sodium azide	0.95 g/L
Reagent 2	
Anti-human HbA1c mouse monoclonal antibody	
Bis-Tris Buffer	
Stabilizers, Sodium azide	< 0.1%

MATERIAL REQUIRED BUT NOT PROVIDED

- Distilled water for the preparation of the hemolysate from samples, e.g.:

REF	Name	Composition	Content
Y04605	H-Reagent	Purified water, < 0.1% S. azide	5 x 100 mL

- Calibrator, e.g.:

REF	Name	Content
Y04603	HbA1c Direct Calibrator 5 Level Series*	5 x 1 mL

- Controls, e.g.:

REF	Name	Content
Y04604	HbA1c Direct Control Set*	4 x 1 mL

- Clinical Laboratory Analyser or Photometer
- General laboratory equipment.

* Hemolysates from packed human erythrocytes lyophilized and stabilized.

REAGENT PREPARATION

Reagent 1 and 2 are ready to use. Mix each bottle gently before use.

STORAGE AND STABILITY

Conditions:	Protect from light. Close immediately after use.
Storage:	at +2 – +8 °C
Stability:	until expiry, at least 24 months after production.

Stability after opening: R1 and R2 are stable for 1 month after opening if stored at +2 – +8 °C and closed tightly

DO NOT FREEZE REAGENTS AND STRICTLY AVOID ANY CONTACT WITH FREEZING TEMPERATURES! DO NOT USE ICE BAGS FOR TRANSPORT!

WARNINGS AND PRECAUTIONS

- For in vitro diagnostic use only.
- This assay should not be used for the diagnosis of diabetes mellitus.
- Freezing temperatures cause agglutination of Reagent R1: Strictly avoid every contact to freezing temperatures including the contact with the wall of your refrigeration system.
- Reagent deterioration: Alterations in the physical appearance (for example agglutination) of the reagents or values of control materials (for example results outside of the manufacturer's acceptable range) may indicate reagent instability. Do not use reagent with any sign of deterioration!
- Use quality controls with each test run and every time in case of an assumed risk of possible impairment to verify the correct functioning of your reagent and measuring system.
- Reagent 2 contains animal material. Handle the product as potentially infectious according to universal precautions and good clinical laboratory practices.
- Each donor unit used in the preparation of the calibrators and controls was found to be negative for the presence of HIV1 and HIV2 antibodies, as well as for the hepatitis B surface antigen and anti-hepatitis C antibodies, using a method approved by the FDA.
- The reagents contain sodium azide (0.95 g/L) as preservative. Do not swallow! Avoid contact with skin and mucous membranes! Sodium azide has been reported to form lead or copper azide in laboratory plumbing, which may explode on percussion. Flush drains with water thoroughly after disposing of fluids containing Sodium azide.
- Please refer to the safety data sheets and take the necessary precautions for the use of laboratory reagents.
- For diagnostic purposes, the results should always be assessed with the patients' medical history, clinical examinations and other findings.
- In the event of an incident related to the device, report it to the manufacturer and your competent authority as required.

SPECIMEN COLLECTION AND STORAGE

Use fresh EDTA whole blood.

To determine HbA1c, a hemolysate must be prepared for each sample:

- Dispense 2 mL distilled water into a test tube.
- Add 20 µL of well mixed EDTA whole blood into the tube and mix.
- Allow to stand for 5 minutes or until complete lysis is evident.

Hemolysate stability: at +2 – +8 °C 72 hours

Discard contaminated specimens.

TEST PROCEDURE

Solve calibrators and controls, do not lyse them. Perform lysis with samples only. Use the **HbA1c Direct Calibrator 5 Level Series** to generate a calibration curve. Assay-Temperature: 37 °C

General procedure, individual setting dependent on your used system:

Pipette into test tubes:	Calibrators	Samples/Controls
HbA1c Reagent 1	180 µL	180 µL
Calibrators/Controls/Samples	5 µL	5 µL
Mix. Incubate for 2 minutes at assay temperature. Then add:		
HbA1c Reagent 2	60 µL	60 µL
Mix and read A1 of calibrators and samples/controls at 660 nm. Incubate for 5 minutes, then read A2 at 660 nm. Calculate: $\Delta A = (A2 - A1)$		

AUTOMATION

Applications for automated systems are available upon request.

INTERPRETATION OF RESULTS

CALCULATION

HbA1c results according to NGSP for the samples and controls are determined using the prepared calibration curve.

For calculating results according to IFCC, use IFCC calibrator values (see insert calibrator series) or use the following equation:

$$IFCC (mmol/mol) = (HbA1c \text{ in } \% \text{ NGSP} - 2.15) / 0.0915$$

QUALITY CONTROL AND CALIBRATION

QUALITY CONTROL

Reliability of results should always be monitored if patient samples are tested. Quality control materials and calibrators must be used. All control material with % HbA1c values determined by this method can be used. We recommend the **DIALAB HbA1c Direct Control Set**. In case of controls of one test run should show abnormal results, patient results should not be valued.

CALIBRATION

The assay requires the use of a HbA1c calibrator. We recommend the **DIALAB HbA1c Direct Calibrator 5 Level Series**.

PERFORMANCE CHARACTERISTICS

LINEARITY

The Hemoglobin A1c assay range is appr. 0 % - 15 % (NGSP).

PRECISION [%CV]

	Low	Medium	High
Intra-Assay	1.14	0.72	1.22
Inter-Assay	1.85	2.05	2.97

ACCURACY [%]

Control	Assigned Value	Measured Value
DIALAB L1	5.8 (4.9 – 6.7)	5.9
DIALAB L2	10.0 (8.5 – 11.5)	9.9
BIORAD L1	5.48 (4.39 – 6.58)	5.58
BIORAD L2	9.42 (7.54 – 11.3)	8.35

METHOD COMPARISON

A comparison study of this Hemoglobin A1c procedure and ROCHE reagent gave the following results: $y = 0.9286x + 0.26 / r = 0.9656$.

ANALYTICAL SENSITIVITY

The change of absorbance units per concentration unit was calculated for the up going part of the nonlinear curve. Result was measured: 0.0458744 Absorbance/Concentration units.

Tests were performed on a Pentra 400 instrument.

TRACEABILITY

DIALAB HbA1c Direct reagent, when used with Dialab HbA1c Direct Calibrator 5 Level Series is traceable against NGSP certification, traceable to the Diabetes Control and Complications Trial (DCCT).

EXPECTED VALUES

According to NGSP: < 6 % for a non-diabetic
 < 7 % for glycemic control of a person with diabetes
 According to IFCC: < 42.0 mmol/mol Hb for a non-diabetic
 < 53.0 mmol/mol Hb for glycemic control of a person with diabetes

This range is for orientation only. Each laboratory establishes its own normal range. In using Hemoglobin A1c to monitor diabetic patients, results should be interpreted individually (the patient should be monitored against him or herself).

LIMITATIONS

No interference of:
 High Lipids: 3000 formazin turbidity units
 Free Bilirubin: 30 mg/dL
 Conjugated Bilirubin: 30 mg/dL

The reagent cannot detect HbF, therefore elevated levels of HbF can lead to underestimation of HbA1c results. Patients suffering from Thalassemia might show lower results due to short lifespan of HbA1c in samples. Due to the relation to red blood cell life, patients suffering from Iron Deficiency Anemia might show falsified results⁵⁻⁶.


Samples of patients treated with Dapsone (diaminodiphenyl sulfone) result in falsely low HbA1c values. For this case other methods of glycemic control, such as fructosamine are recommended⁷.

WASTE MANAGEMENT

Please refer to local legal requirements.

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- Kalibrator, z.B.:

REF	Name	Inhalt
Y04603	HbA1c Direct Calibrator 5 Level Series [*]	5 x 1 mL

- Kontrollen, z.B.:

REF	Name	Inhalt
Y04604	HbA1c Direct Control Set [*]	4 x 1 mL

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