HbA1c Direct Control Set (assayed control hemolysate)

Lyophilized control hemolysate for the use in tests for the quantitative in vitro determination of Hemoglobin A1c in blood on photometric systems

REF	Cont.		
Y04604	4x 1 mL	2x 1 mL Level 1 2x 1 mL Level 2	finished kit
Y04604SV	4x 1 mL	2x 1 mL Level 1 2x 1 mL Level 2	labelled vials

COMPOSITION

Hemolysate prepared from packed human erythrocytes. Stabilizers to maintain hemoglobin in the reduced state for the accurate calibration of the HbA1c procedure. The controls provide two levels of HbA1c, one level in the normal range and the other level in the elevated range.

PREPARATION

- Open the vial carefully, avoiding any loss of the lyophilized material.
- Add exactly 1 mL of deionized water (inaccurate reconstitution of the controls and error in assay technique can cause erroneous results).
- Close the vial carefully and gently mix for 10 minutes, or until all material has dissolved, avoiding the formation of foam. Do not shake!
- Reconstituted controls are now ready for use. Do not lyse controls.

STABILITY AND STORAGE

Storage:	at 2 – 8 °C
Stability:	up to the expiration date if sealed tightly
Stability after reconstitution:	30 days at 2 – 8 °C 3 months at -20 °C

CLOSE IMMEDIATELY AFTER USE. PROTECT FROM LIGHT AND HEAT. Do not use the product if there is visible evidence of microbial growth in the vial. Improper handling and/or storage of the calibrators can affect results.

EXPECTED VALUES

The assayed limits are to be used as a guide in determining the accuracy of the assay procedure. The assay results according to NGSP for the controls should fall within the stated expected range. For control values according to IFCC, please use values reported for IFCC method.

WARNINGS AND PRECAUTIONS

- For in Vitro Diagnostic use.
- Although this product has been tested and found non-reactive for HBsAg, HCV and anti-HIV 1&2, no known test can offer assurance that products derived from human blood will not transmit disease. Therefore all human serum products and patient specimens should be handled in the same manner as an infectious agent.
- Do not pipette by mouth. Avoid contact with skin and mucous membranes.

WASTE MANAGEMENT

Please refer to local legal requirements.

REFERENCES

- Trivelli, L.A., Ranney, H.M., and Lai, H.T., New Eng. J. Med. 284,353 (1971).
- Gonen, B., and Runenstein, A.H., Diabetologia 15, 1 (1978).
- Gabbay, K.H., Hasty, K., Breslow, J.L., Ellison, R.C., Bunn, H.F., and Gallop, P.M., J.Clin. Endocrinol. Metab. 44, 859 (1977).
- Bates, H.M., Lab. Mang., Vol 16 (Jan. 1978)

LOT SPECIFIC ASSAY DATA

Hemoglobin A1c control set with two levels (normal and elevated) of haemoglobin A1c. Values and expiry date are lot specific.

Method: direct immunoturbidimetric latex test.

LOT: 02359		EXP: 2019-06	
Assigned Values according to NGSP [%]			
Level	Target Value %	Range %	
Level 1	5.6	4.8 – 6.4	
Level 2	10.0	8.5 – 11.5	
Assigned Values according to IFCC [mmol/mol Hb] *			
Level	Target Value mmol/mol Hb	Range mmol/mol Hb	
Level 1	38	29 – 46	
Level 2	86	69 – 102	
* HbA1c in mmol/mol according to IFCC = (HbA1c in % NGSP – 2,15) / 0,0915			

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