G6PDH CALIBRATOR

Lyophilized calibration blood for the use in tests for the quantitative in vitro determination of Glucose-6-phosphate dehydrogenase (G6PDH) activity in human blood.

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COMPOSITION

Lyophilized human blood preparation with specified activity of Glucose-6-phosphate dehydrogenase.

CALIBRATOR PREPARATION

- Carefully remove the cap and rubber plug avoiding any loss of contents.
- Add exactly 0.5 mL of deionized water (inaccurate reconstitution of the calibrator can cause erroneous results).
- Close the vial carefully and allow the calibrator to stand for 5 minutes at room temperature
- Swirl the vial carefully and mix thoroughly making sure that all traces of dry material have dissolved. Avoid foaming. Do not shake!
- Again let stand for 5 minutes at room temperature.
- The calibrator should be treated in the same manner as patient samples.

TEST PROCEDURE

Please refer to the reagent package insert for instructions for use.

MATERIALS REQUIRED BUT NOT PROVIDED

Lyse Reagent, G6PDH reagent
 General laboratory equipment.

STABILITY AND STORAGE

Storage: at 2 – 6 °C
 Stability: until date of expiration date if sealed tightly.
 Stability after reconstitution: at 2 – 6 °C 10 days
 Close immediately after use. Avoid contamination!

WARNINGS AND PRECAUTIONS

- For professional use only.
- For in vitro diagnostic use only.
- Handle products with care and precision and with the same precautions used for human blood samples.
- Each individual donation intended for use in manufacture of the G6PDH Calibrator was tested for Hepatitis B surface antigen (HBsAg), Hepatitis C virus (HCV) and anti-HIV by FDA required tests. Only donations with negative findings were used for its manufacture. Nevertheless every product obtained from human blood should be considered potentially infectious.
- Refer to safety data for additional risk information

WASTE MANAGEMENT

Please refer to local legal requirements.

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