"DIAQUICK" COC Dipstick (Cocaine) for urine samples

REF	Content
Z02510CE	- 30 tests individually packed (30 x Ref. No: Z02510B). - 1 package insert
Z02510B	- 1 test individually packed - 1 package insert

**For in vitro diagnostic use only.
 For diagnosis and therapeutic monitoring only.
 For use by medical professionals only.**

GENERAL INFORMATION

Method	competitive immunochromatographic assay
Antibody	Anti-Benzoyllecgonine monoclonal antibody conjugate
Shelf life	24 months from date of production
Storage	2-30°C
Sample	Human Urine
Results	within 5 minutes at room temperature
Sensitivity	300 ng/ml.

INTENDED USE

The "DIAQUICK" COC Dipstick (urine) is a rapid qualitative, competitive binding immunoassay for determination of Cocaine in urine. The test provides only preliminary data, which should be confirmed by other methods (such as GC/MS). The test is intended for in vitro diagnostic use only.

CLINICAL SIGNIFICANCE

Cocaine is used primarily as a CNS stimulant drug. Cocaine is often self-administered by nasal inhalation, intravenous injection and free-base smoking. It is excreted in the urine in a short time primarily as Benzoyllecgonine. Benzoyllecgonine, a major metabolite of Cocaine, has a longer biological half-life (5 - 8 hours) than Cocaine (0.5 - 1.5 hours), and can generally be detected for 24-48 hours after Cocaine exposure. Methods historically used for detecting Cocaine in biological fluids include thin-layer chromatography, gas chromatography, ultraviolet spectroscopy, enzyme immunoassay and radio immunoassay. While confirmation techniques other than GC/MS may be adequate for some drugs of abuse, GC/MS is generally accepted as a vigorous confirmation technique for all drugs since it provides the best level of confidence in the result. The "DIAQUICK" Cocaine Test is an easy, fast and visually read assay without the need for instrumentation. The test system employs unique monoclonal and polyclonal antibodies to selectively identify Cocaine in urine samples with a high degree of sensitivity. The "DIAQUICK" COC Dipstick (urine) yields a positive result when the Cocaine metabolite in urine exceeds 300 ng/ml.

TEST PRINCIPLE

The "DIAQUICK" COC Dipstick (urine) is an immunoassay based on the principle of competitive binding. Drugs which may be present in the urine specimen compete against the drug conjugate for binding sites on the antibody. During testing, a urine specimen migrates upward by capillary action. Benzoyllecgonine, if present in the urine specimen below 300 ng/mL, will not saturate the binding sites of antibody in the test strip. The antibody coated particles will then be captured by immobilized Benzoyllecgonine conjugate and a visible colored line will appear in the test line region. The colored line will not form in the test line region if the Benzoyllecgonine level is above 300 ng/mL because it will saturate all the binding sites of antibodies. A drug-positive urine specimen will not generate a colored line in the test line region because of drug competition, while a drug-negative urine specimen or a specimen containing a drug concentration less than the cut-off will generate a line in the test line region. To serve as a procedural control, a colored line will always appear at the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred.

WARNINGS AND PRECAUTIONS

- Do not use kit beyond the expiration date.
- Samples may be infectious; properly handle and dispose of all used reaction devices in a biohazard container.

STORAGE

Store test kit at 2-30°C; do not freeze. Refer to the expiration date for stability.

SAMPLE COLLECTION AND PREPARATION

Collect a urine sample in clean, dry container, either plastic or glass, without any preservatives. Urine specimens may be refrigerated (2-8°C) and stored up to 2 days. For longer storage, freeze samples (-20°C or below). Bring frozen or refrigerated samples to room temperature before testing. Urine samples exhibiting visible precipitates should be filtered, centrifuged or allowed to settle. Use only clear aliquots for testing.

ASSAY PROCEDURE

Bring an unopened foil pouch and the sample specimen to room temperature.

- Open the foil pouch by tearing along the notch and remove the test strip.
- Dip the stick in the urine specimen for 10-15 seconds submerging only up to the max mark.
- Place the test strip on a non-absorbent flat surface and start the timer
- Read the test result within 5 minutes.


Important: Do not interpret a test result after more than 10 minutes. Waiting longer than five minutes may cause inaccurate interpretation. To avoid confusion, discard the test device after reading the result.

INTERPRETATION OF RESULTS

- Positive Result:**
A colored band is visible in the Control Zone "C" but not in the Test Zone "T". A positive result indicates that the Benzoyllecgonine concentration is at or above the sensitivity of the test.
- Negative Result:**
A colored band is visible in the Control Zone "C", and in the Test Zone "T". A negative result indicates that the Benzoyllecgonine concentration is below the sensitivity of the test.
- Invalid Result:**
There is no distinct color band visible neither in the Test Zone "T" nor in the Control Zone "C". The test should be considered invalid and it is recommended that the sample be retested.

Note: A very faint line on the test region indicates that the Benzoyllecgonine in the sample is near or below the cut-off level for the test. However, if any line is observed at the test region, the sample is NEGATIVE.

QUALITY CONTROL

An internal procedure control has been incorporated into the test to ensure proper kit performance and reliability. The use of an external control is recommended to verify proper kit performance. Quality control samples should be tested according to quality control requirements established by the testing laboratory.

LIMITATIONS

- This product is designed for use with human urine only.
- Although the test is very accurate, there is a possibility false results will occur due to the presence of interfering substances in the urine.
- The test is a qualitative assay and is not for determining quantitative concentration levels or the level of intoxication.
- Adulterants such as bleach or other strong oxidizing agents, when added to urine specimens, may produce erroneous test results regardless of the analysis method used. If adulteration is suspected, obtain another urine specimen and retest.

PERFORMANCE CHARACTERISTICS

Accuracy

A side-by-side comparison was conducted by laboratory personnel using the "DIAQUICK" COC Dipstick (urine) and a commercially available rapid test. Presumptive positive results were confirmed by GC/MS. The following results were tabulated:

Method	Other Rapid Test		Total Results
	Positive	Negative	
"DIAQUICK" COC Dipstick	Positive	0	136
	Negative	157	164
Total Results	143	157	300
% Agreement with this commercial kit	95%	>99%	98%

When compared to GC/MS at the cut-off of 300 ng/mL, the following results were tabulated:

Method	GC/MS	Total Results	
		Positive	Negative
"DIAQUICK" COC Dipstick	Results		
	Positive	119	17
	Negative	5	159
Total Results		124	176
% Agreement with GC/MS Analysis		96%	90%

Eighty (80) of these clinical samples were also run using the "DIAQUICK" COC Dipstick (urine) by an untrained operator at a different site. Based on GC/MS data, the operator obtained a statistically similar Positive Agreement, Negative Agreement and Overall Agreement rate as the laboratory personnel.

A drug-free urine pool was spiked with Benzoyllecgonine at the following concentrations: 0 ng/mL, 150 ng/mL, 225 ng/mL, 300 ng/mL, 375 ng/mL and 450 ng/mL. The result demonstrates 100% accuracy at 50% above and 50% below the cut-off concentration. The data are summarized below:

Benzoyllecgonine Concentration (ng/ml)	Percent of Cut-off	n	Visual Result	
			Negative	Positive
0	0	30	30	0
150	-50%	30	30	0
225	-25%	30	30	0
300	Cut-off	30	4	26
375	+25%	30	0	30
450	+50%	30	0	30

Analytical Specificity

The following table lists compounds that are positively detected in urine by the "DIAQUICK" COC Dipstick (urine) at 5 minutes.

Compound	Concentration (ng/mL)
Benzoyllecgonine	300
Cocaine HCl	780
Cocaethylene	12,500
Ecgonine HCl	32,000

Precision

A study was conducted at three physicians' offices by untrained operators using three different lots of product to demonstrate the within run, between run and between operator precision. An identical panel of coded specimens were blind labeled and tested at each site. The results are given below:

Benzoyllecgonine concentration (ng/ml)	n	Site 1		Site 2		Site 3	
		+	-	+	-	+	-
0	15	0	14*	0	15	0	15
150 ng/mL	15	1	14	0	15	1	14
225 ng/mL	15	11	4	10	5	7	8
375 ng/mL	15	15	0	15	0	15	0
450 ng/mL	15	15	0	15	0	14	1
Non Valid	15	16/16		15/15		15/15	

*Note: Non-valid results were obtained in this treatment. Non-valid tests were provided as part of this study to ensure that readers would accurately identify non-valid test results.

Cross-Reactivity

A study was conducted to determine the cross-reactivity of the test with compounds in either drug-negative urine or Benzoyllecgonine positive urine. The following compounds show no interference when tested with the "DIAQUICK" COC Dipstick (urine) at a concentration of 100 µg/mL.

Non Cross-Reacting Compounds

Acetaminophen	Fenpropfen	Oxymetazoline
Acetophenetidin	Furosemide	Papaverine
N-Acetylprocainamide	Genistic acid	Penicillin-G
Acetylsalicylic acid	Hemoglobin	Pentobarbital
Aminopyrine	Hydralazine	Perphenazine
Amitypyline	Hydrochlorothiazide	Phencyclidine
Amobarbital	Hydrocodone	Phenelzine
Amoxicillin	Hydrocortisone	Phenobarbital
Ampicillin	O-Hydroxyhippuric acid	Phentermine
L-Ascorbic acid	p-Hydroxy-methamphetamine	L-Phenylephrine
DL-Amphetamine Sulfate	3-Hydroxytyramine	β-Phenylethylamine
Apomorphine	Ibuprofen	Phenylpropanolamine
Aspartame	Imipramine	Prednisolone
Atropine	Iproniazid	Prednisone
Benzilic acid	(±) - Isoproterenol	Procaine
Benzoic acid	Isoxsuprine	Promazine
Benzphetamine	Ketamine	Promethazine
Bilirubin	Ketoprofen	DL-Propranolol
(±) - Brompheniramine	Labetalol	D-Propoxyphene
Caffeine	Levorphanol	D-Pseudoephedrine
Cannabidiol	Loperamide	Quinidine
Cannabitol	Maprotiline	Quinine
Chloralhydrate	Meperidine	Ranitidine
Chloramphenicol	Meprobamate	Salicylic acid
Chlordiazepoxide	Methadone	Secobarbital
Chlorothiazide	Methoxyphenamine	Serotonin
(±) - Chlorpheniramine	(±)-3,4-Methylenedioxy-amphetamine	Sulfamethazine
Chlorpromazine	hydrochloride	Sulindac
Chlorquine	(±)-3,4-Methylenedioxy-meth-amphetamine hydrochloride	Temazepam
Cholesterol	Morphine-3-β-D-glucuronide	Tetracycline
Clomipramine	Morphine Sulfate	Tetrahydrocortisone 3-(β-D-glucuronide)
Clonidine	Nalidixic acid	Tetrahydrozoline
Cocaine	Naloxone	Thebaine
Cortisone	Naltrexone	Thiamine
(-) Cotine	Naproxen	Thioridazine
Creatinine	Niacinamide	DL-Tyrosine
Dextroscortosterone	Nifedipine	Tolbutamide
Dextromethorphan	Nifedipine	Triamterene
Diazepam	Norcocaine	Trifluoperazine
Diclofenac	Norethindrone	Trimethoprim
Diffenosal	D-Norpropoxyphene	Trimipramine
Digoxin	Doxylamine	Tryptamine
Diphenhydramine	DL-Octopamine	DL-Tryptophan
Doxylamine	Oxalic acid	Tyramine
Ecgonine methylester	Oxazepam	Uric acid
(-) - Ψ-Ephedrine	Oxolinic acid	Verapamil
Erythromycin	Oxycodone	Zomepirac
β-Estradiol		
Estrone-3-sulfate		
Ethyl-p-aminobenzoate		

REFERENCES

- Stewart DI, T Inoba, M Ducassen, W Kalow. Clin. Pharmacol. Ther. 1979; 25:264
- Ambre J. J. Anal. Toxicol. 1985; 9:241
- Baselt RC. Disposition of Toxic Drugs and Chemicals in Man. 2nd Ed. Biomedical Publ., Davis, CA. 1982: 488
- Hawks RL, CN Chiang. Urine Testing for Drugs of Abuse. National Institute for Drug Abuse (NIDA), Research Monograph 73, 1986

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