"DIAQUICK" MTD Dipstick (Methadone) for urine samples

REF

Content

Z02550CE	- 30 tests individually packed (30 x Ref. No: Z02550B). - 1 package insert
Z02550B	- 1 test individually packed - 1 package insert

**For in vitro diagnostic use only.
 For diagnosis and therapeutic monitoring only.
 For use by medical professionals only.**

GENERAL INFORMATION

Method	competitive immunochromatographic assay
Antigen / Antibodies	monoclonal anti-methadone-antibody-coupled particles Methadone-protein conjugate
Shelf life	24 months from date of production
Storage	2-30°C
Sample	human urine
Results	within 5 minutes at room temperature
Cut-off	300 ng/mL

INTENDED USE

The "DIAQUICK" MTD Dipstick (urine) is a lateral flow chromatographic immunoassay for the detection of Methadone in human urine at a cut-off concentration of 300 ng/mL. This assay provides only a qualitative, preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are used. For in vitro diagnostic use only.

CLINICAL SIGNIFICANCE

Methadone is a narcotic pain reliever for medium to severe pain. It is also used in the treatment of Heroin (Opiate dependence: Vicodin, Percolate, Morphine, etc) addiction. Oral Methadone is very different than the IV Methadone. Oral Methadone is partially stored in the liver for later use. IV Methadone acts more like Heroin. In most states you must go to a pain clinic or a Methadone maintenance clinic to be prescribed Methadone.

Methadone is a long acting pain reliever producing effects that last between twelve to forty-eight hours. Ideally, Methadone frees the client from the pressures of obtaining illegal Heroin, from the dangers of injection, and from the emotional roller coaster that most Opiates produce. Methadone, if taken for long periods and at large doses, can lead to a very long withdrawal period. The withdrawals from Methadone are more prolonged and troublesome than those provoked by heroin cessation, yet the substitution and phased removal of methadone is an acceptable method of detoxification for patients and therapists.

The "DIAQUICK" MTD Dipstick (urine) is an easy to use urine test that can be performed without the use of an instrument. The test utilizes a monoclonal antibody to selectively detect elevated levels of Methadone in urine. The "DIAQUICK" MTD Dipstick (urine) yields a positive result when the Methadone in urine exceeds 300 ng/mL.

TEST PRINCIPLE

The "DIAQUICK" MTD Dipstick (urine) is a rapid chromatographic immunoassay based on the principle of competitive binding. Drugs which may be present in the urine specimen compete against the drug conjugate for binding sites on the antibody. During testing, a urine specimen migrates upward by capillary action. MTD, if present in the urine specimen below 300 ng/mL, will not saturate the binding sites of the antibody coated particles in the test strip. The antibody coated particles will then be captured by immobilized MTD conjugate and a visible colored line will show up in the test line region. The colored line will not form in the test line region if the MTD level is above 300 ng/mL because it will saturate all the binding sites of anti-MTD antibodies. A drug-positive urine specimen will not generate a colored line in the test line region because of drug competition, while a drug-negative urine specimen or a specimen containing a drug concentration less than the cut-off will generate a line in the test line region. To serve as a procedural control, a colored line will always appear at the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred.

MATERIALS REQUIRED BUT NOT PROVIDED

- Specimen collection container
- Timer

WARNINGS AND PRECAUTIONS

- For professional in vitro diagnostic use only. Do not use after the expiration date.
- The test strip should remain in the sealed pouch until use.
- All specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The used test strip should be discarded in a proper

biohazard container after testing.

STORAGE

The kit can be stored at room temperature or refrigerated (2-30°C). The test strip is stable through the expiration date printed on the sealed pouch. The test strip must remain in the sealed pouch until use. DO NOT FREEZE. Do not use beyond the expiration date.

SAMPLE COLLECTION AND PREPARATION

The urine specimen must be collected in a clean and dry container. Urine collected at any time of the day may be used. Urine specimens exhibiting visible precipitates should be centrifuged, filtered, or allowed to settle to obtain a clear supernatant for testing.

SPECIMEN STORAGE

Urine specimens may be stored at 2-8°C for up to 48 hours prior to testing. For prolonged storage, specimens may be frozen and stored below -20°C. Frozen specimens should be thawed and mixed before testing.

ASSAY PROCEDURE

Allow the test strip, urine specimen, and/or controls to equilibrate to room temperature (15-30°C) prior to testing.

- Remove the test strip from the sealed pouch and use it as soon as possible.
- With arrows pointing toward the urine specimen, immerse the test strip vertically in the urine specimen for at least 10-15 seconds. Do not pass the maximum line (MAX) on the test strip when immersing it. See the illustration below.
- Place the test strip on a non-absorbent flat surface, start the timer and wait for the red line(s) to appear. The result should be read at 5 minutes. Do not interpret the result after 10 minutes.


INTERPRETATION OF RESULTS

NEGATIVE:* Two lines appear. One red line should be in the control region (C), and another apparent red or pink line should be in the test region (T). This negative result indicates that the Methadone concentration is below the detectable level of 300 ng/mL.

***NOTE:** The shade of red in the test region (T) may vary, but it should be considered negative whenever there is even a faint pink line.

POSITIVE: One red line appears in the control region (C). No line appears in the test region (T). This positive result indicates that the Methadone concentration is above the detectable level of 300 ng/mL.

INVALID: Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test strip. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

QUALITY CONTROL

An internal procedure control has been incorporated into the test to ensure proper kit performance and reliability. The use of an external control is recommended to verify proper kit performance. Quality control samples should be tested according to quality control requirements established by the testing laboratory.

LIMITATIONS

1. This product is designed for use with urine only.
2. Although the test is very accurate, there is a possibility false results will occur due to the presence of interfering substances in the urine.
3. The test is a qualitative urine assay and is not for determining quantitative concentration levels or the level of intoxication.
4. Adulterants such as bleach or other strong oxidizing agents, when added to urine specimens, may produce erroneous test results regardless of the analysis method used. If adulteration is suspected, obtain another urine specimen and retest.
5. A negative result may not necessarily indicate drug-free urine. Negative results can be obtained when drug is present but below the cut-off level of the test.
6. This test does not distinguish between drugs of abuse and certain

medications.

PERFORMANCE CHARACTERISTICS

Accuracy

A side-by-side comparison was conducted using the "DIAQUICK" MTD Dipstick (urine) and a leading commercially available MTD rapid test. Presumptive positive results were confirmed by GC/MS. The following results were tabulated:

Method	Results	Other MTD Rapid Test		Total Results
		Positive	Negative	
DIAQUICK MTD Dipstick	Positive	132	0	132
	Negative	0	168	168
Total Results		132	168	300
% Agreement with this Rapid Test Kit		>99%	>99%	>99%

When compared to GC/MS at the cut-off of 300 ng/mL, the following results were tabulated:

Method	Results	GC/MS		Total Results
		Positive	Negative	
DIAQUICK MTD Dipstick	Positive	122	10	132
	Negative	1	167	168
Total Results		123	177	300
% Agreement with GC/MS Analysis		99%	94%	96%

Analytical Sensitivity

A drug-free urine pool was spiked with Methadone at the following concentrations: 0 ng/mL, 150 ng/mL, 225 ng/mL, 300 ng/mL, 375 ng/mL and 450 ng/mL. The result demonstrates >99% accuracy at 50% above and 50% below the cut-off concentration. The data are summarized below:

MTD Concentration (ng/mL)	Percent of Cut-off	n	Visual Result	
			Negative	Positive
0	0%	30	30	0
150	-50%	30	29	1
225	-25%	30	24	6
300	Cut-off	30	21	9
375	+25%	30	2	28
450	+50%	30	0	30

Analytical Specificity

The following table lists compounds that are positively detected in urine by the "DIAQUICK" MTD Dipstick (urine) at 5 minutes.

Compound	Concentration (ng/ml)
Methadone	300
Doxylamine	50,000

Precision

A study was conducted at 3 physicians' offices by untrained operators using 3 different lots of product to demonstrate the within run, between run and between operator precision. An identical panel of coded specimens containing no Methadone, 25% Methadone above and below the cut-off and 50% Methadone above and below the 300 ng/mL cut-off was provided to each site. The following results were tabulated:

Methadone conc. (ng/mL)	n	Site A		Site B		Site C	
		-	+	-	+	-	+
0	15	15	0	15	0	15	0
150	15	12	3	15	0	15	0
225	15	8	7	14	1	15	0
375	15	0	15	0	15	1	14
450	15	1	14	0	15	0	15

Effect of Urinary Specific Gravity

Fifteen (15) urine specimens with specific gravity ranging from 1.001 to 1.032 were spiked with 150 ng/mL and 450 ng/mL of Methadone respectively. The "DIAQUICK" MTD Dipstick (urine) was tested in duplicate using the fifteen neat and spiked urine specimens. The results demonstrate that varying ranges of urinary specific gravity do not affect the test results.

Effect of Urinary pH

The pH of an aliquoted negative urine pool was adjusted to a pH range of 5 to 9 in 1 pH unit increments and spiked with Methadone to 150 ng/mL and 450 ng/mL. The spiked, pH-adjusted urine was tested with the "DIAQUICK" Methadone Test Strip (Urine) in duplicate. The results demonstrate that varying ranges of pH do not interfere with the performance of the test.

Cross-Reactivity

A study was conducted to determine the cross-reactivity of the test with compounds in either drug-free urine or Methadone-positive urine. The following compounds show no cross-reactivity when tested with the "DIAQUICK" MTD Dipstick (urine) at a concentration of 100 µg/ml.

Non Cross-Reacting Compounds

Acetaminophen	Estrone-3-sulfate	Oxymetazoline
Acetophenetidin	Ethyl-p-aminobenzoate	Papaverine
N-Acetylprocainamide	Fenoprofen	Penicillin-G
Acetylsalicylic acid	Furosemide	Pentazocine
Aminopyrine	Genistic acid	Pentobarbital
Amtryptiline	Hemoglobin	Perphenazine
Amobarbital	Hydralazine	Phencyclidine
Amoxicillin	Hydrochlorothiazide	Phenelzine
Ampicillin	Hydrocodone	Phenobarbital
L-Ascorbic acid	Hydrocortisone	Phentermine
DL-Amphetamine sulfate	O-Hydroxyhippuric acid	Trans-2-phenylcyclopropylamine hydrochloride
Acopropiline	p-Hydroxyamphetamine	L-Phenylephrine
Aspartame	p-Hydroxy-methamphetamine	β-Phenylethylamine
Atropine	3-Hydroxytyramine	Phenylpropanolamine
Benzilic acid	Ibuprofen	Prednisolone
Benzoic acid	Imipramine	Prednisone
Benzoylcegonine	Iproniazid	Procaine
Benzphetamine	(±) - Isoproterenol	Promazine
Billirubin	(±) - Brompheniramine	Promethazine
Caffeine	Ketamine	DL-Propranolol
Cannabidiol	Ketoprofen	D-Propoxyphene
Cannabinol	Labeltalol	D-Pseudoephedrine
Chloralhydrate	Lorphanol	Quinacrine
Chloramphenicol	Loperamide	Quinidine
Chlorothiazide	Mephentermine	Quinine
(±) - Chlorpheniramine	Isosuxiprine	Ranitidine
Chlorpromazine	Meperidine	Salicylic acid
Chlorquine	Meprobamate	Secobarbital
Cholesterol	Methamphetamine	Serotonin
Clomipramine	Methoxyphenamine	Sulfamethazine
Clonidine	(±) - 3,4-Methylenedioxy-amphetamine hydrochloride	Sulindac
Cocaehtylene	(±) - 3,4-Methylenedioxy-meth-amphetamine hydrochloride	Temazepam
Cocaine hydrochloride	Codeine	Tetracycline
Cortisone	Morphine-3-β-D glucuronide	Tetrahydrocortisone, 3-acetate
(-) Cotinine	Morphine Sulfate	Tetrahydrocortisone 3-(β-D-glucuronide)
Creatinine	Nalidixic acid	Tetrahydrozoline
Deoxycorticosterone	Naloxone	Thebaine
Dextromethorphan	Naltrexone	Thiamine
Diazepam	Nasroxen	Thioridazine
Diclofenac	Niacinamide	DL-Tyrosine
Diffunisal	Nifedipine	Tolbutamide
Digoxin	Norcocaine	Triamterene
Diphenhydramine	Norethindrone	Trifluoperazine
EDDP	D-Norpropoxyphene	Trimethoprim
EMDP	Noscapine	Trimipramine
Ecgonine hydrochloride	DL-Octopamine	Trypamine
Ecgonine methylester	Oxalic acid	DL-Tryptophan
(-) -ψ - Ephedrine	Oxazepam	Tyramine
[1R,2S] (-) Ephedrine	Oxolinic acid	Uric acid
(L) - Epinephrine	Oxycodone	Verapamil
Erythromycin		Zomepirac
β-Estradiol		

REFERENCES

- Glass, IB. The International Handbook of Addiction Behavior. Routledge Publishing, New York, NY. 1991; 216
- Baselt RC. Disposition of Toxic Drugs and Chemicals in Man. 2nd Ed. Biomedical Publ., Davis, CA. 1982; 488
- Hawks RL, CN Chiang. Urine Testing for Drugs of Abuse. National Institute for Drug Abuse (NIDA), Research Monograph 73, 1986

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