“DIAQUICK” CEA Cassette (Carcinoembryonic Antigen) for whole blood, serum or plasma samples

REF

206030

Content

- 30 cassettes individually packed + disposable pipette (30x REF Z03030B)
- 2 vials buffer, sufficient for 30 tests
- 1 package insert

For in vitro diagnostic use only

INTENDED USE

The „DIAQUICK“ CEA Cassette is a rapid, chromatographic immunoassay for the qualitative determination of human carcinoembryonic antigen (CEA) in human whole blood, serum or plasma to aid in monitoring of cancer patients.

CLINICAL SIGNIFICANCE

Carcinoembryonic Antigen (CEA) is a tumour-associated antigen characterized as an oncofetal glycoprotein.¹ CEA is expressed in a variety of malignancies, particularly pulmonary or gastrointestinal tumours (e.g. colon cancer, liver cancer and lung cancer). CEA normally occurs in foetal gut tissue with detectable serum levels essentially disappearing after birth.^{2, 3} Therefore, elevated levels of CEA can be of significant value in the diagnosis of primary carcinomas. In addition to qualitative assessment, CEA testing plays an important role in the monitoring of cancer patients. Clinical evidence indicates that CEA levels can serve as predictive markers in both pre- and post-treatment cancer. Progressive elevation of CEA may signal tumour recurrence 3 – 36 months before clinical evidence of metastasis. Persistent elevation of circulating CEA following treatment is strongly indicative of occult metastatic and residual diseases and deficient therapeutic response.⁴

The “DIAQUICK” CEA Cassette utilizes a combination of anti-CEA antibody coated particles and anti-CEA antibodies to detect elevated levels of CEA in whole blood, serum or plasma. The minimum detection level is 5 ng/mL.

TEST PRINCIPLE

The „DIAQUICK“ CEA Cassette (whole blood/serum/plasma) is a qualitative membrane based immunoassay for the detection of CEA in whole blood, serum or plasma. The membrane is pre-coated with anti-CEA antibodies in the test line region. During testing, the specimen reacts with the particle coated with anti-CEA antibodies. The mixture migrates upward on the membrane chromatographically by capillary action to react with anti-CEA antibodies on the membrane and generate a coloured line. The presence of this coloured line in the test line region indicates a positive result, while its absence indicates a negative result. To serve as a procedural control, a coloured line will always appear in the control line region, indicating that proper volume of specimen has been added and membrane wicking has occurred.

REAGENTS

The test cassette contains anti-CEA antibody coated particles and anti-CEA antibody coated on the membrane.

MATERIALS PROVIDED

- test cassettes
- package insert
- dropper
- buffer

MATERIAL REQUIRED BUT NOT PROVIDED

- specimen collection container
- timer
- centrifuge (for plasma only)

For *fingerstick whole blood*:

- lancets
- heparinized capillary tubes and dispensing bulb

WARNINGS AND PRECAUTIONS

- For in vitro diagnostic use only. Do not use kit beyond the expiration date.
- Do not eat, drink or smoke in the area where the specimens or kits are handled.
- Do not use test if pouch is damaged.
- Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards throughout all procedures and follow the standard procedures for proper disposal of specimens.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are assayed.
- The used test should be discarded according to local regulations
- Humidity and temperature can adversely affect results.

STORAGE

Store as packaged in the sealed pouch either at room temperature or refrigerated (2 – 30 °C). The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. **DO NOT FREEZE.** Do not use after the expiration date.

SAMPLE COLLECTION AND PREPARATION

- The “DIAQUICK” CEA Cassette (whole blood/serum/plasma) can be performed using whole blood (from venipuncture or fingerstick), serum, or plasma.
- To collect **Venipuncture Whole Blood specimens**: Collect anti-coagulated blood specimen (sodium or lithium heparin, potassium or sodium EDTA, sodium oxalate, sodium citrate) following standard laboratory procedures.
- To collect **Fingerstick Whole Blood specimens**:
 - Wash the patient's hand with soap and warm water or clean with an alcohol swab. Allow to dry.
 - Massage the hand without touching the puncture site by rubbing down the hand towards the fingertip of the middle or ring finger.
 - Puncture the skin with a sterile lancet. Wipe away the first sign of blood.

- Gently rub the hand from wrist to palm to finger to form a rounded drop of blood over the puncture site.
- Add the Fingerstick Whole Blood specimen to the test by using **a capillary tube**:
 - Touch the end of the capillary tube to the blood until filled to approx. 50 µL. Avoid air bubbles.
 - Place the bulb onto the top end of the capillary tube, then squeeze the bulb to dispense the whole blood to the specimen well (S) of the test cassette.
- Add the Fingerstick Whole Blood specimen to the test by using **hanging drops**:
 - Position the patient's finger so that the drop of blood is just above the specimen well (S) of the test cassette.
 - Allow 2 hanging drops of fingerstick whole blood to fall into the specimen well (S) of the test cassette, or move the patient's finger so that the hanging drop touches the specimen well (S). Avoid touching the finger directly to the specimen well (S).
- Separate serum or plasma from blood as soon as possible to avoid haemolysis. Use only clear, non-haemolysed specimens.
- Testing should be performed immediately after specimen collection. Do not leave the specimens at room temperature for prolonged periods. Serum and plasma specimens may be stored at 2 – 8 °C for up to 3 days. For long term storage, specimens should be kept below -20 °C. Whole blood collected by venipuncture should be stored at 2 – 8 °C if the test is to be run within 2 days of collection. Do not freeze whole blood specimens. Whole blood collected by fingerstick should be tested immediately.
- Bring specimens to room temperature prior to testing. Frozen specimens must be completely thawed and mixed well prior to testing. Specimens should not be frozen and thawed repeatedly.
- If specimens are to be shipped, they should be packed in compliance with local regulations covering the transportation of etiologic agents.

ASSAY PROCEDURE

Allow the test, specimen, buffer, and/or controls to reach room temperature (15 – 30 °C) prior to testing.


1. Bring the pouch to room temperature before opening it. Remove the test cassette from the sealed pouch and use it as soon as possible.
2. Place the test cassette on a clean and level surface.

For **Serum or Plasma specimens**: Hold the dropper vertically and transfer **1 drop of serum or plasma** (approx. 25 µL) to the specimen well (S) of the test cassette, then **add 1 drop of buffer** (approx. 40 µL) and start the timer.

For **Venipuncture Whole Blood specimens**: Hold the dropper vertically and transfer **2 drops of whole blood** (approx. 50 µL) to the specimen well (S) of the test cassette, then **add 1 drop of buffer** (approx. 40 µL) and start the timer.

For **Fingerstick Whole Blood specimens**:

- To use a capillary tube: Fill the capillary tube and transfer **approx. 50 µL of fingerstick whole blood** to the specimen well (S) of the test cassette, then **add 1 drop of buffer** (approx. 40 µL) and start the timer. See illustration below.
 - To use hanging drops: **Allow 2 hanging drops of fingerstick whole blood** (approx. 50 µL) to fall into the centre of the specimen well (S) on the test cassette, then **add 1 drop of buffer** (approx. 40 µL) and start the timer. See illustration below.
3. Wait for the coloured line(s) to appear. **Read results at 5 minutes.** Do not interpret the result after 20 minutes.

INTERPRETATION OF RESULTS

POSITIVE: Two distinctly coloured lines appear. One line should be in the control line region (C) and another apparent line should be in the test line region (T).

***NOTE:** The intensity of the colour in the test line region (T) will vary depending on the concentration of CEA antibody present in the specimen. Therefore, any shade of colour in the test region (T) should be considered positive.

NEGATIVE: One coloured line appears in the control line region (C). No line appears in the test line region (T).

INVALID: Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

EXPECTED VALUES

The “DIAQUICK” CEA Cassette (whole blood/serum/plasma) has been compared with a leading commercial CEA EIA test. The correlation between these two systems is over 99.1 %.

QUALITY CONTROL

A procedural control is included in the test. A coloured line appearing in the control line region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique. Control standards are not supplied with this kit; however, it is recommended that positive and negative controls be tested as a good laboratory practice to confirm the test procedure and to verify proper test performance.

LIMITATIONS OF THE PROCEDURE

1. The "DIAQUICK" CEA Cassette is for in vitro diagnostic use only. The test should be used for the detection of CEA in whole blood, serum or plasma specimens only. Neither the quantitative value nor the rate of increase in CEA concentration can be determined by this qualitative test.
2. The "DIAQUICK" CEA Cassette will only indicate the presence of CEA in the specimen and should not be used as the sole criteria for the diagnosis of gastrointestinal tract tumours or other cancer.
3. The "DIAQUICK" CEA Cassette cannot detect less than 5 ng/mL of CEA in specimens.
4. As with all diagnostic tests, all results must be interpreted together with other clinical information available to the physician.
5. If the test result is negative and clinical symptoms persist, additional testing using other clinical methods is recommended. A negative result does not at any time preclude the possibility of gastrointestinal tract tumours or other cancer.

PERFORMANCE CHARACTERISTICS

Clinical Sensitivity, Specificity and Accuracy

The „DIAQUICK“ CEA Cassette has correctly identified a panel of specimens and has been compared to a leading commercial CEA EIA test using clinical specimens. The results show that the relative sensitivity of the „DIAQUICK“ CEA Cassette is 98.7 %, and the relative specificity is 99.3 %.

Method	Results	EIA		Total Results
		Positive	Negative	
"DIAQUICK" CEA Cassette	Positive	156	3	159
	Negative	2	424	426
Total Results		158	427	585

Relative Sensitivity: 98.7 % (95.5 % - 99.8 %)*

Relative Specificity: 99.3 % (98.0 % - 99.9 %)*

Accuracy: 99.1 % (98.0 % - 99.7 %)*

* 95% Confidence Interval

Precision

Intra-Assay

Within-run precision has been determined by using 10 replicates of three specimens: a negative, a low positive and a high positive. The negative, low positive and high positive values were correctly identified > 99 % of the time.

Inter-Assay

Between-run precision has been determined by 10 independent assays on the same three specimens: a negative, a low positive and a high positive. Three different lots of the „DIAQUICK“ CEA Cassette have been tested using negative, low positive and high positive specimens. The specimens were correctly identified > 99 % of the time.

Cross-Reactivity

Specimens positive for HCV, HBV, HIV, AFP and rheumatoid factor (RF) have been tested. No cross-reactivity was observed, indicating that the „DIAQUICK“ CEA Cassette (whole blood/serum/plasma) has a high degree of specificity for Carcinoembryonic Antigen.

Interfering Substances

The „DIAQUICK“ CEA Cassette has been tested for possible interference from visibly haemolysed and lipaemic specimens. No interference was observed. In addition, no interference was observed in specimens containing up to 2,000 mg/dL Haemoglobin, 30 mg/dL Bilirubin, 700 mg/dL Triglycerides and 1,700 mg/dL Total Lipids.

BIBLIOGRAPHY

1. Gold P and Freedman SO. Demonstration of Tumor-specific antigens in human colonic carcinomata by immunological tolerance and absorption. *J. Exp. Med.* 121:439, 1965.
2. Banjo C, Gold P, Freedman SO and Krupcy J. Immunologically Active Hetero-saccharides of the Carcinoembryonic Antigen (CEA) of the Human Digestive System, *Nature, New Biol.* 238,183, 1972.
3. Darcy DA, Turberville C, and Janes R. Immunological Study of Carcinoembryonic Antigen (CEA) and Related Glycoprotein, *Br.J.Cancer*, 28,147, 1973.
4. Phil Gold CC, and Goldenberg NA. The Carcinoembryonic Antigen (CEA): Past, Present, and Future. *Perspectives in Colon and Rectal Surgery* 9(2), 1996.

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LITERATUR

- Gold P and Freedman SO. Demonstration of Tumor-specific antigens in human colonic carcinomata by immunological tolerance and absorption. *J. Exp. Med.* 121:439, 1965.
- Banjo C, Gold P, Freedman SO and Krupcy J. Immunologically Active Heterosaccharides of the Carcinoembryonic Antigen (CEA) of the Human Digestive System, *Nature, New Biol*, 238,183, 1972.
- Darcy DA, Turberville C, and Janes R. Immunological Study of Carcinoembryonic Antigen (CEA) and Related Glycoprotein, *Br. J. Cancer*, 28,147, 1973.
- Phil Gold CC, and Goldenberg NA. The Carcinoembryonic Antigen (CEA): Past, Present, and Future. *Perspectives in Colon and Rectal Surgery* 9(2), 1996.