# "DIAQUICK" Syphilis-Cassette

for whole blood, serum and plasma samples

## REF

**Z06903CE** - 30 cassettes individually packed + disposable pipette (30 x Ref. No: Z03903B)  
 - 2 vials buffer, sufficient for 30 tests  
 - 1 package insert

## Content

### For in vitro diagnostic use only

## GENERAL INFORMATION

<b>Method</b>	sandwich type immunochromatographic assay
<b>Antigens</b>	Syphilis antigen coated particles
<b>Antibodies</b>	Anti-human IgG coated membrane
<b>Shelf life</b>	24 months from date of production
<b>Storage</b>	2 - 30°C
<b>Sample</b>	Human whole blood, serum and plasma
<b>Results</b>	after 15 to 20 minutes at room temperature

## INTENDED USE

The "DIAQUICK" Syphilis-Cassette is a rapid chromatographic immunoassay for the qualitative detection of antibodies (IgG and IgM) to *Treponema Pallidum* (TP) in whole blood, serum or plasma to aid in the diagnosis of Syphilis.

## CLINICAL SIGNIFICANCE

*Treponema Pallidum* (TP) is the causative agent of the venereal disease Syphilis. TP is a spirochete bacterium with an outer envelope and a cytoplasmic membrane. Relatively little is known about the organism in comparison with other bacterial pathogens. According to the Center for Disease Control (CDC), the number of cases of Syphilis infection has markedly increased since 1985. Some key factors that have contributed to this rise include the crack cocaine epidemic and the high incidence of prostitution among drug users. One study reported a substantial epidemiological correlation between the acquisition and transmission of the HIV virus and Syphilis.

Multiple clinical stages and long periods of latent, asymptomatic infection are characteristic of Syphilis. Primary Syphilis is defined by the presence of a chancre at the site of inoculation. The antibodies response to the TP bacterium can be detected within 4 to 7 days after the chancre appears. The infection remains detectable until the patient receives adequate treatment.

The "DIAQUICK" Syphilis-Cassette utilizes a double antigen combination of a Syphilis antigen coated particle and Syphilis antigen immobilized on membrane to detect TP antibodies (IgG and IgM) qualitatively and selectively in whole blood, serum or plasma.

## TEST PRINCIPLE

The "DIAQUICK" Syphilis-Cassette (whole blood/serum/plasma) is a qualitative membrane based immunoassay for the detection of TP antibodies in whole blood, serum or plasma specimens. In this test procedure, recombinant Syphilis antigen is immobilized in the test line region of the test. After specimen is added to the specimen well of the cassette, it reacts with Syphilis antigen coated particles in the test. This mixture migrates chromatographically along the length of the test and interacts with the immobilized Syphilis antigens. The double antigen test format can detect both IgG and IgM in specimens. If the specimen contains TP antibodies, a colored line will appear in the test line region indicating a positive result. If the specimen does not contain TP antibodies, a colored line will not appear in this region indicating a negative result. To serve as a procedural control, a colored line will always appear in the control line region, indicating that proper volume of specimen has been added and membrane wicking has occurred.

## REAGENTS

The test contains Syphilis antigen coated particles and Syphilis antigen coated on the membrane.

## WARNINGS AND PRECAUTIONS

- For in vitro diagnostic use only. Do not use kit beyond the expiration date.
- Do not eat, drink or smoke in the area where the specimens or kits are handled.
- Do not use test if pouch is damaged.
- Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards throughout all procedures and follow the standard procedures for proper disposal of specimens.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are assayed.
- The used test should be discarded according to local regulations
- Humidity and temperature can adversely affect results.

## STORAGE

Store as packaged in the sealed pouch either at room temperature or refrigerated (2 - 30 °C). The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. **DO NOT FREEZE.** Do not use after the expiration date.

## SAMPLE COLLECTION AND PREPARATION

- The "DIAQUICK" Syphilis-Cassette (whole blood/serum/plasma) can be performed using whole blood (from venipuncture or fingerstick), serum, or plasma.
- To collect **Venipuncture Whole Blood specimens:**  
 Collect anti-coagulated blood specimen (sodium or lithium heparin, potassium

or sodium EDTA, sodium oxalate, sodium citrate) following standard laboratory procedures.

- To collect **Fingerstick Whole Blood specimens:**

- Wash the patient's hand with soap and warm water or clean with an alcohol swab. Allow to dry.
- Massage the hand without touching the puncture site by rubbing down the hand towards the fingertip of the middle or ring finger.
- Puncture the skin with a sterile lancet. Wipe away the first sign of blood.
- Gently rub the hand from wrist to palm to finger to form a rounded drop of blood over the puncture site.

- Add the Fingerstick Whole Blood specimen to the test by using **a capillary tube:**

- Touch the end of the capillary tube to the blood until filled to approximately 50 µL. Avoid air bubbles.
- Place the bulb onto the top end of the capillary tube, then squeeze the bulb to dispense the whole blood to the specimen well (S) of the test cassette.

- Add the Fingerstick Whole Blood specimen to the test by using **hanging drops:**

- Position the patient's finger so that the drop of blood is just above the specimen well (S) of the test cassette.
- Allow 2 hanging drops of fingerstick whole blood to fall into the specimen well (S) of the test cassette, or move the patient's finger so that the hanging drop touches the specimen well (S). Avoid touching the finger directly to the specimen well (S).

- Separate serum or plasma from blood as soon as possible to avoid hemolysis. Use only clear, non-hemolyzed specimens.

- Testing should be performed immediately after specimen collection. Do not leave the specimens at room temperature for prolonged periods. Serum and plasma specimens may be stored at 2-8°C for up to 3 days. For long term storage, specimens should be kept below -20°C. Whole blood collected by venipuncture should be stored at 2-8°C if the test is to be run within 2 days of collection. Do not freeze whole blood specimens. Whole blood collected by fingerstick should be tested immediately.

- Bring specimens to room temperature prior to testing. Frozen specimens must be completely thawed and mixed well prior to testing. Specimens should not be frozen and thawed repeatedly.

- If specimens are to be shipped, they should be packed in compliance with local regulations covering the transportation of etiologic agents.

## ASSAY PROCEDURE

**Allow the test cassette, specimen, buffer, and/or controls to equilibrate to room temperature (15-30°C) prior to testing.**

1. Bring the pouch to room temperature before opening it. Remove the test cassette from the sealed pouch and use it as soon as possible.
2. Place the test cassette on a clean and level surface.
3. For **Serum or Plasma specimens:** Hold the dropper vertically and transfer **3 drops of serum or plasma** (approximately 75 µL) to the specimen well (S) and start the timer. See illustration below.
4. For **Venipuncture Whole Blood** specimen: Hold the dropper vertically and transfer **2 drops of whole blood** (approximately 50 µL) to the specimen well (S), then **add 1 drop of buffer** (approximately 40 µL), and start the timer. See illustration below.
5. For **Fingerstick Whole Blood** specimen:
  - To use a capillary tube: Fill the capillary tube and **transfer approximately 50 µL of fingerstick whole blood specimen** to the specimen well (S) of the test cassette, then **add 1 drop of buffer** (approximately 40 µL) and start the timer. See illustration below.
  - To use hanging drop: Allow **2 hanging drop of fingerstick whole blood specimen** (approximately 50 µL) to fall into the center of the specimen well (S) on the test cassette, then **add 1 drop of buffer** (approximately 40 µL) and start the timer. See illustration below.
6. Wait for the red line(s) to appear. The result should be read at 10 minutes. Do not interpret results after 30 minutes!

## EXPECTED VALUES

The "DIAQUICK" Syphilis-Cassette (whole blood/serum/plasma) has been compared with a leading commercial TPHA Syphilis test, demonstrating an overall accuracy of 99.7%.

## INTERPRETATION OF RESULTS


(Please refer to the illustration below)

**POSITIVE:** \* **Two lines appear.** One colored line should be in the control line region (C) and another apparent colored line should be in the test line region (T)

\***NOTE:** The intensity of the color in the test line region (T) will vary depending on the concentration of TP antibodies present in the specimen. Therefore, any shade of color in the test line region (T) should be considered positive.

**NEGATIVE:** **One colored line appears in the control line region (C).** No line appears in the test line region (T).

**INVALID: Control line fails to appear.** Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

#### MATERIALS PROVIDED

- Test cassettes
- Droppers
- Buffer (for whole blood only)
- Package insert

#### MATERIALS REQUIRED BUT NOT PROVIDED

- Specimen collection containers
- Lancets (for fingerstick whole blood only)
- Centrifuge (for plasma only)
- Timer
- Heparinized capillary tubes and dispensing bulb (for fingerstick whole blood only)

#### QUALITY CONTROL

A procedural control is included in the test. A colored line appearing in the control line region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

Control standards are not supplied with this kit; however, it is recommended that positive and negative controls be tested as a good laboratory practice to confirm the test procedure and to verify proper test performance.

#### LIMITATIONS

- The "DIAQUICK" Syphilis-Cassette (whole blood/serum/plasma) is for in vitro diagnostic use only. The test should be used for the detection of TP antibodies in whole blood, serum or plasma specimens only. Neither the quantitative value nor the rate of increase in TP antibodies can be determined by this qualitative test.
- The "DIAQUICK" Syphilis-Cassette (whole blood/serum/plasma) will only indicate the presence of TP antibodies in the specimen and should not be used as the sole criteria for the diagnosis of TP infection.
- As with all diagnostic tests, all results must be interpreted together with other clinical information available to the physician.
- If the test result is negative and clinical symptoms persist, additional testing using other clinical methods is recommended. A negative result does not at any time preclude the possibility of TP infection.

#### PERFORMANCE CHARACTERISTICS

##### Clinical Sensitivity, Specificity and Accuracy

The "DIAQUICK" Syphilis-Cassette (whole blood/serum/plasma) has correctly identified specimens of a seroconversion panel and has been compared to a leading commercial TPHA Syphilis test using clinical specimens. The results show that the relative sensitivity of the Syphilis Rapid Test Cassette (whole blood/serum/plasma) is 99.7%, and the relative specificity is 99.6%.

##### DIAQUICK Syphilis-Cassette vs. FTA-ABS

Method	FTA-ABS		Total Results
	Positive	Negative	
"DIAQUICK" Syphilis Ab Cassette	Positive	2	386
	Negative	493	494
<b>Total Results</b>	<b>385</b>	<b>495</b>	<b>880</b>

Relative Sensitivity: 99.7% (98.6-100.0%)\*

Relative Specificity: 99.6% (98.6-100.0%)\*

Accuracy: 99.7% (99.0-99.9%)\*

\* 95% Confidence Interval

##### Precision

###### Intra-Assay

Within-run precision has been determined by testing 10 replicates of four specimens: a negative, a low positive, a medium positive and a high positive. The negative, low positive, medium positive and high positive values were correctly identified >99% of the time.

###### Inter-Assay

Between-run precision has been determined by 10 independent assays on the same four specimens: a negative, a low positive, a medium positive and a high positive. Three different lots of the "DIAQUICK" Syphilis-Cassette (whole blood/serum/plasma) have been tested over a 3-day period using negative, low positive, medium positive and high positive specimens. The specimens were correctly identified >99% of the time.

#### REFERENCES

1. Claire M. Fraser. *Complete genome sequence of Treponema Pallidum, the Syphilis spirochete*, Science (1998); 281 July: 375-381.
2. Center for Disease Control. *Recommendations for diagnosing and treating Syphilis in HIV-infected patients*, MMWR Morb. Mortal Wkly Rep. (1988); 37: 601
3. Aral R. Marx. *Crack, sex and STD*, Sexually Transmitted Diseases, 1991; 18:92-101
4. L.N. Wasserheit. *Epidemiological Synergy: Interrelationships between human immunodeficiency virus infection and other sexually transmitted diseases*, Sexually Transmitted Diseases 1992; 19:61-77
5. Phillip C. Johnson. *Testing for Syphilis*, Dermatologic Clinic (1994); 12 Jan: 9-17.

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#### LITERATUR

1. Claire M. Fraser. *Complete genome sequence of Treponema Pallidum, the Syphilis spirochete*, Science (1998); 281 July: 375-381.
2. Center for Disease Control. *Recommendations for diagnosing and treating Syphilis in HIV-infected patients*, MMWR Morb. Mortal Wkly Rep. (1988); 37: 601
3. Aral R. Marx. *Crack, sex and STD*, Sexually Transmitted Diseases, 1991; 18:92-101
4. L.N. Wasserheit. *Epidemiological Synergy: Interrelationships between human immunodeficiency virus infection and other sexually transmitted diseases*, Sexually Transmitted Diseases 1992; 19:61-77
5. 3.Phillip C. Johnson. *Testing for Syphilis*, Dermatologic Clinic (1994); 12 Jan: 9-17.

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## REFERENCIAS

1. Claire M. Fraser. *Complete genome sequence of Treponema Pallidum, the Syphilis spirochete*, Science (1998); 281 July: 375-381.
2. Center for Disease Control. *Recommendations for diagnosing and treating Syphilis in HIV-infected patients*, MMWR Morb. Mortal Wkly Rep. (1988); 37: 601
3. Aral R. Marx. *Crack, sex and STD, Sexually Transmitted Diseases*, 1991; 18:92-101
4. L.N. Wasserheit. *Epidemiological Synergy: Interrelationships between human immunodeficiency virus infection and other sexually transmitted diseases*, Sexually Transmitted Diseases 1992; 19:61-77
5. Phillip C. Johnson. *Testing for Syphilis*, Dermatologic Clinic (1994); 12 Jan: 9-17.