# "DIAQUICK" TCA Dipstick (Tricyclic Antidepressants) for urine samples

REF	Content
Z05045CE	- 30 tests individually packed (30 x Ref. No: Z05045B). - 1 package insert
Z05045B	- 1 test individually packed - 1 package insert

**For in vitro diagnostic use only.  
 For diagnosis and therapeutic monitoring only.  
 For use by medical professionals only.**

## GENERAL INFORMATION

<b>Method</b>	competitive immunochromatographic assay
<b>Shelf life</b>	24 months from date of production
<b>Storage</b>	2-30°C
<b>Sample</b>	human urine
<b>Results</b>	within 5 minutes at room temperature
<b>Sensitivity</b>	1000 ng/ml

## INTENDED USE

The "DIAQUICK" TCA Dipstick (urine) is a lateral flow chromatographic immunoassay for the detection of Nortriptyline (metabolites of Tricyclic Antidepressants) in human urine at a cut-off concentration of 1000ng/mL. This assay provides only a preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are used.

## CLINICAL SIGNIFICANCE

TCA (Tricyclic Antidepressants) are commonly used for the treatment of depressive disorders. TCA overdoses can result in profound central nervous system depression, cardiotoxicity and anticholinergic effects. TCA overdose is the most common cause of death from prescription drugs. TCAs are taken or sometimes by injection. TCAs are metabolized in the liver, both TCAs and their metabolites are excreted in urine mostly in the form of metabolites for up to ten days. The "DIAQUICK" TCA Dipstick (urine) is an easy to use urine test that can be performed without the use of an instrument. The test utilizes a monoclonal antibody to selectively detect elevated levels of Nortriptyline in urine. The "DIAQUICK" TCA Dipstick (urine) yields a positive result when the Nortriptyline in urine exceeds 1000ng/ml. For in vitro diagnostic use only.

## TEST PRINCIPLE

The "DIAQUICK" TCA Dipstick is a rapid chromatographic immunoassay based on the principle of competitive binding. Drugs which may be present in the urine specimen compete against the drug conjugate for binding sites on the antibody. During testing, a urine specimen migrates upward by capillary action. Tricyclic Antidepressants, if present in the urine specimen below 1000ng/mL, will not saturate the binding sites of the antibody coated particles in the test strip. The antibodies coated particles will then be captured by immobilized Tricyclic Antidepressants conjugate and a visible colored line will show up in the test line region. The colored line will not form in the test line region if the Tricyclic Antidepressants level exceeds 1000ng/mL because it will saturate all the binding sites of anti-Tricyclic Antidepressants antibodies. A drug-positive urine specimen will not generate a colored line in the test line region because of drug competition, while a drug-negative urine specimen will generate a line in the test line region because of the absence of drug competition. To serve as a procedural control, a colored line will always appear at the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred.

## WARNINGS AND PRECAUTIONS

- For in vitro diagnostic use only. Do not use after the expiration date.
- The test strip should remain in the sealed pouch until use.
- All specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The used test strip should be discarded according to federal, state and local regulations.

## STORAGE

Store test kit at 2-30°C; do not freeze. Refer to the expiration date for stability. Do not freeze!

## SAMPLE COLLECTION AND PREPARATION

### Urine Assay

The urine specimen must be collected in a clean and dry container. Urine collected at any time of the day may be used. Urine specimens exhibiting visible particles should be centrifuged, filtered, or allowed to settle to obtain a clear specimen for testing.

### Specimen Storage

Urine specimens may be stored at 2-8°C for up to 48 hours prior to assay. For long-term storage, specimens may be frozen and stored below -20°C. Frozen specimens should be thawed and mixed before testing.

## ASSAY PROCEDURE

Bring an unopened foil pouch and the sample specimen to room temperature.

1. Open the foil pouch by tearing along the notch and remove the test strip.
2. With arrows pointing toward the urine specimen, immerse the test strip vertically in the urine specimen for at least 10-15 seconds. Do not pass the maximum line (MAX) on the test strip when immersing the strip.
3. Place the test strip on a non-absorbent flat surface, start the timer and wait for the red line(s) to appear. The result should be read at 5 minutes. Do not interpret the result after 10 minutes.

**Important:** Do not interpret a test result after more than 10 minutes. Waiting longer than five minutes may cause inaccurate interpretation. To avoid confusion, discard the test device after reading the result.

## INTERPRETATION OF RESULTS

### Positive Result:

A colored band is visible in the control line region "C" but not in the test line region "T". A positive result indicates that the TCA concentration is at or above the sensitivity of the test.


### Negative Result:

A colored band is visible in the control line region "C", and in the test line region "T". A negative result indicates that the TCA concentration is below the sensitivity of the test.

### Invalid Result:

There is no distinct color band visible neither in the test line region "T" nor in the control line region "C". The test should be considered invalid and it is recommended that the sample be retested.

**Note:** A very faint line on the test region indicates that the TCA in the sample is near or below the cut-off level for the test. However, if any line is observed at the test region, the sample is NEGATIVE.

## QUALITY CONTROL

A procedural control is included in the test. A red line appearing in the control region (C) is considered an internal positive procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique. In addition, if the test has been performed properly, the background will be clear to provide a distinctive result.

Control standards are not supplied with this kit; however it is recommended that positive and negative controls be tested as good laboratory testing practice to confirm the test procedure and to verify proper test performance.

## LIMITATIONS

- The "DIAQUICK" TCA Dipstick provides only a preliminary analytical result. A secondary analytical method must be used to obtain a confirmed result. Gas chromatography and mass spectrometry (GC/MS) are the preferred confirmatory methods.
- The "DIAQUICK" TCA Dipstick is a qualitative urine assay and can not determine either the drug concentration in the urine or the level of intoxication.
- It is possible that technical or procedural errors, as well as other interfering substances in the urine specimen may cause erroneous results.
- Adulterants, such as bleach and/or alum, in urine specimens may produce erroneous results regardless of the analytical method used. If adulteration is suspected, the test should be repeated with another urine specimen. A Positive Result indicates presence of the drug or its metabolites but does not indicate level or intoxication, administration route or concentration in urine. A Negative Result may not necessarily indicate drug-free urine. Negative results can be obtained when drug is present but below the cut-off level of the test.

## PERFORMANCE CHARACTERISTICS

### Accuracy

A side-by-side comparison was conducted using the "DIAQUICK" TCA Dipstick and a leading commercially available TCA rapid test. Testing was performed on 300 clinical specimens. Ten percent of the specimens employed were either at -25% or +25% level of the cut-off concentration of Tricyclic Antidepressants. Presumptive positive results were confirmed by GC/MS. The following results were tabulated:

Method	Results	Other TCA Rapid Test		Total Results
		Positive	Negative	
DIAQUICK TCA Dipstick	Positive	140	1	141
	Negative	5	154	159
Total Results		145	155	300
% Agreement with this Rapid Test kit		97%	99%	98%

When compared with GC/MS at a cut-off of 1000ng/ml, the following results were tabulated:

Method	Results	GC/MS		Total Results
		Positive	Negative	
DIAQUICK TCA Dipstick	Positive	135	8	143
	Negative	4	153	157
Total Results		139	161	300
% Agreement with GC/MS Analysis		97%	95%	96%

Eighty (80) of these specimens were also run using the "DIAQUICK" TCA Dipstick by an untrained operator at a different site. Based on GC/MS data, the operator obtained a statistically similar Positive Agreement, Negative Agreement and Overall Agreement rate as the laboratory personnel.

#### Analytical Sensitivity

A drug-free urine pool was spiked with Tricyclic Antidepressants at the following concentrations: 0 ng/mL, 500 ng/mL, 750 ng/mL, 1,000 ng/mL, 1,250 ng/mL and 1,500 ng/mL. The result demonstrates >99% accuracy at 50% above and 50% below the cut-off concentration. The data are summarized below:

Nortriptyline Concentration (ng/mL)	Percent of Cut-off	n	Visual Result	
			Negative	Positive
0	0%	30	30	0
500	-50%	30	30	0
750	-25%	30	22	8
1,000	Cut-off	30	17	13
1,250	+25%	30	5	25
1,500	+50%	30	0	30

#### Analytical Specificity

The following table lists compounds that are positively detected in urine by the "DIAQUICK" TCA Dipstick at 5 minutes.

Compound	Concentration (ng/mL)
Nortriptyline	1,000
Nordoxepine	1,000
Trimipramine	3,000
Amtriptyline	1,500
Promazine	1,500
Desipramine	200
Imipramine	400
Clomipramine	12,500
Doxepine	2,000
Maprotiline	2,000
Promethazine	25,000

#### Precision

A study was conducted at 3 geographically distinct physicians' offices by untrained operators using 3 different lots of product to demonstrate the within run, between run and between operator precision. An identical panel of coded specimens containing, according to GC/MS, no Nortriptyline, 25% Nortriptyline above and below the cut-off and 50% Nortriptyline above and below the 1,000 ng/mL cut-off was provided to each site. For the specimens at -50% cut-off concentration, the 3 sites had a 98% agreement. For the -25% to +25% cut-off specimens, the 3 sites demonstrated a 75% agreement. For specimens at +50% cut-off concentration, the 3 sites demonstrated a 93% agreement. The 3 sites demonstrated an overall agreement of 89%.

#### Effect of Urinary Specific Gravity

Fifteen (15) urine specimens of normal, high, and low specific gravity ranges were spiked with 500 ng/mL and 1,500 ng/mL of Nortriptyline respectively. The "DIAQUICK" TCA Dipstick was tested in duplicate using the fifteen neat and spiked urine specimens. The results demonstrate that varying ranges of urinary specific gravity do not affect the test results.

#### Effect of Urinary pH

The pH of an aliquoted negative urine pool was adjusted to a pH range of 5 to 9 in 1 pH unit increments and spiked with Nortriptyline to 500 ng/mL and 1,500 ng/mL. The spiked, pH-adjusted urine was tested with the "DIAQUICK" TCA Dipstick in duplicate and interpreted according to the package insert. The results demonstrate that varying ranges of pH do not interfere with the performance of the test.

#### Cross-Reactivity

A study was conducted to determine the cross-reactivity of the test with compounds in a drug-free urine pool and a drug-free urine pool spiked to contain a 1,500ng/mL concentration of Nortriptyline. The following compounds show no cross-reactivity when tested with the "DIAQUICK" TCA Dipstick at a concentration of 100 µg/mL.

#### Non Cross-Reacting Compounds

Acetamidophenol	Estrone-3-sulfate	Oxolinic acid
Acetophenetidin	Ethyl-p-aminobenzoate	Oxycodone
N-Acetylprocainamide	Fenfluramine	Oxymetazoline
Acetylsalicylic acid	Fenpropafen	Papaverine
Aminopyrine	Furosemide	Penicillin-G
Amobarbital	Gentisic acid	Pentazocine
Amoxicillin	Hemoglobin	Pentobarbital
Ampicillin	Hydralazine	Perphenazine
L-Ascorbic acid	Hydrochlorothiazide	Phencyclidine
Apomorphine	Hydrocodone	Phenelzine
Aspartame	Hydrocortisone	Phenobarbital
Atropine	p-Hydroxyamphetamine	Phentermine
D,L-Amphetamine	o-Hydroxyhippuric acid	D-Phenylephrine
L-Amphetamine	3-Hydroxytyramine	β-Phenylethylamine
Benzoic acid	p-Hydroxy-methamphetamine	Phenylpropanolamine
Benzoic acid	ibuprofen	Prednisolone
Benzoylgonine	Benzphetamine	Prednisone
Benzphetamine	(±)-Isoproterenol	Procaine
Bilirubin	Isosuprine	D,L-Propranolol
(±)-Brompheniramine	Ketamine	D-Propoxyphene
Caffeine	Ketoprofen	D-Pseudoephedrine
Cannabidiol	Labelalol	Quinidine
Cannabinol	Levorphanol	Quinine
Chloralhydrate	Loperamide	Ranitidine
Chloramphenicol	Meperidine	Salicylic acid
Chloriazepoxide	Meprobamate	Secobarbital
Chlorothiazide	Methadone	Serotonin
(±) Chlorpheniramine	D-methamphetamine	(S-Hydroxytyramine)
Chlorpromazine	(L)-methamphetamine	Sulfamethazine
Chlorquine	Methoxyphenamine	Sulindac
Cholesterol	3,4-Methylenedioxyethyl-amphetamine	Temazepam
Clonidine	(±) 3,4-Methylenedioxy-methamphetamine	Tetracycline
Cocaine hydrochloride	Cocaine	Tetrahydrocortisone, 3 Acetate
Codeine	Methyphenidate	Tetrahydrocortisone 3 (β-D glucuronide)
Cortisone	Morphine-3-β-D-glucuronide	Tetrahydrozoline
(-) Cotinine	D-glucuronide	Thebaine
Creatinine	Nalidixic acid	Thiamine
Deoxycorticosterone	Naloxone	Thioridazine
Dextromethorphan	Naltrexone	Tolbutamide
Diazepam	Naproxen	Triamterene
Diclofenac	Niacinamide	Trifluoperazine
Diffunisal	Nifedipine	Trimethoprim
Digoxin	Norcodein	D, L-Tryptophan
Diphenhydramine	(-)-ψ-Ephedrine	Tyramine
Doxylamine	Egonine hydrochloride	D, L-Tyrosine
Egonine hydrochloride	Egonine methylester	Uric acid
Egonine methylester	(IR,2S)-(-)-Ephedrine	Verapamil
(IR,2S)-(-)-Ephedrine	L-Ephedrine	Oxazepam
L-Ephedrine	Erythromycin	β-Estradiol
Erythromycin	Oxalic acid	Zomepirac
β-Estradiol		

#### REFERENCES

- Rose, J.B., Tricyclic antidepressants toxicity. J. Toxicity Clin. Toxicol., 11,381-402,1977.
- Urine Testing for Drugs of Abuse, National Institute on Drug-Abuse (NIDA), Research Monograph 73, 1986.
- Fed. Register, Department of Health and Human Services, Mandatory Guidelines for Federal Workplace Drug Testing Programs, 53,69,11979,1988.
- McBay, A. J.Clin.Chem.,33,33B-40B,1987.

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Acetamidophenol	Estro-ne-3-sulfate	Oxolinic acid
Acetophenetidin	Ethyl-p-aminobenzoate	Oxycodone
N-Acetylvicprocinamide	Fenfluramine	Oxymetazoline
Acetylsalicylic acid	Fenpropfen	Papaverine
Aminopyrine	Furosemide	Penicillin-G
Amobarbital	Gentisic acid	Pentazocine
Amoxicillin	Hemoglobin	Phenobarbital
Ampicillin	Hydralazine	Perphenazine
L-Ascorbic acid	Hydrochlorothiazide	Phencyclidine
Apomorphine	Hydrocodone	Phenelzine
Aspartame	Hydrocortisone	Phenobarbital
Atropine	p-Hydroxyamphetamine	Phentermine
D,L -Amphetamine	O-Hydroxyhippionic acid	L-Phenylephrine
L-Amphetamine	3-Hydroxytyramine	β-Phenylethylamine
Benzoic acid	p-Hydroxy-	Phenylpropanolamine
Benzoic acid	methamphetamine	Prednisolone
Benzoylecgonine	Ibuprofen	Prednisone
Benzphetamine	(±)-Isoproterenol	Procaine
Bilirubin	Isoxsuprine	D,L-Propranolol
(±)-Brompheniramine	Ketamine	D-Propoxyphene
Caffeine	Ketoprofen	D-Pseudoephedrine
Cannabidiol	Labetalol	Quinidine
Cannabiol	Levorphanol	Quinine
Chloralhydrate	Loperamide	Ranitidine
Chloramphenicol	Meperidine	Salicylic acid
Chlordiazepoxide	Meprobamate	Secobarbital
Chlorothiazide	Methodone	Serotonin
(±) Chlorpheniramine	D-methamphetamine	(S-Hydroxytyramine)
Chlorpromazine	(L)-methamphetamine	Sulfamethazine
Chlorquine	Methoxyphenamine	Sulindac
Cholesterol	3,4-Methylenedioxyethyl-	Temazepam
Clonidine	amphetamine	Tetracycline
Cocaine hydrochloride	(±) 3,4-Methylenedioxy-	Tetrahydrocortisone,
Codeine	-methamphetamine	3 Acetate
Cortisone	Methylphenidate	Tetrahydrocortisone 3 (β-D-
(+) Cotinine	Morphine-3-β-	glucuronide)
Creatinine	D-glucuronide	Tetrahydrozoline
Deoxycorticosterone	Nalidixic acid	Thebaine
Dextromethorphan	Naloxone	Thiamine
Diazepam	Naltrexone	Thioridazine
Diclofenac	Naproxen	Tolbutamide
Diflunisal	Niacinamide	Triamterene
Digoxin	Nifedipine	Trifluoperazine
Diphenhydramine	Norcocaine	Trimethoprim
Doxylamine	Doxylamine	D, L-Tryptophan
Egonine hydrochloride	Norethindrone	Tyramine
Egonine methylester	D-Norpropoxyphene	D, L-Tyrosine
(IR,2S)-(-)-Ephedrine	Noscapine	Uric acid
L-Ephedrine	D,L-Octopamine	Verapamil
Erythromycin	Oxalic acid	Oxazepam
	β-Estradiol	Zomepirac

#### BIBLIOGRAPHIE

- Rose, J.B., Tricyclic antidepressants toxicity. J. Toxicity Clin. Toxicol., 11,381-402, 1977.
- Urine Testing for Drugs of Abuse, National Institute on Drug-Abuse (NIDA), Research Monograph 73, 1986.
- Fed. Register, Department of Health and Human Services, Mandatory Guidelines for Federal Workplace Drug Testing Programs, 53,69,11979, 1988.
- McBay, A. J. Clin. Chem., 33,33B-40B, 1987.

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