“DIAQUICK”

Multi-6 High Sensitivity Touch Test

for surfaces, solid or liquid samples

DOA Multi-6 HS Touch Test AMP, MET, COC, OPI, PCP, THC
 - REF Z05585 Content: - 25 tests, individually packed (25x REF Z09550B)

DOA Multi-6/1 HS Touch Test AMP, MET, COC, OPI, MTD, THC
 - REF Z12586 Content: - 25 tests, individually packed (25x REF Z09553B)

Additionally packed in every kit:
 - 25 collectors
 - 25 collection tubes, filled with Touch Test Buffer
 - 1 package insert
 - 25 security seals
 - 5 plastic pipettes (sterile)

Additionally offered:
 - REF Z09560 THC Extraction Solution, 1x 10 mL
 - REF Z12580 THC Extraction Kit (10 mL Extr. Solution + 60 Extr. Vials)
 - REF Z09570 Skin and Sweat Solvent, 1x 50 mL

For in vitro diagnostic use only. For professional use only.

INTENDED USE

The “DIAQUICK” MULTI-6 HS TOUCH TEST is a lateral flow immunochromatographic assay for the qualitative detection of traces of drugs of abuse (amphetamines, cocaine, marijuana, methamphetamines, opiates and phencyclidine). This test is intended to detect traces of drugs of abuse on surfaces, in solids (e.g. tablets, dusts and powders), in liquids (e.g. urine), on skin and in sweat. The assay provides only a preliminary analytical test result. A more specific alternative method must be used in order to obtain a confirmed analytical result, particularly when preliminary positive test results are indicated. A combination of gas chromatography and mass spectrometry is the established method.

Parameter	Calibrator	Cut-off (ng/mL)
Amphetamine (AMP)	D-Amphetamine	50
Cocaine (COC)	Benzoylcegonine	20
Marihuana (THC)	11-nor- Δ^9 -THC-9 COOH	12
Methamphetamine (MET)	D-Methamphetamine	50
Opiate (OPI)	Morphine	40
Methadone (MTD)	Methadone	30
Phencyclidine (PCP)	Phencyclidine	10

GENERAL INFORMATION

Test Cassette:

Method competitive immunochromatographic assay
Antigen > 2 μ g/test drug-protein-conjugate
Antibody > 2 μ g/test anti-drug monoclonal antibody (mouse)
Shelf life 24 months from date of production
Storage 2-30°C
Results after 10 minutes, do not read after 20 minutes


Touch Test Buffer:

Components bis-tris buffer solvent (pH 7.0-7.2), ionic tenside, organic solvent
Shelf Life 24 months from date of production
Storage 2-30°C

TEST PRINCIPLE

The “DIAQUICK” MULTI-6 HS TOUCH TEST is a chromatographic immunoassay based on the principle of competitive binding. Drugs that might be contained in samples compete against immobilized drug conjugate for a limited amount of binding sites on their respective drug specific antibodies. During testing, the specimen migrates upward to the test zone by capillary action. Dissolved drugs, if present in the tested solution below the cut-off concentration, will not saturate the binding sites of the specific antibodies. In this case the antibodies will react with the drug-protein conjugate and a visible colored line will show up in the test region. This indicates a negative result. Drugs, if present in the specimen above the cut-off concentration, will saturate the binding sites of the respective antibodies. Therefore a colored line will not show up in the test region. This indicates a positive test result. A drug-positive specimen will not generate a colored line in the specific test line region, because of drug competition, while a drug-negative specimen will generate a test line in the test line region, because of the absence of drug competition. To serve as a procedural control a colored line will always appear at the control region, indicating that proper volumes of specimen have been added and membrane wicking has occurred.

TEST SETUP

STORAGE AND STABILITY

Store as packaged in the sealed pouch either at room temperature or refrigerated (2-30°C). The test cassette is stable through the expiration date printed on the sealed pouch. The test cassette must remain in the sealed pouch until use. DO NOT FREEZE. Do not use beyond the expiration date.

WARNINGS AND PRECAUTIONS


- This test is developed for the qualitative detection of drugs of abuse and their metabolites only.
- The kit does not contain any dangerous materials. Still, the test should be used by professionals only. Keep out of children's range.
- Do not use the kit after the expiration date. The test cassette should remain in the sealed pouch until use.
- To avoid any cross contaminations do not reuse any test components.
- The used test components should be discarded according to local regulations.
- The test cassette uses a high sensitivity detection system. Use the test by applying only small amounts of sample.
- Samples may be infectious and dangerous. Please use protective clothing, safety goggles and gloves in order to avoid contamination, carryover and hazards for your health.

SAMPLE COLLECTION AND PREPARATION

If the collection tubes or the test cassettes are stored at low temperature, allow all components used for testing to reach room temperature (15-30 °C) before performing the test. Remove the test cassette from its protective pouch just before performing the test.


FOR POWDER AND SURFACE TESTING, SECURING OF EVIDENCE (EXCEPT THC)

- Remove a collector from the sealed pouch and open the cap of a collection tube.
- Moisten the sponge of the collector with Touch Test Buffer inside the collection tube until the sponge is thoroughly wet (fig. 1) and close the collection tube again.
- The sponge of the collector is dipped into the substance to be tested, or wiped on any suspicious surface (fig. 2). Even the smallest quantities of sample are sufficient for performing the test.
- Open the cap of a collection tube. Press the sponge against the strainer in order to mix the sample with the Touch Test Buffer in the collection tube (figure 3). Snap the cap of the collection tube and mix by shaking vigorously.

FOR SKIN AND SWEAT TESTING (EXCEPT THC)

- Remove a collector from the sealed pouch and apply ~20 drops of “DIAQUICK” Skin and Sweat Solvent on the sponge of the collector (fig. 1).
- Use the sponge of the collector by wiping on skin or absorbing sweat (fig. 2).
- Open the cap of a collection tube. Press the sponge against the strainer in order to mix the sample with the Touch Test Buffer in the collection tube (figure 3). Snap the cap of the collection tube and mix by shaking vigorously.


FOR TESTING SOLID SUBSTANCES (EXCEPT THC)

- Open the cap of a collection tube and transfer a little amount of any suspicious substance (tablets or solids) into the collection tube (fig. 1). Snap the cap of the collection tube and mix by shaking vigorously.


FOR TESTING OF LIQUIDS

- Use a pipette to soak up some drops of the liquid to be tested.
- Open the cap of a collection tube and transfer 3 drops of the sample to be tested into the collection tube (fig. 2). Snap the cap of the collection tube and mix by shaking vigorously. As a result of the very high sensitivity of the test already smallest amounts of the substance to be tested are sufficient to perform the test.

FOR THC DETECTION: MARIHUANA AND HASHISH


- Take 5 flower buds from one Cannabis inflorescence or a small amount of Hashish (in the mg-range) and put them into a suitable container (e.g. Eppendorf tube).
- Add 4 drops of “DIAQUICK” THC Extraction Solution (fig. 1).
- Let it extract for 1 min. (fig. 2).
- Open the cap of a collection tube, add the whole content of the extraction tube (THC Extraction Solution + THC material) and mix. Snap the cap of the collection tube and mix by shaking vigorously to mix the contents with the touch test buffer (fig. 3).


ASSAY PROCEDURE

- Remove a test cassette from the sealed pouch and place it on a clean and level surface. Open the screw cap from the collection tube and transfer 3 drops (approximately 100 µl) into each specimen well of the test cassette (fig. 1) and start the timer. Avoid trapping air bubbles in the specimen well. Wait for the colored lines to appear. Read test results at 5 minutes. Do not interpret test results after more than 20 minutes.
- Screw the cap on the collection tube, secure the collection tube with a security seal and send it to a laboratory for confirmation if necessary (fig. 2).

INTERPRETATION OF RESULTS

NEGATIVE: A colored line in the control line region (C) and a colored line in the test line region (T) for a specific drug indicate a negative result. This indicates that the drug concentration in the oral fluid specimen is below the designated cut-off level for that specific drug.

***NOTE:** The shade of color in the test line region (T) may vary, but it should be considered negative whenever there is even a faint colored line.

POSITIVE: A colored line in the control line region (C) but no line in the test line region (T) for a specific drug indicates a positive result. This indicates that the drug concentration in the oral fluid specimen exceeds the designated cut-off for that specific drug.

INVALID: Control line (C) fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test using a new test cassette. If the problem persists, discontinue using the lot immediately and contact your local distributor.

Note: A faint line appearing in the test region (T) may be caused when the drug concentration is near the cut-off value. Nevertheless the appearance of a line is the indication for a negative test. To get a confirmed result a more sensitive method must be used.

QUALITY CONTROL

A procedural control is included in the test. A colored line appearing in the control region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique. Control standards are not supplied with this kit; however, it is recommended that positive and negative controls be tested as a good laboratory practice to confirm the test procedure and to verify proper test performance.

LIMITATIONS

- The "DIAQUICK" MULTI-6 HS TOUCH TEST provides only a qualitative, preliminary analytical result. A secondary analytical method must be used to obtain a confirmed result. Gas chromatography/mass spectrometry (GC/MS) or gas chromatography/tandem mass spectrometry (GC/MS/MS) is the preferred confirmatory method.
- A positive test result does not indicate the concentration of drug in the specimen or the route of administration.
- A negative result may not necessarily indicate a drug-free specimen. Drug may be present in the specimen below the cut-off level of the test.

PERFORMANCE CHARACTERISTICS

ANALYTICAL SENSITIVITY

A phosphate-buffered saline (PBS) pool was spiked with drugs to target concentrations of ± 50% cut-off and ± 25% cut-off and tested with the "DIAQUICK" MULTI-6 HS TOUCH TEST. The results are summarized below.

Drug Conc.	n	MTD		COC		MET	
		-	+	-	+	-	+
0% Cut-off	30	30	0	30	0	30	0
-50% Cut-off	30	30	0	30	0	30	0
-25% Cut-off	30	24	6	30	0	28	2
Cut-off	30	14	16	20	10	23	7
+25% Cut-off	30	8	22	6	24	7	23
+50% Cut-off	30	0	30	0	30	0	30

Drug Conc.	n	PCP		THC		OPI	
		-	+	-	+	-	+
0% Cut-off	30	30	0	30	0	30	0
-50% Cut-off	30	30	0	30	0	30	0
-25% Cut-off	30	30	0	24	6	26	4
Cut-off	30	22	8	15	15	20	10
+25% Cut-off	30	8	22	11	19	5	25
+50% Cut-off	30	0	30	0	30	0	30

Drug Conc.	n	AMP	
		-	+
0% Cut-off	30	30	0
-50% Cut-off	30	30	0
-25% Cut-off	30	26	4
Cut-off	30	19	11
+25% Cut-off	30	7	23
+50% Cut-off	30	0	30

ANALYTICAL SPECIFICITY

The following table lists the concentration of compounds (ng/mL) above which the "DIAQUICK" DOA SALIVA MULTI-6 TEST (Oral Fluid) identified positive results at 10 minutes.

AMPHETAMINE (AMP)		OPIATE (OPI)	
D-Amphetamine	50	Morphine	40
DL-Amphetamine	125	Codeine	10
β-Phenylethylamine	4,000	Ethylmorphine	24
Tryptamine	1,500	Hydromorphone	100
p-Hydroxyamphetamine	800	Hydrocodone	100
(+)-3,4-Methylenedioxyamphetamine (MDA)	150	Levorphanol	400
L-Amphetamine	4,000	Oxycodone	25,000
COCAINE (COC)		MORPHINE 3-β-D- Glucuronide	
Benzoylcegonine	20	Norcodeine	1,500
Cocaine HCl	20	Normorphine	12,500
Cocaeethylene	25	Nalorphine	10,000
Ecgonine HCl	1,500	Oxymorphone	25,000
Ecgonine methylester	12,500	Thebaine	1,500
MARIJUANA (THC)		DIACETYLMORPHINE (Heroin)	
11-nor-Δ ⁹ -THC-9 COOH	12	6-Monoacetylmorphine	25
Cannabinol	12,500	Bilirubin	3,500
11-nor-Δ ⁹ -THC-9 COOH	2	PHENCYCLIDINE (PCP)	
Δ ⁹ -THC	6,000	Phencyclidine	10
Δ ⁹ -THC	10,000	Tetrahydrozoline	50,000
METHAMPHETAMINE (MET)		METHADONE (MTD)	
D-Methamphetamine	50	Methadone	30
Fenfluramine	60,000	Doxylamine	50,000
p-Hydroxymethamphetamine	400	Estrone-3-Sulfate	50,000
Methoxyphenamine	25,000	Phencyclidine	50,000
3,4-Methylenedioxymethamphetamine (MDMA)	50		
L-Phenylephrine	4,000		
Procaine	2,000		
(1R,2S) - (-) Ephedrine	400		

CROSS-REACTIVITY


A study was conducted to determine the cross-reactivity of the test with compounds spiked into drug-free PBS stock. The following compounds demonstrated no false positive results on the "DIAQUICK" DOA SALIVA MULTI-6 TEST (Oral Fluid) when tested at concentrations up to 100 µg/mL.

Non Cross-Reacting Compounds:

Acetaminophen	Diazepam	MDE	Promazine
Acetophenetidine	Diclofenac	Mehentermine	Promethazine
N-Acetylprocainamide	Dicyclomine	Meperidine	D/L-Propranolol
Acetylsalicylic acid	Diflunisal	Meprobamate	D-Propoxyphene
Aminopyrine	Digoxin	Methadone	D-Pseudoephedrine
Amoxicillin	Diphenhydramine	Methyphenidate	Quinacrine
Ampicillin	Doxylamine	Nalidixic acid	Quinine
Amitypyline	L-ψ-Ephedrine	Naloxone	Quindine
Amobarbital	β-Estradiol	Naltrexone	Ranitidine
Ascorbic acid	Estrone-3-sulfate	Naproxen	Salicylic acid
Apomorphine	Ethyl-p-aminobenzoate	Niacinamide	Secobarbital
Aspartame	L-Epinephrine	Nifedipine	Sulfamethazine
Atropine	Erythromycin	Nimesulide	Sulindac
Benzilic acid	Fenoprofen	Norethindrone	Temazepam
Benzoic acid	Furosemide	D-Norpropoxyphene	Tetracycline
Benzphetamine	Genitic acid	Noscapine	Tetrahydrocortisone
Bupirone	Hemoglobin	D/L-Octopamine	3-acetate
(±)-Brompheniramine	Hydralazine	Oxalic acid	Tetrahydrocortisone
Caffeine	Hydrochlorothiazide	Oxazepam	3 (β-D-glucuronide)
Cannabidiol	Hydrocortisone	Oxolinic acid	Theophylline
Chlordiazepoxide	O-Hydroxyhippuric acid	Oxymetazoline	Thiamine
Chloralhydrate	β-Hydroxynorephedrine	Papaverine	Thioridazine
Chloramphenicol	5-Hydroxytyramine	Penicillin-G	D/L-Tyrosine
Chlorothiazide	(serotonin)	Pentazocine	Tolbutamide
D/L-Chlorpheniramine	3-Hydroxytyramine	hydrochloride	Trazodone
Chlorpromazine	Ibuprofen	Pentobarbital	Triamterene
Chloroquine	Imipramine	Perphenazine	Trifluoperazine
Cholesterol	Iproniazid	Phenelzine	Trimethoprim
Clonidine	(-)Isoproterenol	Trans-2-phenylcyclopropylamine	Trimipramine
Cortisone	Isosuprine	Phentermine	D/L-Tryptophan
L-Cotinine	Ketamine	Phenylpropanolamine	Tyramine
Creatinine	Ketoprofen	Prednisolone	Uric acid
Clomipramine	Labetalol	Phenolbarbital	Verapamil
Deoxycorticosterone	Loperamide	Prednisone	Zomepirac
Dextromethorphan	Maprotiline		

REFERENCES

- Moolchan E, et al. Saliva and Plasma Testing for Drugs of Abuse: Comparison of the Disposition and Pharmacological Effects of Cocaine. Addiction Research Center, IRP, NIDA, NIH, Baltimore, MD. As presented at the SOFT-TIAFT meeting October 1998.
- Schramm W., et al. Drugs of Abuse in Saliva: A Review. J Anal Tox, 16 (1): 1-9, 1992.
- Kim I, et al. Plasma and oral fluid pharmacokinetics and pharmacodynamics after oral codeine administration. Clin Chem, 48 (9):1486-96, 2002.
- McCarron MM, et al. Detection of Phencyclidine Usage by Radioimmunoassay of Saliva. J Anal Tox. 8 (5): 197-201, 1984.
- Kang GI and Abbott FS. Analysis of methadone and metabolites in biological fluids with gas chromatography-mass spectrometry. J Chromatogr. 231(2): 311-319. Sept 1982.

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