# “DIAQUICK” CRP DIPSTICK

## SEMI-QUANTITATIVE DETECTION OF C-REACTIVE PROTEIN

for human whole blood

REF	Content
Z08025CE	- 5 dipsticks individually packed (5 x Ref. No: Z08025B) - 5 buffer vials - 1 plastic rack (for buffer vials) - 1 package insert
Z08026CE	- 20 dipsticks individually packed (20 x Ref. No: Z08025B) - 20 buffer vials - 1 plastic rack (for buffer vials) - 1 package insert

For professional *in vitro* diagnostic use only.

### GENERAL INFORMATION

<b>Method</b>	semi-quantitative, membrane based immunoassay
<b>Shelf life</b>	24 months from date of production
<b>Storage</b>	2-30°C
<b>Sample</b>	whole blood samples
<b>Results</b>	within 5 minutes at room temperature

### INTENDED USE

The „DIAQUICK“ CRP Dipstick (Whole Blood) is a rapid chromatographic immunoassay, based on two specific antibodies against human CRP. The concentration-dependent formation of test lines enables a fast, semi-quantitative detection of CRP in whole blood samples.

### SUMMARY

The acute-phase protein CRP (C-reactive protein) is a classic marker for the diagnosis of an inflammation in the body. The synthesis of CRP happens in the liver and is also triggered by pre-inflammatory cytokines like Interleukine-6. CRP reaches the site of inflammation via plasma and actively contributes to the innate immune response.

Systemic inflammations are normally accompanied by a drastic increase of the CRP concentration in blood. Mainly bacterial infections, sterile tissue damage, e.g. after a surgery, traumas or extreme physical stress, tumors mainly in metastatic stage, some auto-immune diseases and malign systemic diseases (Hodgkin and non-Hodgkin lymphomas) can be the reason for a CRP increase.

In diagnostics CRP is often used to distinguish between viral and bacterial infections, as the amount of CRP increase is notably higher in bacterial infections. Often CRP is also determined to check for the success of a therapy, e.g. an antibiotic therapy.

### PRINCIPLE

The test strip is immersed with the tip into the sample, which is diluted with buffer solution. The sample migrates from bottom to tip along the test strip. If the sample contains CRP, it binds to the first anti-CRP antibody, which is bound to colloidal gold of red color. The red coloured CRP-antibody-gold-complex now migrates together with the sample fluid along the membrane, on which the second anti-CRP antibody is fixed in lines in different concentrations. The CRP-antibody-gold-complex is caught by the immobilized antibodies on the membrane, which results in the formation of red lines. The number of lines depends on the CRP concentration in the sample material. The more CRP is in the sample, the more red lines will appear.

As a functional control, a red control line will appear at the end of the membrane, which states a correct test procedure. The formation of the control line is independent of the CRP concentration of the sample. A failure of the red control line means that the result is invalid. The test should be repeated in this case.

### STORAGE

The “DIAQUICK” CRP Dipstick (Whole Blood) should be stored at room temperature (15-30°C). If the test is stored in a cool area, the salts in the buffer solution could precipitate. Normally this does not have any influence on the product quality.

Unopened test materials can be used until the expiration date.

### PRECAUTIONS

- For professional *in vitro* diagnostic use only.
- Only for one-time use.
- Do not use after expiration date.
- The test must remain in the sealed pouch until use.
- Store and transport the test only under the stated temperature.
- Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards throughout the procedure and follow the standard procedures for proper disposal of specimens.
- Wear protective clothing such as laboratory coats, disposable gloves or eye protection when specimens are being tested.
- The sample buffer contains small amounts of sodium azide below 0,1% as preservative agent.
- Please stick to this package insert and the stated sample volumes and reading times.

### MATERIALS REQUIRED BUT NOT PROVIDED

Sterile lancets	Timer
end-to-end capillary pipettes	

### SPECIMEN COLLECTION

#### Preparation

Make sure that before use all components have reached room temperature. Cold buffer solution or condensation liquid on the test strip can lead to invalid test results.

1. Take one tube with dilution buffer from the kit. Note the name of the patient or the ID. Open the screw cap.

#### Blood Collection

2. Use a sterile lancet to get a drop of blood from the finger tip.
3. Take a drop of blood (10 µl) by using a capillary pipette. The capillary pipette has to be totally filled up to the upper end. Hold the capillary pipette for hygienic reasons with a capillary holder or tweezers. Alternatively the blood can be taken by a micro pipette. Please take care that the exact volume of 10 µl is taken, if micro pipettes or other capillary tubes are used. Please dilute the collected blood immediately to avoid coagulation.

#### Sample Dilution / Sample Stability

4. Put the blood-filled capillary pipette in the prepared tube with the dilution buffer. Alternatively 10 µl of blood can be directly added to the buffer, if it is collected with a micro pipette.
5. Close the tube and mix the sample vigorously by hand for ~10 seconds, so that the blood from the capillary pipette is set free and the sample and dilution buffer are thoroughly mixed.
6. Let the diluted sample stand for ~1 min.
7. The sample can be immediately used or stored for max. 8 h.

NOTE: EDTA-, citrate- or heparin-blood can also be used, which has to be diluted with the provided buffer prior to testing.

### DIRECTIONS FOR USE


After opening the aluminium pouch the test should be immediately performed, as it is sensitive to humidity.

1. Open the aluminium pouch and take out the test strip at the end on which CRP is printed. Avoid touching the white reaction area.
2. Open the bottle with the diluted sample and put the test strip with the other end into the liquid. Make sure, that the test is immersed only to the MAX mark. A direct wetting of the reaction area by spilling the liquid or a too deep immersion should be avoided.
3. Leave the test strip in the sample fluid for ~20 seconds, until the pink liquid front is visible in the reaction area
4. Take out the test strip and put it onto a flat, non-absorbent surface. Alternatively, the test strip can be left in the bottle. Start the timer.
5. Wait for the appearance of the colored lines. Interpret the result after 5 min. Please stick to this timing, as only then a correct semi-quantitative conclusion can be made.

### INTERPRETATION OF RESULTS

Please make sure that the test is interpreted after exactly 5 minutes.

For the interpretation of the results please observe the pattern of test lines, which has formed in the white membrane area.

### Positive Result

If more than one line (control line C) is visible below the handle the test result is positive. This means that CRP is present in the sample.

The formation of test result lines (T1, T2 and T3) is dependent on the concentration. The more lines are visible, the higher is the CRP concentration in the sample.

- 1 red test result line (T1) and 1 red control line (C):  
Above a CRP concentration of ~ 10 mg/l a test line appears below the control line. With rising CRP concentrations the color intensity of the test line increases.
- 2 red test result lines (T1 and T2) and 1 red control line (C):  
Above 40 mg/l a second weaker colored test line appears above the first test result line.
- 3 red test result lines (T1, T2 and T3) and 1 red control line (C):  
Above 80 mg/l a weak third test line appears above the first two test result lines.

The test line formation, which depends on the CRP concentration, is depicted in the diagram above. You can see that the number of test result lines increases with rising CRP concentration. The color intensity of the lines is weak at their first appearance, but later it increases with rising CRP concentrations.

Please note that the three lines have different color intensity. For the depicted concentration range you can see that the color intensity of the lowest line is the strongest. The middle and especially the upper test result line have a clearly lower color intensity.

In comparison to the coloration of the test result lines the coloration of the control line is independent of the CRP concentration and shows a constantly intense red coloring.

### Negative Result

Only the red control line appears. The CRP concentration in the sample is < 10 mg/l.

### Invalid

The intensely red colored control line (C) in the upper area of the membrane did not appear. In this case the test is invalid, even if the test result line(s) are visible. The test should be repeated.

### Note

Please interpret the test only once after 5 minutes. It is normal for immune chromatography and for the kinetics of rapid tests that the intensity of all lines increases over time. The correct, semi-quantitative interpretation of the test results cannot be guaranteed after a longer reading time, as lower concentrations would be classified too high.

### EXPECTED VALUES

The CRP plasma level rises 6-8 hours after the acute event like a bacterial infection or a trauma. They reach the maximum level after about 48 hours after the event.

The level rapidly decreases after a stop of the causing events, with a half life of CRP of 48 hours.

Normally the degree of inflammation and the inflammation activity affect the amount of CRP increase. Values between 10 and 40 mg/l can be found during weak inflammations like local bacterial infections, abscesses, weak traumas, often in malign tumors, most viral diseases and more. Up to 100 mg/l point to severe diseases with inflammation, which normally require instant medical attention. Values above 100 mg/l can be found during bacterial sepsis or surgeries.

### LIMITATIONS

- CRP is no specific marker for a special disease pattern. Like all in-vitro-diagnostic tests the result should never be interpreted isolated, but in combination with other clinical data. Often the CRP increase precedes the actual symptoms, so also the temporal context should be taken into account.
- The intra-individual variation of CRP values is rather high. Normally, values > 10 mg/l are to be considered as increased values.
- The test cannot be used for the risk assessment of coronary heart diseases. For this purpose it is not sensitive enough.

### PERFORMANCE CHARACTERISTICS

The test was calibrated against the international human CRP reference standard WHO 85/506. The detection limit of the "DIAQUICK" CRP Dipstick is calibrated in a way that a concentration of 10 mg/l in undiluted samples results in a positive red line after dilution with the provided buffer.

If the undiluted sample contains more than 40 mg/l CRP, a second red line appears. Concentrations of more than 80 mg/l CRP result in a third red line appearing.

The concentration ranges could be confirmed in comparison with a quantitative reference test.

Hook effects could not be detected up to a concentration of 2.000 mg/l CRP. At this concentration, the test showed results of > 80 mg/l.

### QUALITY CONTROL

Good laboratory practice (GLP) recommends the use of controls to verify the proper test performance.

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